

Triphenyltetrazoliumchlorid p.A. Triphenyltetrazolium chloride (TTC), analytical grade	Kat.Nr./Cat.No. :	37130
--	--------------------------	--------------

Parameter parameter	Methode method	Spezifikation specification
Summenformel formula		C ₁₉ H ₁₅ ClN ₄
Molekulargewicht molecular weight		334.79
Aussehen appearance		cremeweißes Pulver cream-white powder
Reinheit (%) purity	HPLC	min. 99.0
λ_{max} (nm)	0.001 % in H ₂ O	243 – 251
A 1 cm / λ_{max}	1 % in H ₂ O	min. 670
Identität identity	IR	entspricht corresponds
Schmelzpunkt (°C) melting point	Zersetzung decomposition	min. 235
Trocknungsverlust (%) loss on drying		max. 2.0
Lagerung (°C) storage		+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.