

SERVA Unstained SDS PAGE Protein Marker
6.5 – 200 kDa, Liquid Mix

Cat.No. : 39215

| Protein Components | Molecular weight (kDa) |
|--|-----------------------------------|
| Myosin (bovine muscle) | 200 |
| β-Galactosidase (E. coli. recombinant) | 116 |
| Albumin (bovine) | 67 |
| Albumin (egg) | 45 |
| Carbonic anhydrase (bovine erythrocytes) | 29 |
| Trypsin inhibitor (soybean) | 21 |
| Lysozym (egg) | 14.3 |
| Aprotinin (bovine lung) | 6.5 |
| Storage (°C) | -15 to -25 |

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.