

---

<b>Ammoniumsulfat</b> Ammonium sulfate molecular biology grade	<b>Kat.Nr./Cat.No. :</b>	<b>39753</b>
---	--------------------------	--------------

---

<b>Parameter</b> parameter	<b>Methode</b> method	<b>Spezifikation</b> specification
<b>Summenformel</b> formula		(NH <sub>4</sub> ) <sub>2</sub> SO <sub>4</sub>
<b>Molekulargewicht</b> molecular weight		132.2
<b>Aussehen</b> appearance		weiße Kristalle white crystals
<b>Gehalt (%)</b> assay	Titration	min. 99.0
<b>pH</b>	10% in H <sub>2</sub> O	4.5 - 5.5
<b>Identität</b> identity	IR	entspricht corresponds
<b>DNase/ RNase Aktivität</b> DNase/ RNase activity		nicht nachweisbar none detected
<b>Lagerung (°C)</b> storage		+15 bis / to +30

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.