

|                                                                                                                           |                          |              |
|---------------------------------------------------------------------------------------------------------------------------|--------------------------|--------------|
| <b>Glycerin pflanzlich, 87 % (w/w), für Molekularbiologie</b><br>Glycerol from plant, 87 % (w/w), molecular biology grade | <b>Kat.Nr./Cat.No. :</b> | <b>39788</b> |
|---------------------------------------------------------------------------------------------------------------------------|--------------------------|--------------|

| Parameter<br>parameter                                 | Methode<br>method            | Spezifikation<br>specification                             |
|--------------------------------------------------------|------------------------------|------------------------------------------------------------|
| <b>Summenformel</b><br>formula                         |                              | C <sub>3</sub> H <sub>8</sub> O <sub>3</sub>               |
| <b>Molekulargewicht</b><br>molecular weight            |                              | 92.09                                                      |
| <b>Aussehen</b><br>appearance                          |                              | farblose, viskose Flüssigkeit<br>colorless, viscous liquid |
| <b>Brechungsindex</b><br>refractive index              | n <sub>D</sub> <sup>20</sup> | 1.449 – 1.455                                              |
| <b>DNase/ RNase Aktivität</b><br>DNase/ RNase activity |                              | nicht nachweisbar<br>none detected                         |
| <b>Ursprung</b><br>origin                              |                              | pflanzlicher Ursprung<br>plant origin                      |
| <b>Lagerung (°C)</b><br>storage                        |                              | +15 bis / to +30                                           |

|                                                                                                                                                                                |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p><b>Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.</b><br/>We do not guarantee that the product can be used for a special application.</p> |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|