

---

<b>SERVA Tris-Glycin/LDS Probenpuffer (4x)</b> SERVA Tris-Glycine/LDS sample buffer (4x)	<b>Kat.Nr./Cat.No. :</b> 42525
---	--------------------------------

---

<b>Parameter</b> parameter		<b>Spezifikation</b> specification
<b>Aussehen</b> appearance		rotviolette Lösung red violet solution
<b>pH</b>	1x	7.6 – 7.8
<b>Zusammensetzung</b> composition		entspricht corresponds
Glycerin glycerol		40%
Lithiumdodecylsulfat lithium dodecyl sulfate		4%
Ficoll® 400		4%
EDTA		1.89 mmol/ l
Bromphenolblau bromophenol blue		0.025%
Phenolrot phenol red		0.025%
Triethanolamin gepuffert triethanolamine buffered		
<b>Lagerung (°C)</b> storage		+2 bis/ to +8

**Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.**  
We do not guarantee that the product can be used for a special application.