
Towbin Puffer (10x) Towbin Buffer (10x) for Western Blotting	Kat.Nr./Cat.No. :	42558
--	--------------------------	--------------

Parameter parameter	Spezifikation specification
Aussehen appearance	klare farblose Lösung clear colorless solution
pH	8.4 – 8.8
Anwendung application	für Westernblotting for Western blotting
Zusammensetzung composition	
TRIS (# 37180)	0.25M
Glycin / Glycine (# 23390)	1.92M
Arbeitsverdünnung working solution	To 100 ml of concentrate add 200 ml methanol and 700 ml deionized water
	0.025 M Tris 0.192 M Glycine 20 % Methanol
Referenz reference	Towbin, H., Staehelin, T. and Gordon, J. (1979) Electrophoretic transfer of proteins from polyacrylamide gels to nitrocellulose sheets: procedure and some applications. Proc.Nat.Acad.Sci. USA 76 , 4350-4354.
Lagerung (°C) storage	+15 bis/ to +30

Die Eignung des Produktes für spezielle Anwendungszwecke wird nicht zugesichert.
We do not guarantee that the product can be used for a special application.