

# Omni Tissue Master 125

*Value priced, complete kit*



The **Omni Tissue Master 125** is an extremely economical homogenizer that is well suited for customers that are occasionally homogenizing small volume samples with little variation in sample size. The Tissue Master 125 comes complete with a 7 mm stainless steel probe.

The **Omni Tissue Master 125 Homogenizer Complete Kit** includes a 5 mm, 7 mm and 10 mm stainless steel generator probes, specifically designed for processing volumes of <math><200\ \mu\text{L}</math> to 100 mL. These probes can handle a wide variety of different sample sizes, making the **Tissue Master** a great choice for a laboratory with many applications, but a limited budget.

## Features

- > Powerful 125 watt motor
- > Light weight ergonomic design
- > Value priced, complete with 7 mm probe
- > Complete kit includes 5 mm, 7 mm and 10 mm stainless steel generator probes

## VALUE PRICED HOMOGENIZING KIT

WITH 5 MM PROBE    **ORDER NO. TM125-115-05**  
WITH 7 MM PROBE    **ORDER NO. TM125-115**  
WITH 10 MM PROBE    **ORDER NO. TM125-115-10**



**Processing Range:** <math><200\ \mu\text{L}</math> to 100 mL  
**Speed:** Variable from 500 to 35,000 rpm  
**Power Rating:** 125 watt  
**Weight:** 535 g (18.8 oz.)  
**Sound Level:** <math><72\ \text{db}</math>  
**Warranty:** 1 year motor warranty  
**Standards Approval:** CE Approved



## Omni Stainless Steel generator probes

OMNI TISSUE MASTER 125  
HOMOGENIZER COMPLETE KIT

(WITH 5, 7 AND 10 MM PROBES) ALSO AVAILABLE INDIVIDUALLY

ORDER NO. TM125-PKG



## Patented Omni Tip™ plastic disposable generator probes

Proprietary **Omni Tip™ Plastic Probes** are available in both soft tip (deaggregation) and hard tissue (frozen) versions. Engineered from durable plastics, Omni Tips™ can be cleaned and reused, yet are economical enough that they can be disposed of after each use to eliminate cross contamination.

*The Omni Tips' patented design eliminates the hassle of disassembly and maintenance of traditional stainless steel generators.*

### A selection of our best-selling plastic probes:

SOFT TISSUE PLASTIC GENERATOR PROBES		ORDER NO.
Omni Tip™ Plastic Generator Probes, Soft Tissue	25 pack	30750
Omni Tip™ Plastic Generator Probes, Soft Tissue	50 pack	32750
Omni Tip™ Plastic Generator Probes, Soft Tissue	100 pack	34750
HARD TISSUE PLASTIC GENERATOR PROBES		ORDER NO.
Omni Tip™ Plastic Generator Probes, Hard Tissue	25 pack	30750H
Omni Tip™ Plastic Generator Probes, Hard Tissue	50 pack	32750H
Omni Tip™ Plastic Generator Probes, Hard Tissue	100 pack	34750H



Omni International, Inc.  
935C Cobb Place Blvd.  
Kennesaw, GA 30144  
800.776.4431 • 770.421.0058



> For more information — **omni-inc.com**