





VECTASHIELD® Vibrance™ Antifade Mounting Medium with DAPI

Product Images



Short Description

VECTASHIELD Vibrance Antifade Mounting Media are completely new formulations of curing (hard setting) antifade mounting media for immunofluorescence applications. These newest additions to our VECTASHIELD portfolio were developed with customer feedback to improve on parameters such as ease of use and retention of specific staining intensity over time. VECTASHIELD Vibrance Antifade Mounting Media are new tools to help investigators see more and do more with each experiment.

VECTASHIELD Vibrance offers you:

- Superior antifade/anti-photobleaching properties across the spectrum
- Compatibility with commonly used fluorophores
- View sections one hour after mounting
- No tone or autofluorescent background, even after curing
- Room temperature storage of slides with extended archiving time
- Minimal bubble formation, even after several weeks storage
- Curing formulations with choice of counterstain (DAPI) or no counterstain
- 12 month product expiration date

Plus the same features that have made the VECTASHIELD brand successful:

- Ideal refractive index
- Easy to use
- Ready to use, requires no warming
- No sealing of coverslips required (curing formulations)

Additional Information

Unit Size	2 ml, 10 ml
Applications	Immunofluorescence, In situ hybridization, Cellular Imaging
Concentration	0.9 µg/ml
Recommended Storage	Refrigerate
Mounting	Aqueous (Hardening)
Antifade	Yes
Counterstain	DAPI (blue fluorescence)
Cellular Stain	Nucleus


VECTASHIELD® Vibrance™ Antifade Mounting Medium with DAPI

Features:

- Superior antifade/anti-photobleaching properties across the spectrum
- Compatibility with commonly used fluorophores
- View sections one hour after mounting
- No tone or autofluorescent background, even after curing
- Room temperature storage of slides with extended archiving time
- Minimal bubble formation, even after several weeks storage
- Curing formulations with choice of counterstain (DAPI) or no counterstain
- 12 month product expiration date

Plus the same features popular with the VECTASHIELD brand:

- Ideal refractive index
- Easy to use
- Ready to use, requires no warming
- No sealing of coverslips required (curing formulations)



VECTASHIELD® Vibrance™ Antifade Mounting Medium with DAPI

Features:

- Superior antifade/anti-photobleaching properties across the spectrum
- Compatibility with commonly used fluorophores
- View sections one hour after mounting
- No tone or autofluorescent background, even after curing
- Room temperature storage of slides with extended archiving time
- Minimal bubble formation, even after several weeks storage
- Curing formulations with choice of counterstain (DAPI) or no counterstain
- 12 month product expiration date

Plus the same features popular with the VECTASHIELD brand:

- Ideal refractive index
- Easy to use
- Ready to use, requires no warming
- No sealing of coverslips required (curing formulations)



