

VectaMount® Mounting Media & Substrate Compatibility Table

Substrate	Cat. No	VectaMount Express Mounting Medium (Non-Aqueous)	VectaMount Permanent Mounting Medium (Non-Aqueous)	VectaMount AQ Aqueous Mounting Medium
Horseradish Peroxidase (HRP)				
ImmPACT® DAB (brown)	SK-4105	•	•	•
ImmPACT DAB EqV (brown)	SK-4103	•	•	•
DAB (brown)	SK-4100	•	•	•
DAB-Ni (gray-black)	SK-4100	•	•	
ImmPACT AEC (red)	SK-4205			•
ImmPACT AMEC Red (red)	SK-4285			•
AEC (red)	SK-4200			•
TMB (blue)	SK-4400		•	
ImmPACT VIP (purple)	SK-4605	•	•	
Vector® VIP (purple)	SK-4600	•	•	
ImmPACT SG (blue-gray)	SK-4705	•	•	•
Vector SG (blue-gray)	SK-4700	•	•	•
ImmPACT NovaRED® (red)	SK-4805	•	•	
Vector NovaRED (red)	SK-4800	•	•	
Alkaline Phosphatase (AP)				
ImmPACT Vector Red (magenta)	SK-5105	•	•	•
Vector Red (magenta)	SK-5100	•	•	•
Vector Black (brown-black)	SK-5200		•	
Vector Blue (blue)	SK-5300	•	•	
BCIP/NBT (indigo)	SK-5400	•	•	•

For More information please visit:
vectorlabs.com