

Beer Bottle Shaker SM 30 B control



Description

The laboratory shaker SM 30 B control with a horizontal shaking motion up to 240 rpm and a stroke of 50 mm is ideally suited for Gushing Tests. It also has a touch display with programmable control. Up to 6 horizontal bottles (0,33 l or 0,5 l) can be fixed on the shaker with the suitable rack system. With a loading capacity of max. 30 kg it can be supplied with a two- or three-storey top frame.

Basic Equipment

Basic device with shaking plate and rubber mat, **without rack system**
With touch display and programmable control

Loading Capacity

Up to 6 horizontal bottles (0,33 or 0,5l) per storey

1-storey, 2-storey or 3-storey rack system can be selected as accessories

[Loading details](#) against shaking speed.

Options

[USB 2.0 interface](#)

[RS 232 interface](#)

[Control interface](#)

[Safety Package](#)

[Defined stop](#)

Details

| | |
|--------------------------|--|
| Motion: | reciprocating (horizontal) |
| Shaking platform: | 560 x 400 mm |
| Max. load: | 30 kg |
| Shaking speed: | up to 240 rpm, in steps of 5 |
| Stroke: | 50 mm for Gushing Test (No. 0052 216) |
| Runtime: | programmable / continuous |
| Order No.: | 6102 000 |
| Accessory: | 1-storey top frame 2-storey top frame 3-storey top frame |

Technical specifications

| | |
|--------------------------------|---|
| Electrical supply: | 230 V or 115 V, 50 / 60 Hz please indicate in case of order |
| Enclosure protection: | IP 21 |
| Heat emission: | approx. 8 W |
| Ambient temperature: | 5°C to 50°C |
| Relative humidity: | ~ 85 % |
| Dimensions (w x d x h): | 680 x 610 x 160 mm |
| Weight: | 33 kg |