



Infrared Heater IRB 2



Description

The compact laboratory heater IRB 2 with power control and a heating area of 100 x 100 mm is optimally suited for fast heating tasks without open flame up to max. 900°C.

Basic equipment

Basic device incl. integrated power control, support cross

Details

| | |
|----------------------------------|----------------------------------|
| Max. temperature: | 900°C |
| Heating power: | 800 W |
| Heating area: | 100 x 100 mm |
| Power control: | stepless control between 0- 100% |
| Max. ambient temperature: | 30°C |
| Order number: | 6063 000 |
| Accessory: | Support rod for stand connection |

Technical specifications

| | |
|--------------------------------|----------------------------------------------|
| Electrical supply: | 230 V, 50 / 60 Hz (115 V version on request) |
| Dimensions (w x d x h): | 150 x 150 x 170 mm |
| Weight: | 2,5 kg |