




Infrared Heater IRB 6



Description

With the series heater IRB 6 the samples can be heated on 6 separate heating areas of 60 x 60 mm each up to 700°C without open flame. Each heater can be individually used and controlled.

Basic equipment

6 x IRB 1 with support crosses, 1 x power control unit LR 6

Details

Max. temperature:	each 700°C
Heating power:	6 x 200 W
Heating area:	6 x 60 x 60 mm
Power control:	6 stepless controls between 0 - 100 %
Max. ambient temperature:	30°C
Order no.:	9066 000
Accessory:	Support base Horizontal PVC rods for support base Clamps D50 (NS 45) Clamps D32 (NS 29) Clamps D20 (NS 19) Clamps D25 (NS 15)

Technical specifications

Electrical supply:	230 V, 50 / 60 Hz (115 V version on request)
Heating power:	max. 1200 W
Dimensions (w x d x h):	697 x 305 x 122 mm
Weight:	9 kg