


Microplate Shaker TiMix 2



Description

The Microplate Shaker TiMix 2 especially for small samples in microplates, deepwell-plates or other plates in microplate format.

Basic Equipment

Basic device incl. basic platform with holders for 4 microplates

Loading Capacity

4 microplates

(deepwell-plates, or other plates in microplate format)

You will find further loading possibilities [here](#)

Details

Motion:	exact orbital motion
Shaking platform:	200 x 295 mm for 4 microplates
Max. load:	2 kg
Shaking speed:	100 - 1400 min ⁻¹ , steplessly adjustable
Stroke:	3 mm
Timer:	0 - 120 min / continuous
Order No.:	6110 000
Accessory:	Shaking platform KM Clamping strip KM

Technical specifications

Electrical supply:	230 V or 115 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Heat emission:	approx. 6 W
Ambient temperature:	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	220 x 410 x 110 mm
Weight:	8,5kg