


Three-Dimensional Shaker TL 10



Description

The Three-Dimensional Shaker TL 10 works with slow, three-dimensional motion. The angle of tilt is adjustable between 0° and 9° and the shaking speed between 1 and max. 23 rpm. The shaking platform of 340 x 340 mm can be loaded with max. 10 kg.

Basic Equipment

Basic device incl. basic tray and rubber mat

Details

Motion:	three-dimensional orbital motion
Shaking platform:	340 x 340 mm
Max. load:	10 kg
Shaking speed:	1 - 23 rpm, steplessly adjustable
Angle of tilt:	0° to 9° einstellbar
Timer:	0 - 120 min / continuous
Order No.:	6117 000
Accessory:	Multi-storey rack system WS / TL

Technical specifications

Electrical supply:	230 V or 115 V, 50 / 60 Hz please indicate in case of order
Enclosure protection:	IP 21
Heat emission:	approx 6 W
Ambient temperature::	5°C to 50°C
Relative humidity:	~ 85 %
Dimensions (w x d x h):	355 x 450 x 205 mm
Weight:	18,5 kg