

Tube Roller RM 2



Description

Gentle tilting and rotating motion for uniform mixing

Basic Equipment

Basic device with 6 removable plastic rolls and spill tray

Loading Capacity

All standard tube sizes, flasks and universal vessels up to 80 mm in diameter..

Details

Motion:	rotating and tilting motion
Shaking platform:	200 x 300 mm
Max. load:	2 kg
Shaking speed:	5 - 80 rpm, in steps of 5
Angle of tilt:	3° fixed
Stroke at end of rolls:	15 mm
Dimensions of rolls:	300 x Ø 32 mm
Number of rolls:	6 pieces
Runtime:	0 - 99 h 59 min / continuous
Order No.:	6138 000
Accessory:	Overhead Rack RM 2

Technical specifications

Electrical supply:	100 - 240 V, 50 / 60 Hz
Enclosure protection:	IP 21
Heat emission:	max. 0,5 - 2 W
Ambient temperature:	5°C to 45°C
Relative humidity:	~ 85 % (non-condensing)
Dimensions (w x d x h):	244 x 465 x 115 mm
Weight:	7,3 kg