



Product Catalogue Medication Safety



Comprehensive solutions
for preparing, storing
and distributing medicines

Improving patient safety

The WIEGAND® system offers comprehensive and flexible solutions to increase medication and, ultimately, patient safety. Parallel to providing an improvement in the quality of care, the simple, modular system components are applicable to inpatient and ambulatory healthcare environments. Our systems reduce the work load of nurses and pharmacy assistants, resulting in a reduction in medication errors while substantially decreasing costs.



Sustainable know-how

WIEGAND® has been studying medication processes for over forty years. Through close collaboration with care specialists, pharmacists and physicians across Europe, we continually develop our concepts to find new and practical solutions for our customers' ever changing requirements.

WIEGAND®'s concepts increase efficiency and improve the working process of nurses, thus leading to a decrease in medication errors.

The WIEGAND® products are durable and have an outstanding cost-benefit ratio, which effectively facilitates the lowering of health care costs.



Reducing liability risks

The frequency of medication errors occurring between writing a prescription to the administration of medicines is alarmingly high. For years, various statistics from the USA and Europe refer to this harmful and costly problem. Only by integrating these factors into a health care institution's risk management process, the reduction of medication errors will become an achievable strategic goal.

Medication errors are costly


According to different estimates, 60-80 % of all errors in the health care sector can be attributed to human failure. The consequences of medication errors are disturbing. In an extensive report by D.W. Bates and colleagues, the following extrapolations were made:

- Estimated annual costs of Adverse Drug Events (ADE) for a 700-bed hospital: \$ 5.6 million
- ADEs were recorded in 6.5 % of all patients
- Overall costs: \$ 2,595 per ADE
- ADEs lead to 2.2 additional days of hospitalisation
- Preventable ADEs (28 % of all ADE) annual cost: \$ 2.8 million
- Costs for each preventable ADE: \$ 4,685
- Preventable ADEs result in additional 4.6 days of hospitalisation
- 57 % of all ADEs were judged significant (30 % serious, 12 % life-threatening and 1 % fatal)

«These estimates do not include costs of injuries to patients, malpractice costs, or the costs of less serious medication errors or admissions related to ADEs. These results suggest that hospitals can justify devoting additional resources to develop systems that reduce the number of preventable ADEs not only to improve patient care but also to reduce ADE-related expenses.»


Source: Cost of adverse drug events; D.W. Bates et. al., JAMA, Jan 22/29, 1997 – Vol277, No. 4

MEDMARX® (USA) analyzed 21,084 cases with ADEs and highlighted the consequences of the harmful cases.



Allocation of medication errors

The following data from MEDMARX® (USA) are based on 218,347 reported cases and the frequency of the different sources of error was analyzed.



Source: Hicks, R.W., Santell, J.P., Cousins, D.D., and Williams, R.L. (2004). MEDMARXSM 5th Anniversary Data Report: A Chartbook of 2003 Findings and Trends 1999–2003. Rockville, MD: USP Center for the Advancement of Patient Safety.

The Swiss Foundation for Patient Safety carried out a survey in 2006 on patient safety in Swiss hospitals. The six most pertinent sources of medication errors (number of entries x weighting) are listed below:

1. Error while preparing medicines
2. Administration of medicines to the wrong patient
3. Wrong dosage (calculation error in potencies of 10)
4. Bad legibility of drug prescriptions
5. Faulty/missing information of drug therapy with/after discharge
6. Transfer errors when copying prescriptions

Source: Stiftung für Patientensicherheit: Problemfelder (Hot-Spots) in der Patientensicherheit; Schweizerische Ärztezeitung 2008;89: 24

Day medication prepared according to administration times

In order to guarantee a smoothly running day on the ward, the medication must be well planned and prepared in advance. Based on correct care documentation, all drugs for a whole day are put into the appropriate compartments of the medicine dispensers according to the prescribed administration times. Control of the prepared medication and their administration times is thus secured.



Avoiding errors

- Check the patient's identity
- Ensure that dosage calculations are checked independently by another health care professional before the drug is administered
- Make sure that prescription, drug, and patient match
- Make certain that the medication is given at the correct time
- Minimise interruptions during drug rounds

Source: Medication errors, DJP Williams, J R Coll Physicians Edinb 2007; 37:343–346

Reducing errors through the four-eyes principle

Drug therapies (medications) for all patients should be ideally prepared during a calm period. This task requires the highest level of concentration and should never take place at a patient's bedside. Before administering medicines to a patient, drugs prepared in a medication dispenser should be controlled by a second person to aid in reducing medication errors.

Medikament (MediDispenser)	07:00	08:00	12:00	18:00	21:00	Mo	Di	Mi	Do	Fr	Sa	So	Informationen	Bestand
ACTOS 30 Tablette													Tablette, rund, weiss, nicht teilbar Ch: 170492 (1.6.09)	100.0 → 08.06.10
EUTHYROX 125 mg Tablette													Tablette, rund, weiss, teilbar, nüchtern einnehmen Ch: 1468325 (1.6.04)	82.0 → 20.09.10
ISMN 40 retard Heumann Kaps.													Harthaus, rötlich, rot-braun, nicht teilbar Ch: 343514 (1.1.08)	63.0 → 02.08.10
SILIMARTIN Kaps.													Wiedelkaps, Olibon, braun, nicht teilbar Ch: 9910982 (1.4.05)	28.0 → 06.07.10
METOPROLOL 200 retard													Retard-Tablette, rund, weiss, verstellbar Ch: 000432 (1.6.07)	60.0 → 04.08.10
MAGNESIOCARO Sachet													Pulver, in Wasser auflösen Ch: 000350 (1.1.08)	38.0 → 15.06.10
SUPRADYN Tabs.													Braunspachtel, teilbar, in Wasser auflösen Ch: 38200002 (1.8.09)	13.0 → 19.06.10
METFORMIN AL 850													Farmakoblen, oval, weiss, halbkreisförmig Ch: 980182 (1.6.07)	45.0 → 21.07.10

Euthyrox eine halbe Stunde vor Frühstück einnehmen

Wichtig

Verabreichte von:
 Unterschrift:
 Datum:
 Zuständige Abteilungsleiter:
 HEMERAN GbR
 Bedarfsmenge:
 (Porzellan Tablette (1 OP a 20 Tablette)
 max. Dosis: 100mg/Tag

Dr. Pierre Lamaze


Geburtsdatum: 12.03.1926
Haus Sonnenhof
Seestrasse 23
8002 Zürich
Station 3B
Zimmer 204
Arzt
Dr. Pierre Lamaze

gerichtet: 07.06.10
gerichtet von: B. Smeets
gültig: 07.06.10 - 14.06.10
Lohnapothek
Bahnhofstrasse 8
8000 Zürich



WIEGAND® MediDispenser

- The sliding lid opens in one direction only
- This secures the administration sequence
- 4 to 12 flexibly adjustable compartments (1 for each administration)
- Hygienic and economical (semi disposable)
- Tight sliding lid



Adaptable medication processes

Medication processes differ substantially in terms of the logistical framework in which they are applied. A nursing home, operating with centrally or even externally prepared medications has different workflows to the medical or surgical department of a hospital. The more acute the treatments are, the shorter are the processes and operating times and thus the higher are the corresponding needs for quality in care and the costs of an incorrect treatment.



Decentralised ward- and satellite pharmacies – cost-efficient, organised and space-saving

- Provide an excellent overview and fast access
- Low stocks and good availability due to Kanban principle (always 2 units/OP)
- Application of “first in - first out” principle of storage, provides an effective control of expiry dates and lots of the drugs
- The central pharmacy steers the quantity of supplies and determines the permitted pharmaceutical products (cost control)
- Individual drug costs are charged at time of consumption
- Fast turnover of drugs without temporary storage in BIN System



WIEGAND® Medication
Prepare Station for 40 patients

Improving efficiency

When using the Day System for the medication process, the nurses prepare the prescribed drugs – after the doctor's visits – for the next 24 hours. Preparing for a 7-day period (Week System) significantly optimises operational flows in retirement- and nursing homes. To keep patients in a «normal condition», a long-term therapy with a regular, stabilizing drug therapy is applied. Therefore medications are preferably prepared for a whole week.

Professional WIEGAND® Medication Prepare Stations

Empirical research illustrates that when using a professionally furnished Medication Prepare Station, medicines for 20 patients can be prepared in approximately one hour. Decentralised ward- and satellite pharmacies with a Medication Prepare Station are set up to manage between 40 and 100 patients.



WIEGAND® Distribution Tray
8 patients - 24 hours

WIEGAND® Week Tray
1 patient - 7 days

BIN System combined with WIEGAND® Medication Distribution System

When using BIN Systems, a patient's drugs are kept in individual containers. Although these systems need additional storage space, they are quite commonly used and can be easily combined with the WIEGAND® Medication Distribution System in order to be able to prepare medications in advance. In this way, patient safety can be substantially increased.



WIEGAND® medication concepts

- Enhance existing systems, in order to increase patient safety
- Cost efficient investment to reduce the expense of medication errors
- Replace inefficient and outdated systems to reduce liability risks
- Sustainable, modular system leading to greater flexibility in application and investment security

MediDispenser | the «small» ones

Article No. 01/09

Description MediDispenser 01

- Details**
- MediDispenser No. 01 with Sliding-Lid No. 09
 - also available as imprinted Week-Set
 - 165x50x20 mm

Accessories

- 510.14-8xx Label with dispensing periods
- 510.18 Label, extra-small
- 291 Week-Tray 01
- 93T-01 WeekBox 01
- 93-01 Week-Pouch 01
- 295 Sliding-Lid-Opener




Article No. 11/19

Description MediDispenser 11

- Details**
- MediDispenser No. 11 with Sliding-Lid No. 19
 - 225x50x20 mm

Accessories

- 510.14-8xx Label with dispensing periods
- 510.13 Label for patients name and medication
- 510.13Wxx Label for medication with week days
- 510.13A4 Label "No. 13" for imprinting with MS-Word-form
- 591 Week-Tray 11
- 93T-11 WeekBox 11
- 91 Week-Container 11



Article No. 131/139

Description MediDispenser 131

- Details**
- MediDispenser No. 131 with Sliding-Lid No. 139
 - 3 Divisions No. 132, divisible in max. 9 compartments
 - 188x42x37 mm

Accessories

- 510.14-8xx Label with dispensing periods
- 510.18 Label, extra-small
- 510.17xx Label for medication
- 17SD Label holder for Label 510.17xx
- 1191 Tray for MediDispenser 131



MediDispenser | the «medium» ones



Article No. 41/49

Description MediDispenser 41

- Details**
- MediDispenser No. 41 with Sliding-Lid No. 49-69-349
 - also available as imprinted Week-Set
 - 228x81x20 mm

Accessories

- 510.44xx Label with week days
- 510.46 Label, short strip
- 991 Week-Tray 41-61-341
- 193TQ-41-61-341.1bl WeekBox 41-61-341, blue design
- 191 Week-Container 41-61-341



Article No. 61/69

Description MediDispenser 61

Details

- MediDispenser No. 61 with Sliding-Lid No. 49-69-349
- also available as imprinted Week-Set
- 228x81x20 mm

Accessories

510.44xx	Label with week days
510.46	Label, short strip
510.64-7	Label for medication
66.x	Inlay-Dishes, coloured
66.Spez	Singe-Use-Inlay-Dish
67	Inlay-Blister with 4 compartments
991	Week-Tray 41-61-341
193TQ-41-61-341.1bl	WeekBox 41-61-341, blue design
191	Week-Container 41-61-341



Article No. 66.x

Description Inlay-Dishes, coloured

Details

- for MediDispenser No. 61/69
- x to be ordered in different colours:
 - No. 66.1 = Inlay-Dish, transparent-white
 - No. 66.2 = Inlay-Dish, transparent-yellow
 - No. 66.3 = Inlay-Dish, transparent-blue
 - No. 66.4 = Inlay-Dish, transparent-red
- 60x62x18 mm
- [PS]

Accessories

910.66.1 Basket for dish washer for Inlay-Dishes




Article No. 66.Spez

Description Singe-Use-Inlay-Dish

Details

- for MediDispenser No. 61/69
- 60x62x18 mm
- [PET]





Article No. 67

Description Inlay-Blister with 4 compartments

Details

- for MediDispenser No. 61/69
- 219x63x19 mm
- [PET]



Article No. 341/349

Description MediDispenser 341

Details

- MediDispenser No. 341 with Sliding-Lid No. 49-69-349
- 3 Divisions No. 342, divisible in max. 8 compartments
- also available as imprinted Week-Set
- 228x81x20 mm

Accessories

510.44xx	Label with week days
510.46	Label, short strip
991	Week-Tray 41-61-341
193TQ-41-61-341.1bl	WeekBox 41-61-341, blue design
191	Week-Container 41-61-341



MediDispenser | the «large» ones

Article No. 141/149

Description MediDispenser 141

- Details**
- MediDispenser No. 141 with Sliding-Lid No. 149
 - 3 Divisions No. 142, divisible in max. 6 compartments
 - 228x78x23 mm

Accessories

- 510.14-8xx Label with dispensing periods
- 510.54 Label with dispensing times
- 510.13 Label for patients name and medication
- 510.13Wxx Label for medication with week days
- 510.13A4 Label "No. 13" for imprinting with MS-Word form
- 491 Week-Tray 141




Article No. 151/159

Description MediDispenser 151

- Details**
- MediDispenser No. 151 with Sliding-Lid No. 159
 - 3 Divisions No. 152, divisible in max. 7 compartments
 - 221x102x45 mm

Accessories

- 510.14-8xx Label with dispensing periods
- 510.54 Label with dispensing times
- 510.18 Label, extra-small
- 891 Storage Tray for 8 MediDispensers 151



Article No. 171/179

Description MediDispenser 171

- Details**
- MediDispenser No. 171 with Sliding-Lid No. 179-189
 - 3 Divisions No. 172, divisible in max. 8 compartments
 - 260x98x23 mm

Accessories

- 510.44xx Label, transparent (week days)
- 510.48 Label, long strip
- 193TQ-171.2gn WeekBox 171, green design



Article No. 181/189

Description MediDispenser 181

- Details**
- MediDispenser No. 181 with Sliding-Lid No. 179-189
 - 3 Divisions No. 182, divisible in max. 8 compartments
 - 260x98x39 mm

Accessories

- 510.44xx Label, transparent (week days)
- 510.48 Label, long strip



Article No. 161/169

Description MediDispenser 161

- Details**
- MediDispenser No. 161 with Sliding-Lid No. 169-269
 - 4 Divisions No. 162, divisible in max. 12 compartments
 - 363x77x23 mm

- Accessories**
- 510.54 Label with dispensing times
 - 510.13 Label for patients name and medication
 - 510.13Wxx Label for medication with week days
 - 510.13A4 Label "No. 13" for imprinting with MS-Word form
 - 691 Week-Tray 161-261




Article No. 261/269

Description MediDispenser 261

- Details**
- MediDispenser No. 261 with Sliding-Lid No. 169-269
 - 4 Divisions No. 262, divisible in max. 12 compartments
 - 363x77x45 mm

- Accessories**
- 510.54 Label with dispensing times
 - 510.13 Label for patients name and medication
 - 510.13Wxx Label for medication with week days
 - 510.13A4 Label "No. 13" for imprinting with MS-Word form
 - 691 Week-Tray 161-261




Article No. 271/279

Description MediDispenser 271

- Details**
- MediDispenser No. 271 with Sliding-Lid No. 279-289
 - 4 Divisions No. 272, divisible in max. 12 compartments
 - 338x67x23 mm

- Accessories**
- 510.54 Label with dispensing times
 - 510.20 Label for week days and patients name
 - 1091 Storage Tray for 8 MediDispensers 271-281



Article No. 281/289

Description MediDispenser 281

- Details**
- MediDispenser No. 281 with Sliding-Lid No. 279-289
 - 4 Divisions No. 282, divisible in max. 12 compartments
 - 338x67x45 mm


- Accessories**
- 510.54 Label with dispensing times
 - 510.20 Label for week days and patients name
 - 1091 Storage Tray for 8 MediDispensers 271-281




MediDispenser | Inscription

Article No. 510.13**Description** Label for patient's name and medication

- Details**
- for MediDispenser No. 11, 141, 161, 261
 - bundle: 350 labels
 - 75x48 mm
 - adhesive, removable

**Article No.** 510.13Wxx**Description** Label for medication with week days

- Details**
- with imprint: «MONDAY..., SUNDAY»
 - sheet for 1 Week (bundle for 50 weeks)
 - for MediDispenser No. 11, 141, 161, 261
 - 75x48 mm
 - adhesive, removable
 - Replace **xx** in the article number with the desired language code:
510.13Wde = Deutsch, 510.13Wfr = Français,
510.13Wit = Italiano, 510.13Wen = English,
510.13Wnl = Nederlands

**Article No.** 510.13A4**Description** Label "No. 13" for imprinting with MS-Word form

- Details**
- A4-sheet: 10 labels (box with 2500 labels)
 - for MediDispenser No. 11, 141, 161, 261
 - 75x48 mm
 - MS-Word form for imprinting: www.wiegand.ch
 - adhesive, removable

**Article No.** 15/15SD**Description** Label holder

- Details**
- with Sliding-Lid No. 15SD
 - suitable for Label No. 510.13, No. 510.13Wxx and No. 510.13A4
 - to be plugged on a Distribution-Tray
 - 77x55 mm

**Article No.** 64.738**Description** Dispenser Marker


- Details**
- to mark MediDispenser, MediCups, Week-Sets etc.
 - with integrated eraser
 - dish washer proof ink



Article No. 510.46

Description Label, short strip

- Details**
- sheet with 11 labels (bundle of 20 sheets)
 - for MediDispenser No. 41, 61, 341
 - 191x11 mm
 - adhesive, removable



Article Nr 510.64-7

Description Label for medication, insert

- Details**
- for MediDispenser No. 61/69
 - 50 sheets with 7 labels each
 - paper: US-Legal form
 - 206x46 mm



Article No. 510.48

Description Label, long strip

- Details**
- sheet with 10 labels (bundle with 20 sheets)
 - for MediDispenser No. 171, 181
 - 235x12 mm
 - adhesive, removable



Article No. 510.17xx

Description Label for medication

- Details**
- A4 sheet with 8 labels
 - to be placed on Distribution Trays or in Storage Tray No. 1191
 - MS-Word form: www.wiegand.ch
 - 140x42 mm
 - Replace **xx** in the article number with the desired language code:
510.17de = Deutsch, 510.17en = English, 510.17fr = Français

Accessories

17SD Label holder for label 510.17xx



Article No. 17SD

Description Label holder for label 510.17xx

- Details**
- to be placed on Distribution Trays or in Storage Tray No. 1191
 - 150x45 mm
 - [PET]


Accessories

510.17xx Label for medication




Article No. 510.44xx**Description** Label with week days

- Details**
- with imprint: «MONDAY..., SUNDAY»
 - sheet of 2x7 labels, transparent
 - 41x8 mm
 - adhesive, removable
 - replace **xx** in the article number with the desired language code:
510.44de = Deutsch, 510.44fr = Français,
510.44it = Italiano, 510.44en = English,
510.44nl = Netherlands

**Article No.** 510.14-8xx**Description** Label with distribution periods

- Details**
- with imprint: «MORNING, NOON, EVENING, NIGHT, RESERVE»
 - sheet of 8x5 labels, transparent
 - 35x4 mm
 - adhesive, removable
 - replace **xx** in the article number with the desired language code:
510.14-8de = Deutsch, 510.14-8fr = Français,
510.14-8it = Italiano, 510.14-8en = English,
510.14-8nl = Netherlands

**Article No.** 510.54**Description** Label with distribution times

- Details**
- with imprint: «01:00, 02:00..., 24:00»
 - sheet of 35 labels, transparent
 - 13x13 mm
 - ideal for divisible MediDispenser
 - adhesive, removable

**Article No.** 510.18**Description** Label extra-small

- Details**
- A4 sheets of 96 labels
 - for MediDispenser or MediCups
 - 40x10 mm
 - MS-Word form for imprinting:
www.wiegand.ch
 - adhesive, removable

**Article No.** 510.20**Description** Label small

- Details**
- A4 sheets of 39 labels
 - for MediDispenser No. 271, 281
 - 60x20 mm
 - MS-Word form for imprinting:
www.wiegand.ch
 - adhesive, removable



Article No. 31

Description Distribution Tray simple

- Details**
- with 19 grooves
 - for MediDispensers
No. 01/09, No. 11/19, No. 131/139,
No. 41/49, No. 61/69 or No. 341/349
 - for 19 MediCups No. 71.x
 - 450x340 mm

Accessories
16 Universal Mounting Plate
for Distribution Trays

Note also available with
enforced base, No. 31B




Article No. 51

Description Distribution Tray extended

- Details**
- with 17 grooves
 - for MediDispensers
No. 01/09, No. 11/19, No. 131/139,
No. 41/49, No. 61/69 or No. 341/349
 - for 35 MediCups No. 71.x
 - 450x340 mm

Accessories
16 Universal Mounting Plate
for Distribution Trays

Note also available with
enforced base, No. 51B




Article No. 81

Description Distribution-Tray Universal,
small

- Details**
- with 10 grooves
 - for all MediDispenser
 - for 50 MediCups No. 71.x
 - 300x400 mm

Accessories
16 Universal Mounting Plate
for Distribution Trays




Article No. 241

Description Distribution-Tray Universal,
medium

- Details**
- with 16 grooves
 - for all MediDispenser
 - for 80 MediCups No. 71.x
 - 450x400mm

Accessories
16 Universal Mounting Plate
for Distribution Trays




Article No. 251

Description Distribution-Tray Universal, large

- Details**
- with 23 grooves
 - for all MediDispenser
 - for 115 MediCups No. 71.x
 - 600x400 mm

Accessories
16 Universal Mounting Plate for Distribution Trays



Trays | Storage Trays for MediDispensers

Article No. 891

Description Storage Tray for MediDispenser 151

- Details**
- for 8 MediDispensers No. 151/159
 - 450x450x47 mm



Article No. 1091

Description Storage Tray for MediDispenser 271-281

- Details**
- for 8 MediDispensers No. 271/279 or No. 281/289
 - 600x400x25 mm



Article No. 1191

Description Storage Tray for MediDispenser 131

- Details**
- for 8 MediDispensers No. 131/139 or 8x4 MediCups No. 71.x
 - with plug-in for MediCups
 - with 8 inscription grooves
 - 400x204x45 mm

Accessories
510.17xx Label for medication
17SD Label holder for label 510.17xx
299-1199 Sliding-Lid for Trays 291 and 1191



Week-Systems | Working aids

Article No. 295

Description Sliding-Lid-Opener

- Details**
- for Week-Tray No. 291 with MediDispenser No. 01/09
 - consists of plug-in unit No. 296 and pusher No. 297

Accessories
01/09 MediDispenser 01
291 Week-Tray 01




Week-Systems | Week-Trays with Sliding-Lids

Article No. 291

Description Week tray 01

Details - for 7 MediDispenser No. 01/09
- 400x204 mm

Accessories
299-1199 Sliding-Lid for Tray 291 and 1191
510.46 Label, short strip




Article No. 591

Description Week-Tray 11

Details - for 7 MediDispenser No. 11/19
- 400x245 mm

Accessories
599 Sliding-Lid for Week-Tray 11
510.46 Label, short strip




Week-Systems | Week-Trays

Article No. 991

Description Week-Tray 41-61-341

Details - for 7 MediDispensers No. 41/49,
No. 61/69 or No. 341/349
- 600x260 mm

Accessories
510.46 Label, short strip



Article No. 491

Description Week-Tray 141

Details - for 7 MediDispensers No. 141/149
- 600x260 mm

Accessories
510.46 Label, short strip



Article No. 691

Description Week-Tray 161-261

Details - for 7 MediDispensers No. 161/169
or No. 261/269
- 600x400 mm

Accessories
510.46 Label, short strip



Week-Systems | Week Container

Article No. 91

Description Week Container 11

Details - for 7 MediDispensers No. 11/19
- 265x177x58 mm

Accessories

510.46 Label, short strip
99 Sliding-Lid for Week-Container No. 91



Article No. 191

Description Week Container 41-61-341

Details - for 7 MediDispensers No. 41/49,
No. 61/69 or No. 341/349
- 265x177x87 mm

Accessories

510.46 Label, short strip
199 Sliding-Lid for Week-Container No. 191



Week-Systems | WeekBox

Article No. 193TQ-171.2gn

Description WeekBox 171, green design

Details - Week container for
7 MediDispensers No. 171/179
- foldable, washable
- translucent white, green design
- 267x160x102 mm

Accessories

510.95 Label for WeekBox




Article No. 193TQ-41-61-341.1bl

Description WeekBox 41-61-341, blue design

Details - Week container for 7 MediDispenser
No. 41/49, No. 61/69 or No. 341/349
- foldable, washable
- translucent white, blue design
- 138x84x229 mm

Accessories

510.95 Label for WeekBox
7-41en/49 Week set 41 with week days
7-61en/69 Week set 61 with week days
7-341en/349 Week set 341 with week days



MediDispenser Special | BlisterBox

Article No. 311

Description Dosys Dispenser

Details - Dispenser for blistered
unit dose drugs
for 7-14 days
- 82x150x150 mm



Article No.	311T		
Description	BlisterBox Large		
Details	<ul style="list-style-type: none">- Dispenser for blistered unit dose drugs for 7-14 days- foldable, washable- translucent white- 78x162x139 mm		
Article No.	312T		
Description	BlisterBox Compact		
Details	<ul style="list-style-type: none">- Dispenser for blistered unit dose drugs for 7 days- foldable, washable- translucent white- 78x128x107 mm		
Note	also available in white with imprint, No. 312W.1bl		
Article No.	500.240XO.075		
Description	Holder for labels and blister bags		
Details	<ul style="list-style-type: none">- for blister bags- to be plugged on a Distribution-Tray- can be used as Label holder on SwissModul Trays and Baskets- 75x40 mm		

Storage | PharmaBox

Article No.	111/119	
Description	PharmaBox small	
Details	<ul style="list-style-type: none">- Container to store drugs in defined quantities- Sliding-Lid No. 119- 130x72x30 mm	
Accessories	<ul style="list-style-type: none">112e Label, mini-card112R Mini-cards on rolls, white/yellow112a Label holder for mini-cards119weiss Sliding Lid, UV, for PharmaBox smallH121-111 Support plate for PharmaBoxes small on shelves	

Article No.	411/419	
Description	PharmaBox large	
Details	<ul style="list-style-type: none">- Container to store drugs in defined quantities- Sliding-Lid No. 419- 191x99x44 mm	
Accessories	<ul style="list-style-type: none">112e Label, mini-card112R Mini-cards on rolls, white/yellow112a Label holder for mini-cards419w Sliding Lid, UV, for PharmaBox largeH121-411 Support plate for PharmaBoxes large on shelves	

19

Tablet processing

Article No. 97

Description PillCrusher WIEGAND®

- Details**
- mortar for tablets
 - 2-pieces, Ø 70 mm, height 105 mm
 - with 1 MediCup, 25 ml

Accessories
71.x MediCup



Article No. 96T

Description MediSplitter WIEGAND®

- Details**
- tablet cutter
 - 85x35x25 mm
 - blade: stainless steel
 - with deposit compartment



Article No. 98T.1

Description MediCrusher WIEGAND®

- Details**
- mortar for tablets
 - ergonomically designed
 - with deposit compartment in the cover



Article No. 75.01

Description Dispenser Bag

- Article Info**
- for partitioned tablets or effervescent tabs
 - 45x60 mm



Drug storage | TRESTON

Description TRESTON Container, transparent

- Details**
- incl. Front Label and Label Protection
 - ordering:
 - No. 3010-1 = 92x300x82 mm
 - No. 3015-1 = 132x300x100 mm
 - No. 3020-1 = 186x300x82 mm
 - No. 4010-1 = 92x400x82 mm
 - No. 4015-1 = 132x400x100 mm
 - No. 4020-1 = 186x400x82 mm

- Accessories**
- 10-D Division matching with No. 3010, No. 4010
 - 15-D Division matching with No. 3015, No. 4015
 - 20-D Division matching with No. 3020, No. 4020
 - 10-E Front-Label matching with No. 3010, No. 4010
 - 15-E Front-Label matching with No. 3015, No. 4015
 - 20-E Front-Label matching with No. 3020, No. 4020
 - 06-E Label for Divisions



MediCups and accessories

Article No. 200.77.x

Description MediCups light 25 ml

- Details**
- medicine cup with graduation to 25 ml, CE 0413
 - x to be ordered in different colours:
 No. 200.77.1 = MediCup transparent-white
 No. 200.77.2 = MediCup transparent-yellow
 No. 200.77.3 = MediCup transparent-blue
 No. 200.77.4 = MediCup transparent-red
 No. 200.77.5 = MediCup transparent-green

Accessories
200.78 Cover for MediCups



Article No. 71.x

Description MediCups 25 ml

- Details**
- medicine cup with graduation to 25 ml, CE 0543
 - x to be ordered in different colours:
 No. 71.1 = MediCup transparent-white
 No. 71.2 = MediCup transparent-yellow
 No. 71.3 = MediCup transparent-blue
 No. 71.4 = MediCup transparent-red
 No. 71.5 = MediCup transparent-green

Accessories
72 Cover for MediCups
910.71.1 Dish Washer Basket with insert
64.746 Dispenser for MediCups



Article No. 64.746

Description Dispenser for MediCups

- Details**
- for 50 MediCups, 25 ml
 - to be fixed on the wall

Accessories
71.x MediCups



Article No. 64.747

Description Dispenser for MediCup covers

- Details**
- Pic-Box Container, transparent, with Inox dispenser holder
 - to be fixed on the wall

Accessories
72 Cover for MediCups



Accessories for injections

Article No. 64.679-2.0

Description Risk-Waste-Container large

- Details**
- E-SAFE, Universal-A-cover, yellow, 2.0 litre
 - with cannula-, scalpel and Insulin pen removal
 - flab against over-filling
 - closure for disposal
 - extra thick and secure material

Accessories

64.678 Holder for Risk-Waste-Container large for FlexLine Trolleys

Movix2.0 Movix, Vacuum-adhaesion-plate

E-TRAY45101 Syringe-Tray with adapter for small and large Risk-Waste-Containers

- Note**
- 0.5 litre container (No. 64.679-0.5) available
 - Vacuum-adhaesion-plate for 0.5 litre container (No. Movix0.5) available
 - Syringe-Tray with adapter for small Risk-Waste-Containers (E-TRAY45041)



Article No. E-Tray45101

Description Syringe-Tray

- Details**
- white, with Adapter for 0.5litre and for 2.0 litre Risk-Waste-Containers
 - with steril-unsteril-division
 - 350x250x40 mm




Baskets for dish washer

Article No. 910.71.1

Description Dish washer basket with MediCup insert


- Details**
- dish washer basket with lid
 - insert for 22 MediCups
 - 415x90x70 mm



Article No. 910.66.1

Description Dish washer basket with Inlay-Dish insert

- Details**
- dish washer basket with lid
 - insert for 18 Inlay-Dishes
 - 415x90x70 mm





Transparent Sliding-Lids

Sliding-Lid production is a well guarded secret. The covers close the MediDispenser so firmly that no humidity can penetrate. A series of scientific tests in England showed outstanding results.

Material

PET: Polyethylene terephthalate is a polyester plastic (suitable for food). Contains no latex and no softeners. PET is UV-B, UV-C-FUV, UV-C-VUF, EUV and XUV resistant.

Recycling

- Emission free recycling
- Energy recycling (fuel) no reservations
- Very good ecological balance

Cleaning instructions

Washing	Clean with luke-warm water and mild detergent. (Hydrolysis)
Chemicals	Not allowed
Dishwasher	Not suitable
Sterilisation	Gas, radio-isotope (no steam!)

From approx. 45° C, WIEGAND® Sliding-Lids begin to open up their fold, due to the special manufacturing procedure. Thus the strength of the clamp reduces and a hermetic closure is no longer ensured.

Blue MediDispenser and Trays

All blue WIEGAND® products are manufactured deliberately in elaborate deep-drawing procedures. Material formed in this way becomes very resistant, and even withstands dropping. The "blue line" is light sensitive. Avoid UV-exposure.

Material

Polystyrene (food foil). Contains no latex, fluorine, cadmium, lead, bromine or chlorine.

Recycling

- Emission free recycling
- Energy recycling (fuel) no reservations
- Very good ecological balance

Cleaning reference

Chemicals	with frequent use of caustic solutions, alcohol and acids, the plastic becomes brittle
Dishwasher	at 55° C (eco-program) without any problem, at 65° C maximum 45 min. running time at 95° C (industrial rinsing machines) 2 min. running time
Sterilisation	gas, radioisotopes, (no steam!)

WIEGAND® Dispensers have been certified by various European states.



Welcome to our Wiegand Online Shop

Register and profit from our Online Shop:

- ▶ The clear representation of our MVS products helps you when searching for products.
- ▶ You are provided with a quick access to your personal volume-based prices.
- ▶ Move interesting products in to your shopping list, so you can retrieve these products at any time from there and move them to your shopping cart.
- ▶ Read further information on products that fit the product of your choice.
- ▶ Clearly structured article lists and details help you with your order around the clock.

You can sign up below:

You can find the link to the Online Shop on our Website www.wiegand.ch.

Wiegand AG
Schlosserstrasse 5
CH-8180 Bülach
Tel.: +41 44 872 74 74
Fax: +41 44 872 74 75
wiegand@wiegand.ch
www.wiegand.ch

WIEGAND®
Medication Safety