


# UNLOADER®

## Clinically proven to relieve OA knee pain

Unloader braces unload the affected, painful side of the knee using a 3-Point Leverage System. The thigh and calf shells account for two points of leverage, while the Dynamic Force Strap (the diagonal strap above the knee) provides the third. This system “unloads” the pressure from the affected area, providing a reduction in pain.



Knee OA without bracing  
(bone-on-bone contact)



The 3-Point Leverage  
System



Knee OA with bracing  
(space created between bones)



## UNLOADER ONE®

Clinically proven to reduce the pain associated with unicompartmental osteoarthritis of the knee, this brace is designed for average to high-impact activities, indoors as well as outdoors.



## UNLOADER ONE® LITE

At just 10 oz, this is the lightest functional brace available for OA knee pain. Perfect for lean, athletic legs, Unloader One® Lite helps to alleviate pain and improve function.

To learn more, visit [ossur.com/unloader-braces](https://ossur.com/unloader-braces)

FOLLOW ÖSSUR ON



**ÖSSUR**®  
LIFE WITHOUT LIMITATIONS