

BRACING & SUPPORT SOLUTIONS CATALOG



BRACING & SUPPORT SOLUTIONS CATALOG

TABLE OF CONTENTS

ÖSSUR OA SOLUTIONS.....	4	ELASTIC ABDOMINAL BINDER.....	50
UNLOADER® HIP	6	BACK SUPPORT SUSPENDERS.....	50
UNLOADER® HIP.....	7	UPPER EXTREMITY	51
UNLOADER ONE® LITE	8	INNOVATOR X®	53
UNLOADER ONE® LITE	9	FORMFIT® SHOULDER BRACE	54
UNLOADER ONE®	10	SHOULDER ABDUCTION SLING	55
UNLOADER ONE®	11	PREMIUM SHOULDER IMMOBILIZER.....	55
UNLOADER ONE® CUSTOM	13	PADDED STRAP SHOULDER IMMOBILIZER.....	56
UNLOADER ONE® PLUS	16	PREMIUM CONTACT CLOSURE ARM SLING	56
UNLOADER® CUSTOM	17	BUCKLE CLOSURE ARM SLING.....	57
CTI® OA.....	19	BERGSCHULTZ UNIVERSAL ARM SLING.....	57
FORMFIT® PRO KNEE OA	21	PREMIUM PADDED ARM SLING	58
FORMFIT® OA WRAPAROUND.....	22	FRONT CLOSURE CLAVICLE SUPPORT.....	58
ÖSSUR INJURY SOLUTIONS	24	FIGURE-8 CLAVICLE SPLINT	59
INTRODUCTION TO INJURY SOLUTIONS.....	25	SLING AND SWATHE	59
MIAMI LSO™	26	PEDIATRIC UNIVERSAL ARM SLING	60
MIAMI TLSO 464	28	TENNIS ELBOW SUPPORT WITH	
MIAMI TLSO 456	30	HOT/COLD GEL THERAPY	60
MIAMI LUMBAR® POSTEO	32	FORMFIT® TENNIS ELBOW	61
ÖSSUR FORMFIT® PRO BACK	33	AIRFORM® TENNIS ELBOW SUPPORT	61
MIAMI J® SELECT	34	EXOFORM® WRIST	63
MIAMI J®	35	EXOFORM® CARPAL TUNNEL WRIST	63
MIAMI JTO® THORACIC EXTENSION	37	FORMFIT® WRIST / FORMFIT® WRIST & FOREARM.....	64
OCCIAN® BACK & OCCIAN® BACK JR	38	6" WRIST	65
MIAMI Jr®	39	FORMFIT® WRIST & THUMB, UNIVERSAL.....	66
PAPOOSE®	40	FORMFIT® THUMB SPICA.....	67
PHILADELPHIA® ADJUSTABLE		6" THUMB.....	68
TRACHEOTOMY COLLAR	41	UNIVERSAL THUMB SPLINT.....	69
PHILADELPHIA® TRACHEOTOMY COLLAR.....	42	HIP.....	70
UNIVERSAL CERVICAL COLLAR	43	REBOUND® HIP	71
FOAM CERVICAL COLLAR.....	43	KNEE.....	72
NECLOC®	44	REBOUND® CARTILAGE	73
NECLOC® KIDS	44	REBOUND® PCL.....	76
PATRIOT® COLLAR.....	45	CTI® CUSTOM	78
OAM RIGID LUMBAR	46	CTI®	81
FORMFIT® BACK SUPPORT, AIR	47	REBOUND® DUAL	84
FORMFIT® BACK SUPPORT	48	REBOUND® DUAL ST	86
UNIVERSAL ELASTIC RIB BELT	49	REBOUND® DUAL ST CUSTOM.....	88
UNIVERSAL ABDOMINAL BINDER.....	49	PARADIGM®	90

BRACING & SUPPORT SOLUTIONS CATALOG

EXTREME®	91	AIRFORM® INFLATABLE ANKLE STIRRUP	130
FLEX™	92	AIRFORM® PRE-INFLATED ANKLE STIRRUP	130
FORMFIT® MCL	94	FORMFIT® HONEYCOMB ANKLE STIRRUP	131
C180™ ROCKET	95	FORMFIT® FOAM ANKLE BRACE.....	131
REBOUND® POST OP KNEE.....	96	GEL ANKLE STIRRUP	132
FORMFIT® POST-OP KNEE.....	97	EXOFORM® DORSAL NIGHT SPLINT	132
EXOFORM® KNEE IMMOBILIZER.....	98	FORMFIT® NIGHT SPLINT	133
FORMFIT® KNEE IMMOBILIZER.....	99	AIRFORM® NIGHT SPLINT	134
FX™ PATELLA STABILIZER	100	SQUARE TOE POST-OP SHOE	135
FORMFIT® HINGED LATERAL J.....	101	SOFT TOP POST-OP SHOE.....	135
FORMFIT® KNEE, HINGED.....	102	MESH TOP POST-OP SHOE	136
FORMFIT® KNEE ROM	103	CANVAS ROCKER BOTTOM CAST SHOE	136
REBOUND® KNEE BRACE	104	SPECIALITY CARE,	
NEOPRENE HINGED KNEE SUPPORT	105	CASTING & TRACTION	137
NEOPRENE WRAPAROUND		COLD RUSH®	138
HINGED KNEE SUPPORT	105	COLD RUSH® COMPACT	138
NEOPRENE KNEE SUPPORT		COLD RUSH® PADS.....	139
WITH STABILIZED PATELLA	106	ELASTIC BANDAGE®	140
NEOPRENE KNEE SLEEVE	108	CAST PADDING.....	140
PREMIUM ELASTIC KNEE SUPPORT	109	STOCKINETTE	140
ELASTIC KNEE SUPPORT	109	ALOETOUCH® EXAM GLOVES.....	141
FOOT AND ANKLE.....	110	TECHFORM® PREMIUM	141
REBOUND® DIABETIC WALKER.....	111	TECHFORM® ORTHO-ROLL®	142
REBOUND® AIR WALKER	113	TECHFORM® PRE-CUT SPLINTS.....	142
AFO DYNAMIC	115		
AFO LEAF SPRING	116		
AFO LIGHT.....	117		
REBOUND® FOOT-UP®	118		
FORMFIT® WALKER, AIR	119		
FORMFIT® WALKER	120		
ROM WALKER	121		
DH OFFLOADING WALKER®	122		
DH OFFLOADING POST-OP STYLE SHOE	122		
PEDIATRIC WALKER	123		
METAL UPRIGHT AIR WALKER	124		
REBOUND® ANKLE BRACE.....	126		
FORMFIT® ANKLE BRACE.....	127		
FORMFIT® ANKLE WITH SPEEDLACE.....	128		
AIRFORM® UNIVERSAL INFLATABLE			
ANKLE STIRRUP.....	129		

ÖSSUR OA SOLUTIONS


Enhance quality of life and mobility for people living with osteoarthritis.

These product solutions have been built on decades of experience and are designed to enhance the quality of life of people living with osteoarthritis by reducing pain and improving mobility. Through innovation leadership, OA expertise and partnerships with medical professionals, we offer comprehensive, proven and effective treatment options.

UNLOADER®

The Össur Unloader® brand name has become synonymous with Osteoarthritis bracing, representing OA Solutions' technology as the industry standard, with more than 15 clinical studies demonstrating its efficacy.

Össur strives to maintain its position as the innovation leader for OA Solutions by continuing to develop simple, yet effective clinically proven treatment solutions that increase people's quality of life.





HIP OA



UNLOADER® HIP

PAGE 6

KNEE OA



UNLOADER ONE® LITE

PAGE 8



UNLOADER ONE®

PAGE 10

UNLOADER[®] HIP

Unloader Hip is designed to reduce pain by optimizing load dispersion for patients suffering from mild and moderate osteoarthritis of the hip, thus contributing to hip stability and improved mobility.




INDICATIONS FOR USE

- Mild and moderate osteoarthritis of the hip
- Conditions benefitting from repositioning of the hip and/or compression around the greater trochanters

FEATURES AND BENEFITS

- The Dynamic Rotation Straps with SmartDosing dial allows individual adjustment of the external rotation force, to achieve optimum pain relief for each patient
- Easy donning and doffing thanks to pocket wrap closure system and zipper
- Pulley system for improved stability
- Premium quality lycra for breathability and a snug fit



REQUIRED MEASUREMENTS

- Measure the hip circumference at the anterior superior iliac spine.



UNLOADER[®] HIP

ORDERING INFORMATION

UNLOADER HIP

Size	Hip Circumference	Woman	Men
Small	31.5 – 35" (80 – 88 cm)	B-465350030	B-465350031
Medium	35 – 38" (88 – 96 cm)	B-465350032	B-465350033
Large	38 – 41" (96 – 104 cm)	B-465350034	B-465350035
X-Large	41 – 45" (104 – 114 cm)	B-465350036	B-465350037
XX-Large	45 – 49" (114 – 124 cm)	B-465350038	B-465350039

ACCESSORIES/ OPTIONS

REPLACEMENT TROUSERS

Size	Measurements	Woman	Men
Small	31.5 – 35" (80 – 88 cm)	B-465350077	B-465350078
Medium	35 – 38" (88 – 96 cm)	B-465350079	B-465350080
Large	38 – 41" (96 – 104 cm)	B-465350081	B-465350082
X-Large	41 – 45" (104 – 114 cm)	B-465350083	B-465350084
XX-Large	45 – 49" (114 – 124 cm)	B-465350085	B-465350086

OTHER ACCESSORIES

Description	Part#
Rotation strap	B-465350064
Washbag	B-465350065
Trochanter pad	B-465350066
Thigh strap S	B-465350067
Thigh strap M	B-465350068
Thigh strap L	B-465350069
Thigh strap XL	B-465350070
Thigh strap XXL	B-465350071
Bilateral kit* S	B-465350072
Bilateral kit* M	B-465350073
Bilateral kit* L	B-465350074
Bilateral kit* XL	B-465350075
Bilateral kit* XXL	B-465350076

*The bilateral kit consists of a rotation control strap, a thigh strap and a trochanter pad

UNLOADER ONE® LITE

Unloader One Lite is a slim and low profile knee brace that relieves pain from knee osteoarthritis and degenerative meniscal tears by means of the clinically proven Dynamic Force System that provides uni-compartmental unloading of the knee joint.

With an integrated Lycra® sleeve and shells with overmold, the brace is an optimal combination of comfort and sleek look.



INDICATIONS FOR USE


- Mild to moderate unicompartmental osteoarthritis
- Degenerative meniscus
- Other knee conditions that may benefit from unicompartmental unloading

FEATURES AND BENEFITS

- 3-Point Leverage System provides clinically relevant pain relief of affected compartment
- Single SmartDosing dial allows fine-tuned pain management
- Helps keep patients mobile during daily activities
- Discreet hinge and integrated sleeve ensure a sleek look on leg and snug fit underneath clothing
- Lightweight and easy to use

REQUIRED MEASUREMENTS

- Circumference measurement: 6"/15 cm above and below mid-patella



UNLOADER ONE® LITE

ORDERING INFORMATION

UNLOADER ONE LITE

Size	Circumference 6" (15 cm) above mid patella	Circumference 6" (15 cm) below mid-patella	Left medial	Right medial	Left lateral	Right lateral
X-Small	12 - 15" (30 - 38 cm)	10 - 12" (25 - 30 cm)	UOL02300	UOL01300	UOL02400	UOL01400
Small	15 - 17.75" (38 - 45 cm)	12 - 13.75" (30 - 35 cm)	UOL02301	UOL01301	UOL02401	UOL01401
Medium	17.75 - 22" (45 - 56 cm)	13.75 - 16.5" (35 - 42 cm)	UOL02302	UOL01302	UOL02402	UOL01402
Large	22 - 26" (56 - 66 cm)	16.5 - 20" (42 - 51 cm)	UOL02303	UOL01303	UOL02403	UOL01403
X-Large	26 - 29" (66 - 73 cm)	20 - 24" (51 - 61 cm)	UOL02304	UOL01304	UOL02404	UOL01404

Hinge is placed on the side of the affected compartment

ACCESSORIES

Size	Part#
UNLOADER ONE LITE DFS KIT XS	UNF01302
UNLOADER ONE LITE DFS KIT SM	UNF01304
UNLOADER ONE LITE DFS KIT LG	UNF01308
UNLOADER ONE LITE DFS KIT XL	UNF01310
UNLOADER ONE LITE ALLIGATOR DFS KIT	UNF01502
UNLOADER ONE LITE CALF STRAP KIT XS	UNF01402
UNLOADER ONE LITE CALF STRAP KIT SM	UNF01404
UNLOADER ONE LITE CALF STRAP KIT MD	UNF01406
UNLOADER ONE LITE CALF STRAP KIT LG	UNF01408
UNLOADER ONE LITE CALF STRAP KIT XL	UNF01410
STOP UNLR ONE FLEX/EXT PK GRAY	B-301040046
UNLOADER ONE LITE HINGE PAD (Qty 10)	UNF00950

UNLOADER ONE®

The Unloader One is a comfortable, lightweight, low-profile brace that provides excellent suspension and pain relief. This lightweight brace preserves the clinically proven Össur Unloader 3-Points of Leverage design by incorporating Dual Dynamic Force Straps that are easily fine-tuned with the SmartDosing® dial. Color-coded Quick Fit Buckles simplify application and removal, while flexible shells and breathable tibial liners, coated with Össur Sensil® Silicone, improve comfort and reduce migration, helping to ensure patient compliance. Recent research demonstrates that the Unloader One reduces pain, improves function and decreases the use of pain medication.



INDICATIONS FOR USE

- Moderate to severe unicompartmental osteoarthritis

FEATURES AND BENEFITS

- Dual Dynamic Force Straps designed to reduce pressure on sensitive bony surfaces of the knee, resulting in reduction of pain
- Tool-free SmartDosing dial enables fine-tuned pain management and can be adjusted even when wearing the brace under clothing
- Color coded quick fit buckles enable easy application and removal
- Proprietary Sensil® silicone liners designed to reduce migration
- Flexible, low-profile upright contours to patient anatomy and allows use underneath clothing
- Quick Connect Strapping for tool-free removal and replacement
- Lightweight frame is perforated for flexibility and ventilation with custom paint options for the Custom version
- Unloader One is available in pre-fabricated and Custom version

REQUIRED MEASUREMENTS

- Circumference measurement: 6" (15 cm) below mid-patella

CUSTOM STANDARD COLORS



CUSTOM PAINTED DESIGNS



No charge for custom standard colors. Additional fees required for other custom colors, custom painted designs and custom printed graphics.



UNLOADER ONE®

ORDERING INFORMATION

UNLOADER ONE WITH SMARTDOSING

Size	Calf circumference	Left medial	Right medial	Left lateral	Right lateral
Small	12.2 – 14" (31 – 36 cm)	B-240618712	B-240518712	B-240628712	B-240528712
Medium	14 – 16.75" (36 – 43 cm)	B-240618713	B-240518713	B-240628713	B-240528713
Large	16.75 – 20" (43 – 51 cm)	B-240618714	B-240518714	B-240628714	B-240528714
X-Large	20 – 24" (51 – 61 cm)	B-240618715	B-240518715	B-240628715	B-240528715

UNLOADER ONE SHORT WITH SMARTDOSING

Size	Calf circumference	Left medial	Right medial	Left lateral	Right lateral
X-Small	9.8 – 12.2" (25 – 31 cm)	B-241618711	B-241518711	B-241628711	B-241528711
Small	12.2 – 14" (31 – 36 cm)	B-241618712	B-241518712	B-241628712	B-241528712
Medium	14 – 16.75" (36 – 43 cm)	B-241618713	B-241518713	B-241628713	B-241528713
Large	16.75 – 20" (43 – 51 cm)	B-241618714	B-241518714	B-241628714	B-241528714
X-Large	20 – 24" (51 – 61 cm)	B-241618715	B-241518715	B-241628715	B-241528715

UNLOADER ONE WITH RATCHET

Size	Calf circumference	Left medial	Right medial	Left lateral	Right lateral
Small	12.2 – 14" (31 – 36 cm)	B-240619712	B-240519712	B-240629712	B-240529712
Medium	14 – 16.75" (36 – 43 cm)	B-240619713	B-240519713	B-240629713	B-240529713
Large	16.75 – 20" (43 – 51 cm)	B-240619714	B-240519714	B-240629714	B-240529714
X-Large	20 – 24" (51 – 61 cm)	B-240619715	B-240519715	B-240629715	B-240529715

UNLOADER ONE SHORT WITH RATCHET

Size	Calf circumference	Left medial	Right medial	Left lateral	Right lateral
X-Small	9.8 – 12.2" (25 – 31 cm)	B-241619711	B-241519711	B-241629711	B-241529711
Small	12.2 – 14" (31 – 36 cm)	B-241619712	B-241519712	B-241629712	B-241529712
Medium	14 – 16.75" (36 – 43 cm)	B-241619713	B-241519713	B-241629713	B-241529713
Large	16.75 – 20" (43 – 51 cm)	B-241619714	B-241519714	B-241629714	B-241529714
X-Large	20 – 24" (51 – 61 cm)	B-241619715	B-241519715	B-241629715	B-241529715

Standard version is recommended for patients between 5'5" and 6'4" (1,65 and 1.92 meters tall).

Short version is recommended for patients under 5'5" (1,65 meters tall).

Hinge is placed on the side of the affected compartment

ACCESSORIES

Accessory Kits	SMARTDOSING GREY		RATCHET GREY	
	RM/LL	RL/LM	RM/LL	RL/LM
Strap & Strap Pad & Buckle Clasp Kit	B-705277003	B-705278003	B-701217008	B-701218008
Strap & Strap Pad & Buckle Clasp Kit, Short	B-705277004	B-705277005	-	-
Strap & Strap Pad Kit	B-705270003		B-701210006	
Strap & Strap Pad Kit, Short	B-705277006		B-701240004	
DFS Strap Kit, Standard	B-701240002		-	-
DFS Strap Kit, Short	B-701240012		-	-
Pads Kit	-	-	B-702210003	

UNLOADER ONE®

Description	Right medial/Left lateral	Right lateral/Left medial
Buckle Clasp Kit	B-708217003	B-708218003
Alligator Tab (Qty. 4)	B-704210003	
Alligator Tab	B-310005334	
Extension Stop Kit	B-301040046	
Velcro Twist-Fit Attachment Kit	B-721240003	
Hinge Screws	B-301010076	
Post	B-301010075	

ÖSSUR KNEE UNDERSLEEVE

Size	Thigh Circumference	Part#
X-Small	12 – 15" (30.5 – 38 cm)	B-250150201
Small	15 – 18" (38 – 45.5 cm)	B-250150202
Medium	18 – 22" (45.5 – 56 cm)	B-250150203
Large	22 – 26" (55 – 66 cm)	B-250150204
X-Large	26 – 29" (66 – 73.5 cm)	B-250150205
XX-Large	29 – 33" (73.5 – 84 cm)	B-250150206

SUSPENSION STRAP KIT

Description	Part#
Suspension Strap Kit	UNL00212

PADS

Description	Part#
Hinge Pads	B-304120301
Popliteal Pad/ strap Grey 6in	B-304010103

When using an Össur Knee Undersleeve we advise using the suspension strap kit to minimize migration

SENSIL LINERS

Size	Calf		Thigh	
	Right medial/Left lateral	Right lateral/Left medial	Right medial/Left lateral	Right lateral/Left medial
X-Small	B-419700311	B-419800311	B-409700321	B-409800321
Small	B-419700312	B-419800312	B-409700322	B-409800322
Medium	B-419700313	B-419800313	B-409700323	B-409800323
Large	B-419700314	B-419800314	B-409700324	B-409800324
X-Large	B-419700315	B-419800315	B-409700325	B-409800325

DOESKIN LINERS

Size	Calf		Thigh	
	Right medial/Left lateral	Right lateral/Left medial	Right medial/Left lateral	Right lateral/Left medial
X-Small	B-424703721	B-424803721	B-412701721	B-412801721
Small	B-424703722	B-424803722	B-412701722	B-412801722
Medium	B-424703723	B-424803723	B-412701723	B-412801723
Large	B-424703724	B-424803724	B-412701724	B-412801724
X-Large	B-424703725	B-424803725	B-412701725	B-412801725

When using a doeskin calf liner we advise using a suspension strap kit to minimize migration

SENSIL AMS (ANTI-MIGRATION SYSTEM) WRAPS

Size	Right medial/Left lateral	Right lateral/Left medial
X-Small	B-715217031	B-715218031
Small	B-715217032	B-715218032
Medium	B-715217033	B-715218033
Large	B-715217034	B-715218034
X-Large	B-715217035	B-715218035

UNLOADER ONE® CUSTOM



SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.

Description	Part#
SmartMeasure Measuring Kit	B-312000005
Constant Force Caliper	B-312000000

Scans, CMS (CTI only), CCS and casting also accepted. Knee M/L required

ORDERING INFORMATION

UNLOADER ONE CUSTOM WITH SMARTDOSING

Color	Left	Right
Black	B-115600112	B-115500112
Grey	B-115601712	B-115501712

UNLOADER ONE CUSTOM

Color	Left	Right
Black	B-115600111	B-115500111
Grey	B-115601711	B-115501711

ACCESSORIES

ÖSSUR KNEE UNDERSLEEVE

Size	Circumference	Grey	Black
X-Small	12 – 15" (30.5 – 38 cm)	B-250150201	B-250150101
Small	15 – 18" (38 – 45.5 cm)	B-250150202	B-250150102
Medium	18 – 22" (45.5 – 56 cm)	B-250150203	B-250150103
Large	22 – 26" (56 – 66 cm)	B-250150204	B-250150104
X-Large	26 – 29" (66 – 73.5 cm)	B-250150205	B-250150105
XX-Large	29 – 33" (73.5 – 84 cm)	B-250150206	B-250150106

When using an Össur Knee Undersleeve we advise using the suspension strap kit to minimize migration

SENSIL LINERS

Size	Color	Calf		Thigh	
		Right medial/Left lateral	Right lateral/Left medial	Right medial/Left lateral	Right lateral/Left medial
X-Small	Grey	B-419700311	B-419800311	B-409700321	B-409800321
Small	Grey	B-419700312	B-419800312	B-409700322	B-409800322
Medium	Grey	B-419700313	B-419800313	B-409700323	B-409800323
Large	Grey	B-419700314	B-419800314	B-409700324	B-409800324
X-Large	Grey	B-419700315	B-419800315	B-409700325	B-409800325
X-Small	Black	B-419700111	B-419800111	B-409700121	B-409800121
Small	Black	B-419700112	B-419800112	B-409700122	B-409800122
Medium	Black	B-419700113	B-419800113	B-409700123	B-409800123
Large	Black	B-419700114	B-419800114	B-409700124	B-409800124
X-Large	Black	B-419700115	B-419800115	B-409700125	B-409800125

UNLOADER ONE® CUSTOM

ACCESSORIES

DOESKIN LINERS

Size	Color	Calf		Thigh	
		Right medial/Left lateral	Right lateral/Left medial	Right medial/Left lateral	Right lateral/Left medial
X-Small	Grey	B-424703721	B-424803721	B-412701721	B-412801721
Small	Grey	B-424703722	B-424803722	B-412701722	B-412801722
Medium	Grey	B-424703723	B-424803723	B-412701723	B-412801723
Large	Grey	B-424703724	B-424803724	B-412701724	B-412801724
X-Large	Grey	B-424703725	B-424803725	B-412701725	B-412801725
X-Small	Black	B-424702121	B-424802121	B-412700121	B-412800121
Small	Black	B-424702122	B-424802122	B-412700122	B-412800122
Medium	Black	B-424702123	B-424802123	B-412700123	B-412800123
Large	Black	B-424702124	B-424802124	B-412700124	B-412800124
X-Large	Black	B-424702125	B-424802125	B-412700125	B-412800125

When using a doeskin calf liner we advise using a suspension strap kit to minimize migration. For replacement liners, please contact Customer Service.

SENSIL AMS (ANTI-MIGRATION SYSTEM) WRAPS

Size	Grey		Black	
	Right medial/Left lateral	Right lateral/Left medial	Right medial/Left lateral	Right lateral/Left medial
X-Small	B-715217031	B-715218031	B-715217011	B-715218011
Small	B-715217032	B-715218032	B-715217012	B-715218012
Medium	B-715217033	B-715218033	B-715217013	B-715218013
Large	B-715217034	B-715218034	B-715217014	B-715218014
X-Large	B-715217035	B-715218035	B-715217015	B-715218015

SUSPENSION STRAP KIT

Description	Grey	Black
Suspension Strap Kit	UNL00212	UNL00211

SMARTDOSING ACCESSORY KITS

Description	Grey RM/LL	Grey RL/LM	Black RM/LL	Black RL/LM
SmartDosing Strap & Strap Pad & Buckle Clasp Kit	B-705277003	B-705278003	B-705277001	B-705278001
Strap & Strap Pad Kit	B-705270003		B-705270001	
Strap Kit	B-701270013		B-701270011	

UNLOADER ONE® CUSTOM

ACCESSORIES

RATCHET ACCESSORY KITS

Description	Grey RM/LL	Grey RL/LM	Black RM/LL	Black RL/LM
Strap & Strap Pad & Buckle Clasp Kit	B-701217008	B-701218008	B-701217009	B-701218009
Strap & Strap Pad Kit	B-701210006		B-701210007	
Strap Kit	B-701210004		B-701210005	
Pad Kit	B-702210003		B-702210001	
Ratchet Complete Strap Set	B-701210004			

ACCESSORY KITS

Description	Grey RM/LL	Grey RL/LM	Black RM/LL	Black RL/LM
Buckle Clasp Kit	B-708217003	B-708218003	B-708217001	B-708218001
Alligator Clip Kit (included in all strap kits)	B-704210003		B-704210001	
Alligator Tab	B-310005334		B-31005335	
Extension Stop Kit	B-301040046		B-301040047	
Velcro Twist-Fit Attachment Kit	B-721240003		B-721240001	

UNLOADER ONE® PLUS

ORDERING INFORMATION

UNLOADER ONE PLUS WITH RATCHET

Size	Thigh Circumference	Calf circumference	Left medial	Right medial	Left lateral	Right lateral
Large	21 – 25" (53 – 64 cm)	14 – 16.75" (36 – 43 cm)	B-245610004	B-245510004	B-245620004	B-245520004
X-Large	25 – 29" (64 – 74 cm)	16.75 – 20" (43 – 51 cm)	B-245610005	B-245510005	B-245620005	B-245520005

UNLOADER ONE PLUS WITH RATCHET, LONG VERSION

Size	Thigh Circumference	Calf circumference	Left medial	Right medial	Left lateral	Right lateral
Large	21 – 25" (53 – 64 cm)	14 – 16.75" (36 – 43 cm)	B-244610004	B-244510004	B-244620004	B-244520004
X-Large	25 – 29" (64 – 74 cm)	16.75 – 20" (43 – 51 cm)	B-244610005	B-244510005	B-244620005	B-244520005

ACCESSORIES

SENSIL AMS (ANTI-MIGRATION SYSTEM) WRAPS

Size	Right medial/Left lateral	Right lateral/Left medial
Large	B-715237034	B-715238034
X-large	B-715237035	B-715238035

LINERS

Size	Calf, sensil		Thigh, doeskin	
	Right medial/Left lateral	Right lateral/Left medial	Right medial/Left lateral	Right lateral/Left medial
Large	B-413700314	B-413800314	B-414700324	B-414800324
X-Large	B-413700315	B-413800315	B-414700325	B-414800325

ACCESSORY KITS

Description	Right medial/Left lateral	Right lateral/Left medial
Strap & Strap Pad & Buckle Clasp Kit	B-701237008	B-701238008
Strap & Strap Pad Kit	B-701230006	
Strap Kit	B-701230004	
Pad Kit	B-702230003	
Extra Pads	B-702230013	
Buckle Clasp Kit	B-708237003	B-708238003
Alligator Clip Kit (included in all strap kits)	B-704210003	
Gastroc Kit	B-705230001	
Elastic Thigh Strap 8"	B-701237083	B-701238083
Elastic Thigh Strap 10"	B-701237103	B-701238103

UNLOADER[®] CUSTOM

With a semi-rigid frame and an adjustable hinge, the Unloader Custom accommodates a wide variety of leg shapes and sizes. The Adjustable Dynamic Joint enables changes to the hinge and additional leverage for optimal pain management.



INDICATIONS FOR USE

- Mild to severe unicompartmental osteoarthritis

FEATURES AND BENEFITS

- Custom molded semi-flexible shells are durable, adaptable to leg shape changes, yet rigid enough to contain tissue while force is applied at the knee
- Patented Adjustable Dynamic Joint™ (ADJ) high-carbon steel hinge is strong enough to unload the pressure from the compartment of the knee. Adjustments are made easily
- The Dynamic Force Strap® designed to reduce pressure on sensitive bony surfaces of the knee, resulting in reduction of pain
- Available in standard and in SmartDosing version. Tool-free SmartDosing dial option enables fine-tuned pain management and can be adjusted even when wearing the brace under clothing
- SmartDosing version available with Sensil® silicone calf liner to reduce migration
- Available in 12 standard colors in matte or glossy finish or custom painted


REQUIRED MEASUREMENTS

- Custom-to-measurement

CUSTOM STANDARD COLORS



CUSTOM PAINTED DESIGNS



No charge for custom standard colors. Additional fees required for other custom colors, custom painted designs and custom printed graphics.



UNLOADER[®] CUSTOM



SMARTMEASURE[®]

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.

Description	Part#
SmartMeasure Measuring Kit	B-312000005
Constant Force Caliper	B-312000000

Scans, CMS (CTI only), CCS and casting also accepted. Knee M/L required

ORDERING INFORMATION

UNLOADER CUSTOM

Description	Part#
Right	B-122500020
Left	B-122600020

Description	Part#
Proximal Non-Elastic Thigh Strap	B-507010000
AMS Wrap (custom made)	B-426000150
AMS Sleeve (custom made)	B-426000250
ACL Band Kit	B-704000003

ACCESSORIES

Description	Part#
Knee Shield	B-250110000
Hinge cover	B-251030000
Water Ski Strap Lock	B-303010084
Distal Non-Elastic Thigh Strap	B-507010000

ELASTIC CALF SUSPENSION STRAP

Size	Part#	Size	Part#	Size	Part#
14"	B-504020114	18"	B-504020018	22"	B-504020022
15"	B-504020115	19"	B-504020019	23"	B-504020023
16"	B-504020116	20"	B-504020020	24"	B-504020024
17"	B-504020017	21"	B-504020021		

SLEEVES

Size	Circumference	Sofsleeve	Neoprene undersleeve	Neoprene oversleeve
Small	< 15" (38 cm)	20643	20251	20051
Medium	15 – 18" (38 – 46 cm)	20644	20252	20052
Large	18 – 19.5" (46 – 50 cm)	20644	20253	20053
X-Large	19.5 – 21" (50 – 55 cm)	20645	20254	20054
XX-Large	21.5 – 24.5" (55 – 62 cm)	20645	20255	20055
XXX-Large	24.5 – 28" (62 – 71 cm)	20646	20256	20056
XXXX-Large	28 – 34" (71 – 86 cm)	20647	20257	20057

Size	Measuring	Össur knee undersleeve
X-Small	12 – 15" (30.5 – 38 cm)	B-250150101
Small	15 – 18" (38 – 45.5 cm)	B-250150102
Medium	18 – 22" (45.5 – 56 cm)	B-250150103

Size	Measuring	Össur knee undersleeve
Large	22 – 26" (56 – 66 cm)	B-250150104
X-Large	26 – 29" (66 – 73.5 cm)	B-250150105
XX-Large	29 – 33" (73.5 – 84 cm)	B-250150106

CTI® OA

CTI® OA's composite shells can be constructed with up to 7° of varus or valgus correction, offering a high level of protection and stability for clinicians interested in treating osteoarthritis and ligament instability equally. Key features include a Total Support System, Accutrac® hinges that mimic the natural motion of the knee joint and breathable liners, coated with Össur® Sensil® Silicone to improve efficacy and virtually eliminate migration.



INDICATIONS FOR USE

- Mild to severe unicompartmental osteoarthritis
- ACL, MCL, LCL, PCL, rotary and combined instabilities
- Demanding activities of daily living

FEATURES AND BENEFITS

- Anatomically-correct Accutrac® hinges with extension stops
- Proprietary Sensil® silicone liners designed to reduce migration
- Total Support System™ for ACL, MCL, LCL, PCL, rotary and combined instabilities
- Carbon composite, rigid frame construction
- Custom-fabricated to each patient's unique measurements with up to 7° varus/valgus correction built into the frame
- Non-corrosive materials, ideal for water sports
- Standard model for medium to high impact levels. Vapor model for low to medium impact levels. Pro Sport model for highest impact levels

REQUIRED MEASUREMENTS

- Custom-made to measure



CTI® OA

ORDERING INFORMATION

CTI OA

Description	Left	Right
Standard	B-116630010	B-116530010
Pro Sport	B-117630010	B-117530010
Vapor	B-118630010	B-118530010



SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.

Description	Part#
SmartMeasure Measuring Kit	B-312000005
Constant Force Caliper	B-312000000

Scans, CMS (CTI only), CCS and casting also accepted. Knee M/L required

CUSTOM STANDARD COLORS



CUSTOM PAINTED DESIGNS



No charge for custom standard colors. Additional fees required for other custom colors, custom painted designs and custom printed graphics.

ACCESSORIES

Description	Part#
OA Condyle Pad Kit 1/4"	25021
OA Condyle Pad Kit 3/8"	25644

For all other Cti OA Accessories, see CTi Custom Accessories. Page 76.

FORMFIT® PRO KNEE OA

Össur Formfit® Pro Knee OA is a unique, compressive knee sleeve with a Dynamic Force Strap system, designed to provide pain relief from the early stages of osteoarthritis (OA). The dynamic, 3-D MotionTech™ knit is supportive and meets certified medical grade compression requirements, and has been proven to be superior to the leading brands at moisture wicking and breathability, keeping your knee cool during wear. CoolVent™ technology in the popliteal area enhances breathability and comfort, while the distinct patella graphics help with intuitive placing of the sleeve. The integrated StayBilizer™ hinge, constructed from lightweight, flexible EVA material, is suitable for sports, such as football.



MotionTech™



CoolVent™



StayBilizer™



FoldLock™

INDICATIONS FOR USE

For osteoarthritis and other knee conditions which may benefit from unloading and/or compression, such as:

- Medial or lateral knee pain
- Early or mild unicompartmental osteoarthritis of the knee
- Symptomatic, unicompartmental meniscus degeneration

CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated.

FEATURES AND BENEFITS

- Patented DFS system.
- Certified medical grade compression (20-30 mmHg).
- Lightweight CoolVent™ knit in popliteal area for added breathability and comfort.
- Patented FoldLock™ system for easier and more secure application.
- Sports-compatible StayBilizer™ hinge made from light, flexible EVA material for optimum shape and fit.
- Patella graphics for intuitive placement.
- Available in 5 standard sizes and 3 comfort sizes.

REQUIRED MEASUREMENTS



Measure calf circumference 15 cm (6") below mid-patella. Leg in extension.



Measure thigh circumference 15 cm (6") above mid-patella. Leg in extension.

Size	Calf circumference 15 cm (6") below mid patella		Thigh circumference 15 cm (6") above mid patella		Right Medial, Left Lateral	Right Lateral, Left Medial
	cm	inch	cm	inch		
XS	26 - 31	10¼ - 12¼	36 - 41	14½ - 16½	FP110401	FP110301
S	31 - 36	12¼ - 14½	41 - 46	16½ - 18½	FP110402	FP110302
M	36 - 41	14½ - 16½	46 - 51	18½ - 20½	FP110403	FP110303
M Comfort	36 - 41	14½ - 16½	51 - 56	20½ - 22	FP110413	FP110313
L	41 - 46	16½ - 18½	51 - 56	20½ - 22	FP110404	FP110304
L Comfort	41 - 46	16½ - 18½	56 - 60	22 - 23½	FP110414	FP110314
XL	46 - 50	18½ - 19¾	56 - 60	22 - 23½	FP110405	FP110305
XL Comfort	46 - 50	18½ - 19¾	60 - 64	23½ - 25¼	FP110415	FP110315

FORMFIT® OA WRAPAROUND

The Formfit OA Wraparound is a cost effective, low-profile option for patients with mild to moderate unicompartmental knee osteoarthritis. The Dynamic Force Strap tensions as the leg extends, providing pain relief in a single pull motion. The slim polycentric hinge design and breathable material offers a discreet and comfortable fit. Dot-coded color wrap closures and non-slippery rubber tabs allows for a simple donning procedure for the user.



INDICATIONS FOR USE

- Mild to moderate unicompartmental knee osteoarthritis

FEATURES AND BENEFITS

- Easy wraparound design
- Dot coded closure wraps
- Non-slippery rubber tab
- Dynamic Force Strap
- Low-profile design & slim hinge
- Breathable lycra

REQUIRED MEASUREMENTS

- Thigh Circumference: 15 cm (6") above mid-patella
- Calf Circumference: 15 cm (6") below mid-patella

ORDERING INFORMATION – FORMFIT® OA WRAPAROUND

FORMFIT OA WRAPAROUND

Size	Thigh Circumference	Calf Circumference	Left medial	Right medial	Left lateral	Right lateral
X-Small	14 – 16" (36 – 41 cm)	11 – 12.5" (28 – 32 cm)	B-254901008	B-254901001	B-254901022	B-254901015
Small	16 – 18" (41 – 46 cm)	12.5 – 14" (32 – 36 cm)	B-254901009	B-254901002	B-254901023	B-254901016
Medium	18 – 21" (46 – 53 cm)	14 – 16" (36 – 41 cm)	B-254901010	B-254901003	B-254901024	B-254901017
Large	21 – 23.5" (53 – 60 cm)	16 – 18" (41 – 46 cm)	B-254901011	B-254901004	B-254901025	B-254901018
X-Large	23.5 – 26.5" (60 – 67 cm)	18 – 20" (46 – 51 cm)	B-254901012	B-254901005	B-254901026	B-254901019
2X-Large	26.5 – 29.5" (67 – 75 cm)	20 – 22" (51 – 56 cm)	B-254901013	B-254901006	B-254901027	B-254901020
3X-Large	29.5 – 32" (75 – 81 cm)	22 – 24" (56 – 61 cm)	B-254901014	B-254901007	B-254901028	B-254901021

ACCESSORIES

SUSPENSION STRAP

Description	Part#
Suspension Strap	OAW00007

ÖSSUR INJURY SOLUTIONS

For people recovering from fractures, ligament injuries or in need of post-operative treatment solutions.

These product solutions have been built on decades of experience and are designed to enhance the healing process of bone and soft tissue injuries. Through innovation leadership, proven biomechanical approaches, and partnerships with leading medical professionals, we strive to offer comprehensive and effective treatment options.

Össur Injury Solutions are globally recognized, primarily under the brand names CTi®, Rebound® and Miami®, all of which are developed to assist healing and increase people's quality of life.

Our expertise in Injury Solutions has driven us to take the next step forward by developing indication-specific products which are designed to optimize healing while maintaining function and mobility. These Functional Healing products incorporate Össur proprietary technology and are designed to enhance clinically relevant outcomes.



FUNCTIONAL HEALING –
WHERE MOBILITY MEANS RECOVERY.



SPINE



MIAMI J®
PAGE 35



FUNCTIONAL HEALING



REBOUND® CARTILAGE
PAGE 73



REBOUND® PCL
PAGE 76

KNEES



CTi®
PAGE 81



REBOUND® DUAL
PAGE 84

FEET



REBOUND® AIR WALKER
PAGE 113

MIAMI LSO™

Miami LSO is a spinal orthosis designed for pain relief, protection of injured ligaments or muscles and post-surgical immobilization. It provides low-profile, anatomical immobilization, easy fit and wear, options for universal or modular support and an advanced pulley system with a length-adjustable pulley cord.



INDICATIONS FOR USE

Conditions requiring gross immobilization of the trunk in the lumbar region. Such conditions may include:

- Stable, non-displaced spinal fractures (L1-L5)
- Spinal stenosis
- Herniated discs
- Degenerative spinal pathologies
- Spondylolisthesis
- Spondylolysis

Contra-indication: Unstable, displaced spinal fractures

FEATURES AND BENEFITS

- Slim, breathable material
- Proven Miami double-pulley system with length-adjustable pulley cord
- Intuitive fitting without trimming with fold-over belt arms & Precise-Size™ landmarks
- Optional fit tool for perfect fitting directly on the patient
- Universal belt sizing accommodates waists up to 178 cm/70"
- New Universal Posterior panel to accommodate more patients with less inventory
- Contoured posterior panels sized S-XL for a precise, anatomical fit
- Quick-Fit pocket enclosure and blue highlights for patient ease of use
- Immobilization from L1 – L5 with coverage from S1 – T9


REQUIRED MEASUREMENTS



- Waist circumference






MIAMI LSO™ ORDERING INFORMATION


OPTION 1: Order braces for inventory

MIAMI LSO BELT		
		
SIZE	WAIST CIRC.	PART #
6" Universal (Slim)	26-50" (66-127 cm)	MLSO-1
8" Universal (Standard)	26-50" (66-127 cm)	MLSO-2
6" (Slim) Belt Extension	26-70" (66-177 cm)	MLSO-1EXT
8" (Standard) Belt Extension	26-70" (66-177 cm)	MLSO-2EXT


CONTOURED POSTERIOR PANEL			CONTOURED POSTERIOR PANEL WITH LATERAL SUPPORTS		
					
SIZE	WAIST CIRC.	PART #	SIZE	WAIST CIRC.	PART #
Small (S)	26-34" (66-86 cm)	MLSO-P11	Small (S)	26-34" (66-86 cm)	MLSO-P12
Medium (M)	34-38" (86-96 cm)	MLSO-P21	Medium (M)	34-38" (86-96 cm)	MLSO-P22
Large (L)	38-50" (96-127 cm)	MLSO-P31	Large (L)	38-50" (96-127 cm)	MLSO-P32
Extra Large (XL)	50-66" (127-167 cm)	MLSO-P41	Extra Large (XL)	50-66" (127-167 cm)	MLSO-P42

STANDARD POSTERIOR PANEL		NEW UNIVERSAL POSTERIOR		ANTERIOR PANEL	
					
SIZE	PART #	TYPE	PART #	SIZE	PART #
Standard	MLSO-P56	with Lateral	MLSO-P55	6" Panel (Slim)	MLSO-PA6
		without Lateral	MLSO-P57	8" Panel (Standard)	MLSO-PA8

OPTION 2: Order a brace in a kit.

WITH NEW UNIVERSAL POSTERIOR +/- ANTERIOR	
	
Belt Width (X)	X = 1 (6" Slim) or 2 (8" Standard)
Posterior Panel (Y)	Y = 5 (Universal)
Posterior Panel Type (Z)	Z = 5 (Universal w/ lats) or 7 (Universal no lats)
Anterior Option (A)	A = 1 (Anterior panel) or 2 (No Anterior)

ÖSSUR FIT TOOL	
	
PART #	B-04010000

WITH CONTOURED OR STANDARD POSTERIOR +/- ANTERIOR	
	
Belt Width (X)	X = 1 (6" Slim) or 2 (8" Standard)
Posterior Panel Size (Y)	Y = 1 (S), 2 (M), 3 (L), 4 (XL) or 5 (Standard)
Posterior Panel Type (Z)	Z = 1 (Contoured), 2 (Contoured w/ lats), or 6 (Standard)
Anterior Option (A)	A = 1 (Anterior Panel) or 2 (No Anterior Panel)

Part number =
 MLSO - -
 X Y Z A

X = Belt
 Y = Posterior Panel Size
 Z = Posterior Panel Type
 A = Anterior Panel Option

Note: The standard posterior panel is always recommended when no other posterior panel is being used.

MIAMI TLSO 464

The Miami TLSO from Össur directly addresses the needs of patients and practitioners, simplifying preparation without compromising on fit or performance. This easy-to-use pre-assembled Miami TLSO brace offers three configurations in one. Compliance: Easily applied in over the shoulder or underarm configuration, donning and doffing is further simplified with the conveniently located quick release buckles. Adjustability: The universal belt designs allows for setting the belt length on the patient for a more accurate sizing experience. Comfort: Truly shaped to fit the anatomy, the anatomically contoured design better immobilizes the spine while maximizing patient comfort. Ease of Use: The Miami TLSO brace requires minimal preparation and is designed to provide optimal patient comfort.



INDICATIONS FOR USE


Conditions requiring gross immobilization of the thoracic and lumbar spine. Such conditions may include:

- Post-surgical immobilization
- Spinal Support TLSO attachment is up to T7
- Spinal stenosis
- Herniated disc
- Degenerative spinal fractures
- Spondylolithesis
- Spondylolysis

CONTRAINDICATIONS

- Unstable, displaced fractures

REQUIRED MEASUREMENTS



Measure circumference of the trunk at a level midway between the lowest ribs and the upper iliac crest. Standing.

Size	Waist circumference	
	cm	inch
Universal	66 - 177	26 - 70
S	66 - 86	26 - 34
M	86 - 107	34 - 42
L	107 - 127	42 - 50
XL	127 - 177	50 - 70

*Belt extensions are required for waist circumferences larger than 127 cm / 50"

FEATURES AND BENEFITS

- Significantly superior comfort and ease of donning/doffing.
- Fully assembled for fitting with ease.
- Accommodates heights 155–193 cm (5'1"–6'5")
- Posterior thoracic extension option with ergonomic axillary strapping.
- Step-down to an LSO as one progresses through rehabilitation.
- One-touch telescoping anterior-thoracic extension option with single or double-hinge.
- The convenience of a one-size Universal or individualize with a uniquely-sized anatomical design.
- Improve comfort and immobilization with multiple strapping configurations and removable anchors.

MIAMI TLSO 464

ORDERING INFORMATION

ATE/PTE, LORDOTIC

Size	Part#
Universal	MTLSO-2351
S	MTLSO-2111
M	MTLSO-2121
L	MTLSO-2131
XL	MTLSO-2141

ATE/PTE, NON-LORDOTIC

Size	Part#
S	MTLSO-2211
M	MTLSO-2221
L	MTLSO-2231
XL	MTLSO-2241

ATE/PTE, LORDOTIC, SLIM

Size	Part#
Universal	MTLSO-1351
S	MTLSO-1111
M	MTLSO-1121
L	MTLSO-1131
XL	MTLSO-1141

ATE/PTE, NON-LORDOTIC, SLIM

Size	Part#
S	MTLSO-1211
M	MTLSO-1221
L	MTLSO-1231
XL	MTLSO-1241

OTHER ACCESSORIES

Description	Part#
5MM Hex Long Wrench	MTLSO-0004
Axillary strap accessory kit	MTLSO-0003
Quick fit tool	B-040100000
Belt extension, standard (Recommended for waist circumferences 127-177 cm (50"-70"). Fits size XL.)	MLSO-2EXT
Belt extension, slim (Recommended for waist circumferences 127-177 cm (50"-70"). Fits size XL.)	MLSO-1EXT

MIAMI TLSO 456

The Miami TLSO 456 from Össur directly addresses the needs of patients and practitioners, simplifying preparation without compromising on fit or performance. This easy-to-use pre-assembled Miami TLSO brace offers the convenience of a one-size Universal or a uniquely-sized anatomical design.

Compliance: Blue is You™ patient touch-points help send patients home with confidence.

Adjustability: The universal belt design allows for setting the belt length on the patient for a more accurate sizing experience. Heat moldable posterior panels allow a customized fit, while the modular design allows step-down support to an LSO.

Comfort: Featuring the Miami LSO™ Advanced Compression pulley system with lumbar compression you can feel, and an anatomically contoured design that better immobilizes the spine while maximizing patient comfort.

Ease of Use: Easily applied with over-the-shoulder straps, the Miami TLSO brace requires minimal preparation and is designed to provide optimal patient comfort.



INDICATIONS FOR USE

Conditions requiring gross immobilization of the thoracic and lumbar spine. Such conditions may include:

- Post-surgical immobilization
- Spinal Support TLSO attachment is up to T7
- Spinal stenosis
- Herniated disc
- Degenerative spinal fractures
- Spondylolithesis
- Spondylolysis


FEATURES AND BENEFITS

- Significantly superior comfort and ease of donning/ doffing (when compared to Miami Lumbar TLSO).
- Fully assembled for fitting with ease.
- Accommodates heights 155–193 cm (5'1"–6'5")
- Step-down to an LSO as one progresses through rehabilitation.
- The convenience of a one-size Universal or individualize with a uniquely-sized anatomical design.

CONTRAINDICATIONS

- Unstable, displaced fractures

REQUIRED MEASUREMENTS



Measure circumference of the trunk at a level midway between the lowest ribs and the upper iliac crest. Standing.

Size	Waist circumference*	
	cm	inch
Universal	66 - 177	26 - 70
S	66 - 86	26 - 34
M	86 - 107	34 - 42
L	107 - 127	42 - 50
XL	127 - 177	50 - 70

*Belt extensions are required for waist circumferences larger than 127 cm / 50"

MIAMI TLSO 456

ORDERING INFORMATION

PTE LORDOTIC

Size	Part#
Universal	MTLSO-2352
S	MTLSO-2112
M	MTLSO-2122
L	MTLSO-2132
XL	MTLSO-2142

PTE NON-LORDOTIC, SLIM

Size	Part#
S	MTLSO-1212
M	MTLSO-1222
L	MTLSO-1232
XL	MTLSO-1242

PTE NON-LORDOTIC

Size	Part#
S	MTLSO-2212
M	MTLSO-2222
L	MTLSO-2232
XL	MTLSO-2242

PTE LORDOTIC, SLIM

Size	Part#
Universal	MTLSO-1352
S	MTLSO-1112
M	MTLSO-1122
L	MTLSO-1132
XL	MTLSO-1142

OTHER ACCESSORIES

Description	Part#
Quick fit tool	B-040100000
Belt extension, standard (Recommended for waist circumferences 127-177 cm (50"-70"). Fits size XL.)	MLSO-2EXT
Belt extension, slim (Recommended for waist circumferences 127-177 cm (50"-70"). Fits size XL.)	MLSO-1EXT
Axillary strap accessory kit	MTLSO-0003
5MM Hex Long Wrench	MTLSO-0004

MIAMI LUMBAR® POSTEO

Designed with patient comfort and compliance in mind, the Miami Lumbar Posteo provides functional relief of vertebral fracture pain, and facilitates return to mobility. The lumbar compression system offers both pain relief via active vertebral offloading and better grip of patient's hips for greater ease of use and compliance. The Miami Lumbar Posteo includes an innovative approach to ease of use and consistent donning.



INDICATIONS FOR USE

- Back pain caused by stable vertebral compression fractures (VCF)

FEATURES AND BENEFITS

- Active offloading via thoracic extension and lumbar compression
- Integrated pulley system for customized fit
- Trimmable axillary straps with adjustable tension for optimal stabilization
- Cold-moldable aluminum PTE to fit different anatomies
- Extra cushioning around the aluminum strut for minimized skin pressure

REQUIRED MEASUREMENTS

- Back length measured from the top of the sacrum up to C7 vertebrae, following the spinal curvature

Please note that the Miami Lumbar Posteo includes a universal lumbar belt which can be trimmed down to the appropriate size. Calculation of the appropriate back length size, measured from the tip of the sacrum to the C7 vertebra, following the curvature of the spine.

ORDERING INFORMATION

MIAMI LUMBAR POSTEO

Size	Back Length	Part#
X-Small	14.5" – 17" (37 – 44 cm)	LMB-POSTEO-XS
Small	17" – 20" (44 – 51 cm)	LMB-POSTEO-SM
Medium	19" – 22" (49 – 56 cm)	LMB-POSTEO-MD
Large	21" – 24" (55 – 61 cm)	LMB-POSTEO-LG

ACCESSORIES

Description	Part#
Strap Pad	LMB-POSTEO ST PD



ÖSSUR FORMFIT® PRO BACK

Össur Formfit® Pro Back provides relief for patients suffering from pains and strains in the lower back. Manufactured with our dynamic, 3D-Motion Tech™ knitting technology, our belts feature superior moisture wicking capabilities and are highly breathable. The patented back panel design utilizes our StayBilizer™ technology to provide optimum shape and fit, while conveniently placed finger loops make donning easy. The SmartPad™ provides a targeted support and cushioning and can be removed if desired.



SPINE SOLUTIONS

INDICATIONS FOR USE

For lower back conditions that may benefit from compression and gross motion restriction, such as:

- Lower back pain (Lumbago)
- Lumbosacral sprain/strain
- Degenerative joint/disc disease


CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated.

FEATURES AND BENEFITS

- MotionTech™ for compression and dynamic fit.
- Removable lumbar SmartPad™ for targeted support and cushioning.
- Lightweight, malleable aluminum stays are shaped to the anatomy of the lower back.

REQUIRED MEASUREMENTS



Measure circumference of the trunk at a level midway between the lowest ribs and the upper iliac crest. Standing.

Size	Waist circumference	
	cm	inch
XS	70 - 80	27½ - 31½
S	80 - 90	31½ - 35¾
M	90 - 100	35¾ - 39¾
L	100 - 115	39¾ - 45¼
XL	115 - 130	45¼ - 51¾

ORDERING INFORMATION

BLACK

Size	Part#
XS	FP500001
S	FP500002
M	FP500003
L	FP500004
XL	FP500005

MIAMI J® SELECT

The Miami J Select cervical collar by Össur offers the widest range of height adjustability in the market with a single, universal design that replicates the comfort, immobilization and quality you expect from a Miami J product. Featuring a patented, first-to-market locking mechanism to encourage patient compliance and anti-microbial Sorbatex padding, Miami J Select is designed to support successful patient outcomes.



INDICATIONS FOR USE

Situations requiring gross immobilization of the cervical spine, these may include:

- Pre- and post c-spine surgery
- C-spine precaution for trauma patients

FEATURES AND BENEFITS


- Easy height adjustability with a tamper-proof locking mechanism – ‘Left to Lock’.
- Anti-microbial & biocompatible Sorbatex™ padding prevents bacterial build-up.
- Universal design that doesn’t compromise immobilization or comfort.
- 7 precise fitting sizes with opportunity to fine-tune.
- Miami J® inspired chin tray.
- One-piece, anatomically inspired collar back with Flex Edge® designed to reduce occipital pressure.
- Easy to don and doff with patented Reproducible Fit option.
- Intuitive blue Patient Touch-Points to send patients home with confidence – Blue is You™.
- Simplified comfort adjustments with integrated and patented Sternal-Relief Dial.
- Large tracheal opening – improved function for ICU patients.
- Easy to clean and replace multi-colored padding.
- X-ray and CT lucent, MR-safe.
- Compatible with Occian™ Back & Miami JTO® accessories.

ORDERING INFORMATION

Description	Part#	Part#
Miami J Select Collar Set	One-size adjustable (Universal)	MJSR-101
Miami J Select Collar	One-size adjustable (Universal)	MJS-101
Miami J Select Replacement Pads	One-size adjustable (Universal)	MJSP-101
Miami J Select Reproducible Fit Tabs (10 pack)	One-size adjustable (Universal)	MJST-101

MIAMI J®

Known as the superior c-spine immobilizer through multiple independent studies, the Miami J's patented design meets the specialized needs of all those afflicted with cervical injuries. The Miami J offers features that enhance compliance while providing immobilization, inhibition of skin breakdown and patient comfort.



SPINE SOLUTIONS

INDICATIONS FOR USE

- Post-trauma
- Stable cervical fractures
- Cervical spondylitis
- Cervical spondylosis
- RA and OA cervical spine
- Herniated cervical disc
- Post-operative
- Motor Neurone Disease

FEATURES AND BENEFITS

- Sorbatex™ padding is antibacterial and clinically proven to inhibit microbial growth, which enhances comfort and helps assure patient compliance
- Bioengineered to minimize pressure points in key known areas: Chin, occiput, trapezius and clavicle
- Angulation buttons for occiput and sternal pad is removable when supine, lying and swallowing
- MRI, CT and X-ray lucent, MR safe

REQUIRED MEASUREMENTS

- Sizing based on 6 phenotypes to ensure correct fit and optimum immobilization

Phenotype-driven sizing system makes size selection simple while drastically reducing margin of error in fitting a patient. Each collar has a bioengineered front with a specific back designed to work together for superior immobilization. For pediatric sizing, please see Miami Jr in this Catalog.

Specialty Sizes (40% of population)

Standard Sizes (60%)

Kyphotic (10%)	XSmall (5%)	Stout (15%)	Tall (10%)	Regular (30%)	Short (30%)
<ul style="list-style-type: none"> • kyphotic • chin-on-chest • geriatric • osteoporitic and/or ankylosing spondylitic 	<ul style="list-style-type: none"> • Very short, thin neck • Principally female • More prevalent in the Asian population 	<ul style="list-style-type: none"> • Very large neck circumference • Obese • "No neck" • Massive shoulders 	<ul style="list-style-type: none"> • Long, tall neck • "Swan neck" • Young Women • Adolescents 	<ul style="list-style-type: none"> • Female • Standard size for mature women • Thin, mature male 	<ul style="list-style-type: none"> • Standard size for mature men • Short-necked women

*For individuals with neck circumference > 25", the Bariatric collar is recommended.

MIAMI J®

ORDERING INFORMATION

MIAMI J

Size	Collar set*	Collar	Collar Back	Collar Front
Super Short	MJR-200S	MJ-200S	MJ-400BK	MJ-200SFRT
Stout	MJR-200L	MJ-200L	MJ-200LBK	MJ-200LFRT
XS	MJR-250	MJ-250	MJ-250BK	MJ-250FRT
Short	MJR-300	MJ-300	MJ-300BK	MJ-300FRT
Regular	MJR-400	MJ-400	MJ-400BK	MJ-400FRT
Tall	MJR-500	MJ-500	MJ-500BK	MJ-500FRT
Bariatric	MJR-200B	MJ-200B		

*Collar sets come with one additional set of replacement pads.

Miami J Front and Back Panels may be interchanged to accommodate patient anatomy, e.g. a Short Front may be used with a Tall Back. Choose collar size based on patient neck height, not patient body size. Super Short collars are packaged with Regular Back Panels

ACCESSORIES

REPLACEMENT PADS

Size	Part#
Super Short	MJP-100
Stout	MJP-200L
XS	MJP-250
Short	MJP-100
Regular	MJP-100
Tall	MJP-100
Bariatric	MJP-200B

MIAMI J PATIENT PACK

Description	Part#
Miami J Patient Pack	MJ-PTPK2

MIAMI JTO® THORACIC EXTENSION

The Miami JTO Thoracic Extension is designed in conjunction with leading spine surgeons, orthotists and nurses. The orthosis is used exclusively with all Miami J® collars when even greater stability at the cervical junctions is indicated.



SPINE SOLUTIONS

INDICATIONS FOR USE

- | | |
|--|--|
| <p>Lower C-Spine disruptions</p> <ul style="list-style-type: none"> • Level one instabilities • Isolated compartment fracture • Lateral mass fracture • Spinous Process fracture | <p>Upper C-Spine disruptions</p> <ul style="list-style-type: none"> • For increased stability over the collar alone • Stable Jefferson or Hangman’s Fracture • Type I Odontoid Fracture |
|--|--|

FEATURES AND BENEFITS

- Provides more control than a collar alone can provide
- Innovative device designed to treat cervical and high-thoracic injuries
- Patented design with ‘floating back’ permits maximum immobilization without the risk of compensatory motion
- Snap on and off without tools
- Sorbatex™ II padding is antibacterial and clinically proven to inhibit microbial growth, which enhances comfort and helps assure patient compliance
- MRI, CT and X-ray lucent, MR compatible

REQUIRED MEASUREMENTS

- Universal, one size

ORDERING INFORMATION

MIAMI JTO THORACIC EXTENSION

Description	Part#
Miami JTO	MT-900
Miami JTO with replacement softgoods	MTR-900

Compatible with all sizes of the Miami J®, level 1 through 5 of the Miami J® Select, and the P3 Miami Jr®. Any model of Miami J collar is not included with Miami JTO Thoracic Extension and is purchased separately.

ACCESSORIES

Description	Part#
Replacement softgoods	MTP-950

OCCIAN® BACK & OCCIAN® BACK JR

The Occian Back is designed exclusively for the specialized needs of the multi-system trauma patient in the ICU. Universally sized, specifically engineered and designed with advanced, Intuitech™ pressure-relieving memory foam, the Occian Back helps caregivers manage c-spine patients who are at high risk of occipital breakdown.



INDICATIONS FOR USE

- Patients not awake and alert within 24 – 48 hours after admission
- Patients with pre-existing skin injuries and/or lacerations to the head/neck
- Burn patients who must wear a collar
- Elective surgery patients desiring night time post-operative comfort
- Patients who remain in a supine position for long periods of time related to respiratory, neurological, or haemodynamic instability

FEATURES AND BENEFITS

- Intuitech pressure relieving memory foam virtually eliminates occipital breakdown
- Maintains C spine stability
- Provides unique pressure-relieving comfort by replacing the back section of a Miami J Collar.
- Clinically proven to virtually eliminate occipital breakdown without sacrificing c-spine stability.
- Easy to clean
- MRI, CT & X-ray lucent, MR safe

REQUIRED MEASUREMENTS

- Universal, one size

ORDERING INFORMATION

Description	Part#
Occian Back*	ACB-100
Occian Back Junior**	ACB-200

*Compatible with Miami Jr® size P3 and all sizes of Miami J® and Miami J® Select.

** Compatible with Miami Jr® size P2.

MIAMI Jr®

The Miami Jr is the breakthrough collar designed exclusively for children of 12 years and under. The Miami Jr line, which spans the growth cycle, is sized according to age. Össur is the only company that has a pediatric line of collars based on the anatomical differences between adults and children. We realize that children are not just “little adults” – they require innovative design.



SPINE SOLUTIONS

INDICATIONS FOR USE

- Post-trauma
- Stable cervical fractures
- Cervical spondylitis
- Cervical spondylosis
- RA and OA cervical spine
- Herniated cervical disc
- Post-operative
- Motor Neurone Disease

FEATURES AND BENEFITS

- Innovative collar designed exclusively for children 12 years and under anatomy pediatric
- Sorbatex™ II pads make the collar “skin friendly” and protect skin integrity during extended wear
- MRI, CT and X-ray lucent, MR safe

REQUIRED MEASUREMENTS

- Age

ORDERING INFORMATION

MIAMI JR

Size	Collar set*	Collar
0 – 6 months	MJR-P0	MJ-P0
6 – 24 months	MJR-P1	MJ-P1
2 – 6 years	MJR-P2	MJ-P2
6 – 12 years	MJR-P3	MJ-P3

*Collar sets come with one additional set of replacement pads.

ACCESSORIES

REPLACEMENT PADS

Size	Part#
0 – 6 months	MJP-P0
6 – 24 months	MJP-P1
2 – 6 years	MJP-P2
6 – 12 years	MJP-P3

PAPOOSE®

The Papoose is an innovative design for the youngest and most challenging patients – newborn infants. The Papoose cradles the baby’s body securely while gently immobilizing the head and spine.



INDICATIONS FOR USE

- Immobilization of the spine
- Spinal fractures

FEATURES AND BENEFITS

- True cervical spine immobilizer for infants up to 3 months
- Allows for continued intimacy between caregiver and child
- MRI, CT and X-ray lucent, MR compatible

REQUIRED MEASUREMENTS

- Universal, one size

ORDERING INFORMATION

PAPOOSE

Description	Part#
Papoose without front collar	PO-100
Papoose without front with replacement pads	POR-100
Papoose with front collar	PA-100
Papoose with front with replacement pads	PAR-100

ACCESSORIES

REPLACEMENT PADS

Description	Part#
Papoose replacement pads (with front)	PAP-100
Papoose replacement pads (without front)	POP-100

PHILADELPHIA® ADJUSTABLE TRACHEOTOMY COLLAR

INDICATIONS FOR USE

- Post-trauma
- Herniated cervical disc
- Cervical spondylitis
- Post-operative
- Cervical spondylosis
- RA and OA cervical spine

FEATURES AND BENEFITS

- Circumference and height adjustable
- Permanent lock on height adjustment
- Anatomically shaped Plastazote® foam (Not made with natural rubber latex, non-toxic and hypoallergenic to reduce skin irritation), leaving clearance for cervical veins
- Total cervical arch support maintains neutral alignment of the cervical spine
- Water-resistant for bathing and hydro-therapy.
- Well ventilated
- X-ray and CT lucent. MR conditional.



SPINE

REQUIRED MEASUREMENTS

- Neck height: 1 inch (or 25mm) of neck length adjustability accommodates varying patient neck length.

ORDERING INFORMATION

PHILADELPHIA ADJUSTABLE TRACHEOTOMY COLLAR

Size	Neck circumference
Small (2 tears)	10 – 13" (25.4 – 33 cm)
Medium (1 tear)	13 – 16" (33 – 40.6 cm)
Large (no tear)	16 – 19" (40.6 – 48.3 cm)

NOTE: The table above represents common sizing only. Sizing must be verified by the healthcare provider for each patient individually.

Part#	Size / Description
PHP-AA	Conforms to sizes 2, 3 and 4 of the original Philadelphia Tracheotomy Collar

PHILADELPHIA® TRACHEOTOMY COLLAR

INDICATIONS FOR USE

- Post-trauma
- Herniated cervical disc
- Cervical spondylitis
- Post-operative
- Cervical spondylosis
- Motor Neurone Disease
- RA and OA cervical spine

FEATURES AND BENEFITS

- Not made with natural rubber latex, non-toxic and hypoallergenic Plastazote® foam reduces skin irritation
- Total cervical arch support maintains neutral alignment of the cervical spine
- Water-resistant for bathing and hydro-therapy.



REQUIRED MEASUREMENTS

- Neck circumference and height from tip of chin to sternal notch

ORDERING INFORMATION

Size	Neck circumference	6 cm / 2.25" height	8 cm / 3.25" height	11 cm / 4.25" height	13 cm / 5.25" height
X-Small	up to 10" (< 25 cm)	PHP-T2XS			
Small	10 – 13" (25 – 33 cm)	PHP-T2S	PHP-T3S	PHP-T4S	PHP-T5S
Medium	13-16" (33 – 41 cm)	PHP-T2M	PHP-T3M	PHP-T4M	PHP-T5M
Large	16-19" (41 – 48 cm)	PHP-T2L	PHP-T3L	PHP-T4L	PHP-T5L
X-Large	19" up (48 cm)	PHP-T2XL	PHP-T3XL	PHP-T4XL	PHP-T5XL

Size	Neck circumference	Height	Part#
Infant	6-8" (15 – 20 cm)	1.25" (3 cm)	PHP-TI
Pediatric	8-11" (20 – 28 cm)	1.75" (4.5 cm)	PHP-TP

ACCESSORIES

PHILADELPHIA COOLMAX COTTON LINER

Size	Neck circumference	Part#
Infant	6-8" (15 – 20 cm)	PHP-LCI
X-Small	up to 10" (< 25 cm)	PHP-LCXS
Small	10-13" (25 – 33 cm)	PHP-LCS
Medium	13-16" (33 – 41 cm)	PHP-LCM
Large	16-19" (41 – 48 cm)	PHP-LCL
X-Large	19" up (> 48 cm)	PHP-LCX

PHILADELPHIA TRACHEOTOMY ATTACHMENTS

Size	Stabilizer	Halo
Pediatric	PHP-OSP	PHP-OHP
Adult	PHP-OS	PHP-OH

UNIVERSAL CERVICAL COLLAR

INDICATIONS FOR USE

- Neck extension injuries (Whiplash)
- Cervical spondylitis
- Cervical spondylosis
- RA and OA cervical spine
- Herniated cervical disc

FEATURES AND BENEFITS

- Medium density foam
- Contact closure for easy fit

REQUIRED MEASUREMENTS

- Mandible to clavicle height



ORDERING INFORMATION

Height	Part#
2" (5 cm)	203020
3" (7.5 cm)	203000
4" (10 cm)	203010

FOAM CERVICAL COLLAR

INDICATIONS FOR USE

- Neck extension injuries (Whiplash)
- Cervical spondylitis
- Cervical spondylosis
- RA and OA cervical spine
- Herniated cervical disc

FEATURES AND BENEFITS

- Easy to apply

REQUIRED MEASUREMENTS

- Neck circumference and mandible to clavicle height
- Ordering Information



Size	Circumference	Neck size	Part#
Small	3" x 18.5" (7 cm x 47 cm)	11-15"	203013
Small/Long	3" x 22" (7 cm x 56 cm)	15-20"	203014
Medium	3.5" x 20" (9 cm x 51 cm)	13-18"	203015
Medium/Long	3.5" x 23" (9 cm x 58 cm)	16-21"	203016
Large	4" x 23" (10 cm x 58 cm)	15-20"	203017
Large/Long	4" x 25" (10 cm x 64 cm)	16-22"	203018

NECLOC®

INDICATIONS FOR USE

- Post-trauma

FEATURES AND BENEFITS

- Ergonomic two-piece clamshell design and retaining strap secure flexion during application
- Five unique phenotype-driven sizes make it easy to select the correct size
- Minimal manipulation of patient during application
- X-ray and CT lucent. MR conditional.



ORDERING INFORMATION

Size	Part#
Stout	NL-200EL
X-Small	NL-250E
Small	NL-300E
Medium	NL-400E
Large	NL-500E

NECLOC® KIDS

INDICATIONS FOR USE

- Post-trauma

FEATURES AND BENEFITS

- Ergonomic two-piece clamshell design and retaining strap secure flexion during application
- Three pediatric sizes ensure optimal airway and cervical spine immobilization
- Minimal manipulation of patient during application to children
- Meets the unique developing anatomical requirements of children
- X-ray and CT lucent. MR conditional



REQUIRED MEASUREMENTS

- Age

ORDERING INFORMATION

Age	Part#
0 – 2 years	NL-P1
2 – 6 years	NL-P2
6 – 12 years	NL-P3

PATRIOT® COLLAR

INDICATIONS FOR USE

- Post-trauma

FEATURES AND BENEFITS

- Multi-height adjustment – Adult (4 positions) and Pediatric (3 positions)
- Open rear ventilation panel for fluid drainage, airflow and cervical spine palpation
- Universal sizing optimizes inventory requirements

REQUIRED MEASUREMENTS

- Neck circumference



SPINE

ORDERING INFORMATION

Size	Neck circumference	Blue
Child	8-18" (20 – 46 cm)	PHP-PP
Adult	11-23" (28 – 58 cm)	PHP-PA

OAM RIGID LUMBAR

From the makers of the Miami Lumbar,[®] the OAM Rigid Lumbar offers you a discrete look with endless options.

First, you choose the belt.

Slim or regular, 10 pre-cut sizes or universal.

Then, you choose the panels.

Different panels for different control, modular to allow you to remove as the patient heals.

Or, you choose the kit.

We make it even simpler by providing our most popular sizes pre-configured for you.



TWO SIMPLE OPTIONS FOR ORDERING

OPTION 1 TO ORDER BRACES BY COMPONENTS

OAM RIGID LUMBAR BELT			
SIZE	WAIST CIRC.	REGULAR PART #	SLIM PART #
XS	26-30" (66-76cm)	LMBGEN-001	LMBSLGEN-001
S	30-34" (76-86cm)	LMBGEN-002	LMBSLGEN-002
M	34-38" (86-96cm)	LMBGEN-003	LMBSLGEN-003
L	38-42" (96-106cm)	LMBGEN-004	LMBSLGEN-004
XL	42-46" (106-116cm)	LMBGEN-005	LMBSLGEN-005
2XL	46-50" (116-127cm)	LMBGEN-006	LMBSLGEN-006
3XL	50-54" (127-137cm)	LMBGEN-007	LMBSLGEN-007
4XL	54-58" (137-147cm)	LMBGEN-008	LMBSLGEN-008
5XL	58-62" (147-157cm)	LMBGEN-009	LMBSLGEN-009
6XL	62-66" (157-167cm)	LMBGEN-010	LMBSLGEN-010

ANTERIOR PANELS		
SIZE	WAIST CIRC.	PART #
Standard	<46" <116cm	LMBGEN-ST
Slim	<46" <116cm	LMBGEN-SLST
Pendulous	>46" >116cm	LMBGEN-PN

POSTERIOR PANELS			
SIZE	WAIST CIRC.	WITHOUT LATERAL	WITH LATERAL
S	26-34" (66-86cm)	LMBGEN-NS	LMBGEN-LS
M	34-38" (86-96cm)	LMBGEN-NM	LMBGEN-LM
L	38-50" (96-127cm)	LMBGEN-NL	LMBGEN-LL
XL	50-66" (127-167cm)	LMBGEN-NX	LMBGEN-LX

OPTION 2 TO ORDER BRACES BY KIT

FILL IN SPACES TO COMPLETE PART NUMBER

SLIM CONFIGURATION
PART NUMBER = OAM

STANDARD CONFIGURATION
PART NUMBER = OAMST

XXX = Panel Configuration
 YY = Panel Size/ Belt Length

1. CHOOSE PANEL CONFIGURATION X X X

L0627	L0631	L0637
XXX= 627	XXX= 631	XXX= 637

2. CHOOSE A BELT LENGTH & PANEL SIZE

Use waist circ. for panel/belt size		Y Y	Y Y			
BELT LENGTH	WAIST CIRC.	ANTERIOR PANEL (627 ONLY)	POSTERIOR PANEL SIZE (631 & 637)			
			SM	MD	LG	XL
Universal	34-50" (86-127cm)	01	01	06	16	18
Universal +	34-66" (86-167cm)	06	02	10	20	22

Note: 631 & 637 configurations are kitted with their respective anterior panel by default.

FORMFIT® BACK SUPPORT, AIR

INDICATIONS FOR USE

- Muscular low back pain
- Lumbosacral Strains
- Chronic low back Pain
- Sciatica

FEATURES AND BENEFITS

- Intuitive integrated pneumatic system for more customized compression and support
- Anti-Migration system reduces slippage
- Double wrap design for improved compression
- Finger loops for easy application
- Hot/Cold Gel pad available separately



SPINE

REQUIRED MEASUREMENTS

- Waist circumference

ORDERING INFORMATION

Size	Waist circumference	Standard
Small	25-30" (63 – 76 cm)	209133
Medium	30-34" (76 – 86 cm)	209135
Large	34-38" (86 – 97 cm)	209137
X-Large	38-42" (97 – 107 cm)	209138
XX-Large	42-56" (107 – 142 cm)	209139

ACCESSORIES

Description	Part#
Gel pad	GB19028

Description	Part#
Back Support Suspenders Clip On	10912

FORMFIT® BACK SUPPORT

INDICATIONS FOR USE


- Lumbar strains
- General posture support
- Lumbosacral strains
- Sciatica

FEATURES AND BENEFITS

- Muscular low back pain
- Lumbosacral strains
- Chronic low back pain
- Sciatica

REQUIRED MEASUREMENTS

- Waist circumference



ORDERING INFORMATION

Size	Waist Circumference	Standard
Small	25 – 30" (63 – 76 cm)	209143
Medium	30 – 34" (76 – 86 cm)	209145
Large	34 – 38" (86 – 97 cm)	209147
X-Large	38 – 42" (97 – 107 cm)	209148
XX-Large	42 – 56" (107 – 142 cm)	209149
XXX-Large	56 – 60" (142 – 152 cm)	209151

ACCESSORIES

Description	Part#
Back Support Suspenders Clip On	10912
Hot/Cold Therapy Gel Pad	GB19028

UNIVERSAL ELASTIC RIB BELT

INDICATIONS FOR USE

- Pre-or post operative use in the thorax area.

FEATURES AND BENEFITS

- Flannel-lined
- Contact closure

REQUIRED MEASUREMENTS

- Measure circumference at waist



ORDERING INFORMATION

Description	Part#
Male – 6.25" (16 cm) wide x 38.75" (98.5 cm) long	209150
Female – 6.25" (16 cm) wide x 38.75" (98.5 cm) long	209160

UNIVERSAL ABDOMINAL BINDER

INDICATIONS FOR USE

- Sacroiliac support
- Low back pain

FEATURES AND BENEFITS

- Elastic belt has Velcro closures for a perfectly contoured fit
- Available in unisex universal size

REQUIRED MEASUREMENTS

- Measure circumference at waist

ORDERING INFORMATION

Description	Fits	Part#
10" Universal binder	26 – 50" (66 – 127 cm)	SS1100

ELASTIC ABDOMINAL BINDER

INDICATIONS FOR USE

- Sacroiliac support
- Low back pain

FEATURES AND BENEFITS

- Durable Elastic
- Abdominal, post-natal and post-surgical application
- Easy adjustment with contact closure
- Provides compression and support for strains and weakness of the abdominal area
- 3 panel and 4 panel options available
- Does not contain natural rubber latex



REQUIRED MEASUREMENTS

- Measure waist circumference

ORDERING INFORMATION

Waist size	3-Panel	4-Panel
30 – 42" (76 – 107 cm)	5437691	54376121
42 – 55" (107 – 140 cm)	5437692	54376122
55 -69" (140 – 175 cm)	5437693	54376123
69 -84" (175 – 213 cm)	5437694	54376124
84 – 100" (213 – 254 cm)	5437695	54376125

BACK SUPPORT SUSPENDERS

FEATURES AND BENEFITS

- Suspenders can be added to each support upon request.

REQUIRED MEASUREMENTS

- Please specify sewn-in or removable suspender design.



ORDERING INFORMATION

Description	Part#
Clip-on	10912



ÖSSUR[®]
LIFE WITHOUT LIMITATIONS

UPPER EXTREMITY



FORMFIT® PRO ELBOW

Össur Formfit® Pro Elbow provides relief from elbow pain. This 3D compressive elbow sleeve is manufactured using our MotionTech™ technology and features two integrated epicondyle SmartPads™ for more targeted support and cushioning. The RadialStrap™ provides additional compression and support.



MotionTech™



CoolVent™



SmartPad™



FoldLock™

UPPER EXTREMITY SOLUTIONS

INDICATIONS FOR USE

For elbow conditions that may benefit from compression, such as:

- Medial and lateral epicondylitis
- Elbow pain and swelling

CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated.

FEATURES AND BENEFITS

- 3D MotionTech™ for compression and dynamic fit keep.
- Epicondyle SmartPads™ for targeted support and cushioning.
- Patented FoldLock™ system for easier and more secure application.
- RadialStrap™ for increased compression and support.

REQUIRED MEASUREMENTS



Measure arm circumference exactly at the elbow. Elbow in full extension.



Size	Elbow circumference		Blue	Black
	cm	inch	Part#	Part#
XS	17 - 20	6¾ - 7⅞	FP301001	FP300001
S	20 - 23	7⅞ - 9	FP301002	FP300002
M	23 - 26	9 - 10¼	FP301003	FP300003
L	26 - 29	10¼ - 11¾	FP301004	FP300004
XL	29 - 32	11¾ - 12⅝	FP301005	FP300005

INNOVATOR X®

Benefiting from our proven, patented technologies, Innovator X, with extending struts, can help improve patient outcomes. With superior comfort, clean aesthetics and the most user-friendly hinge in post-op bracing, Innovator X is designed to increase patient satisfaction and compliance.



INDICATIONS FOR USE

- Controlled elbow flexion/extension following surgery or trauma

FEATURES AND BENEFITS

- Easy-to-use Innovator hinge
- Slide-to-Size™ struts lead to a more customized fit
- Sensil® silicone padding on the strapping system reduces slipping and malleable aluminium cuffs allow for customization
- Patented Formfit® 3-dimensional padding and lycra lining provide superior comfort and compliance
- Comfort pad at the neck and shoulder
- Flexion range can be quickly adjusted to 0, 10, 20, 30, 45, 60, 75, 90, 105 and 120. Extension range can be adjusted to 0, 10, 20, 30, 45, 60, 75 and 90.
- Optional anti-supination/pronation arm bar kit is available

REQUIRED MEASUREMENTS

- Universal, one size

ORDERING INFORMATION

INNOVATOR X

Side	Part#
Left	W-823101
Right	W-823100

ACCESSORIES

Description	Part#
Innovator X Bar Kit	W-A82310
Replacement Softgoods Left	W-823111
Replacement Softgoods Right	W-823110

FORMFIT® SHOULDER BRACE

The Formfit Shoulder Brace is designed to provide immobilization following shoulder injury/surgery. With its streamlined fit, sporty look and breathability it offers excellent comfort and ease of use. Same as its predecessor, the SmartSling®, the Formfit Shoulder Brace also provides abduction of the shoulder (when the abduction pillow is attached), or external rotation of the shoulder (when the external rotation pillow is attached).

Optional axilla pillow helps distribute weight over a larger surface area and decreases pressure on the ulnar nerve.

INDICATIONS FOR USE

- Immobilization for rotator cuff repairs
- Bankart Repairs
- Capsular Shift
- Gleno humeral dislocation/subluxations

FEATURES AND BENEFITS

- Breathable material
- Quick-release buckles
- Offloading strap
- Extra padding
- Universal thumb loop

REQUIRED MEASUREMENTS

- Olecranon to end of closed fist



Axilla Pillow



External rotation pillow



Airplane Accessory



ORDERING INFORMATION

Size	Fits Length	Description	Part
Small	10.5-13" (up to 29 cm)	Formfit Shoulder Brace - sling only	SSS-SO
Medium	13.25-15.75" (29-40 cm)	Formfit Shoulder Brace - sling only	SSR-SO
Large	16-17.5" (40-44.5 cm)	Formfit Shoulder Brace - sling only	SSL-SO
Small	10.5-13" (up to 29 cm)	Formfit Shoulder Brace with Abduction	SSS
Medium	13.25-15.75" (29-40 cm)	Formfit Shoulder Brace with Abduction	SSR
Large	16-17.5" (40-44.5 cm)	Formfit Shoulder Brace with Abduction	SSL
Small	10.5-13" (up to 29 cm)	Formfit Shoulder Brace with External Rotation	SSS-ER
Medium	13.25-15.75" (29-40 cm)	Formfit Shoulder Brace with External Rotation	SSR-ER
Large	16-17.5" (40-44.5 cm)	Formfit Shoulder Brace with External Rotation	SSL-ER
Small	10.5-13" (up to 29 cm)	Formfit Shoulder Brace with Airplane Accessory	SSS-AK
Medium	13.25-15.75" (29-40 cm)	Formfit Shoulder Brace with Airplane Accessory	SSR-AK
Large	16-17.5" (40-44.5 cm)	Formfit Shoulder Brace with Airplane Accessory	SSL-AK

ACCESSORIES

Description	Part
Axilla Pillow	SSAXP
External Rotation Pillow	SSE
Airplane Accessory	SSAK

SHOULDER ABDUCTION SLING

INDICATIONS FOR USE

- Shoulder injuries
- Post-operative

FEATURES AND BENEFITS

- Provides shoulder abduction at 10° with support
- Cushion rests on waist to support weight of arm
- Breathable cotton provides a comfortable fit
- Fits left and right

REQUIRED MEASUREMENTS

- Measure elbow crease to fingertip



ORDERING INFORMATION

Size	Elbow crease to fingertip	Part#
Small	up to 12.5" (< 32 cm)	204053
Medium	12.5-15.5" (32 – 39 cm)	204055
Large	15.5" and up (39 cm +)	204057


PREMIUM SHOULDER IMMOBILIZER

INDICATIONS FOR USE

- Post-trauma
- Post-operative

FEATURES AND BENEFITS

- Comfortable fit
- Easily adjusts
- One size fits all (max circumference 45")



ORDERING INFORMATION

Description	Part#
Universal	204100

PADDED STRAP SHOULDER IMMOBILIZER

FEATURES AND BENEFITS

- Padded pressure sensitive straps improve comfort
- Waist strap is easily adjusted or can be removed
- Deep pocket design
- Two adjustment points for each strap allows for better fit
- Hook and loop closure make it easy to adjust



ORDERING INFORMATION

Size	Length	Height	U/M	Part#
X-Small	11" (28 cm)	7.5" (44 cm)	1 ea	204120
Small	13.5" (34 cm)	7.5" (44 cm)	1 ea	204122
Medium	15.5" (39 cm)	7.5" (44 cm)	1 ea	204123
Large	17.5" (44 cm)	7.5" (44 cm)	1 ea	204125
X-Large	20" (51 cm)	7.5" (44 cm)	1 ea	204127

PREMIUM CONTACT CLOSURE ARM SLING

INDICATIONS FOR USE

- Post upper limb trauma
- Post upper limb surgery

FEATURES AND BENEFITS

- Envelope style with thumb support
- Individually boxed
- Includes double "D" ring and slide buckle adjustment
- Canvas Material



ORDERING INFORMATION

Size	Length	Depth	Strap length	Standard	With Comfort Pad
X-Small	12" (30 cm)	4" (10 cm)	28" (71 cm)	204112	P204112
Small	14" (35 cm)	5" (13 cm)	30" (76 cm)	204113	P204113
Medium	17" (43 cm)	6.75" (17 cm)	35" (89 cm)	204115	P204115
Large	19" (48 cm)	6.75" (17 cm)	38" (97 cm)	204117	P204117
X-Large	19.5" (50 cm)	8" (20 cm)	43" (109 cm)	204118	P204118

BUCKLE CLOSURE ARM SLING

INDICATIONS FOR USE

- Post-trauma
- Post-surgery

FEATURES AND BENEFITS

- Contact closure shoulder strap with clip buckle swathe
- Available in 6 pack only



ORDERING INFORMATION

Size	Length	Depth	Strap length	Part#
Child	10.5" (27 cm)	5" (13 cm)	24" (61 cm)	204020
X-Small	12.25" (31 cm)	4" (10 cm)	32" (81 cm)	204022
Small	14.5" (37 cm)	5.25" (13 cm)	35" (89 cm)	204023
Medium	17.5" (45 cm)	6.5" (17 cm)	39.5" (100 cm)	204025
Large	19" (48 cm)	7.5" (19 cm)	43.5" (110 cm)	204027

BERGSCHULTZ UNIVERSAL ARM SLING

INDICATIONS FOR USE

- Post-trauma
- Post-operative

FEATURES AND BENEFITS

- Designed by Jim Bergschultz. R.O.T.
- Blue heavy denim material with contact closure

REQUIRED MEASUREMENTS

- Universal, one size



ORDERING INFORMATION

Description	Part#
Bergschulz Universal Arm Sling	204070

PREMIUM PADDED ARM SLING

INDICATIONS FOR USE

- Post-trauma
- Post-surgery

FEATURES AND BENEFITS

- Blue envelope style
- Canvas Material



ORDERING INFORMATION

Size	Length	Depth	Strap length	Part#
Small	13" (33 cm)	7" (18 cm)	31" (79 cm)	204313
Medium	16" (41 cm)	8" (20 cm)	35" (89 cm)	204315
Large	19" (48 cm)	9.5" (24 cm)	36" (91 cm)	204317
X-Large	22" (56 cm)	9.5" (24 cm)	40" (102 cm)	204318

FRONT CLOSURE CLAVICLE SUPPORT

INDICATIONS FOR USE


- Clavicle fractures and dislocations
- Postural support
- Soft tissue injuries

FEATURES AND BENEFITS

- Features front closure design for convenience in application
- Comfortable padding over clavicle area

REQUIRED MEASUREMENTS

- Chest circumference at the sternum area



ORDERING INFORMATION

Size	Circumference	Part#
X-Small	20 – 24" (51 – 61 cm)	205032
Small	24 – 30" (61 – 76 cm)	205033
Medium	30 – 36" (76 – 91 cm)	205035
Large	36 – 42" (91 – 107 cm)	205037
X-Large	42 – 48" (107 – 122 cm)	205038
XX-Large	48 – 54" (122 – 137 cm)	205039

FIGURE-8 CLAVICLE SPLINT

INDICATIONS FOR USE

- Clavicle fractures and dislocations
- Postural support
- Soft tissue injuries

FEATURES AND BENEFITS

- Prong buckle closure ensures correct fit without slippage

REQUIRED MEASUREMENTS

- Chest circumference



ORDERING INFORMATION

Size	Circumference	Part#
X-Small	20 – 24" (51 – 61 cm)	205002
Small	24 – 30" (61 – 76 cm)	205003
Medium	30 – 36" (76 – 91 cm)	205005
Large	36 – 42" (91 – 107 cm)	205007
X-Large	42 – 48" (107 – 122 cm)	205008

SLING AND SWATHE

INDICATIONS FOR USE

- Shoulder injuries
- Post-operative

FEATURES AND BENEFITS

- Comfortable position for patients having the arm next to the body
- One size fits all



ORDERING INFORMATION

Description	Part#
Sling and Swathe	13074

PEDIATRIC UNIVERSAL ARM SLING

INDICATIONS FOR USE

- Post upper limb trauma
- Post upper limb surgery

FEATURES AND BENEFITS

- Prong buckle closure correct fit without slippage

REQUIRED MEASUREMENTS

- Universal, one size



ORDERING INFORMATION

Description	Part#
Pediatric Universal Arm Sling Packed in boxes of 6	204020

TENNIS ELBOW SUPPORT WITH HOT/COLD GEL THERAPY

INDICATIONS FOR USE

- Epicondylitis

FEATURES AND BENEFITS

- Universal 2" / 5 cm wide with loop-lock Velcro closure
- Loop-lock contact closure
- Twin blue cushions provide counter force, reducing traction on the tendons and inflammation

REQUIRED MEASUREMENTS

- Circumference of the elbow



ORDERING INFORMATION

Description	Size	Part#
Tennis Elbow Support with hot/cold gel therapy	Universal (fits elbow circumference 8-14"/20 - 36 cm)	E310000

FORMFIT® TENNIS ELBOW

INDICATIONS FOR USE

- Epicondylitis

FEATURES AND BENEFITS

- Patented ergonomically-shaped plastic eliminates tourniquet effect associated with traditional designs.
- Twin blue cushions provide counterforce, reducing traction on the tendons and thus inflammation.
- Ribbed surfaces reduce vibration to the tendon, increasing comfort.
- Durable design and materials enable everyday use for both work and athletic activities.



REQUIRED MEASUREMENTS

- Circumference of the elbow

ORDERING INFORMATION

Size	Circumference	Part#
Small	8.7-10.4" (22-26 cm)	B-247000002
Medium	10.4-12.8" (27-32 cm)	B-247000003
Large	12.8-15" (33-38 cm)	B-247000004

AIRFORM® TENNIS ELBOW SUPPORT

INDICATIONS FOR USE

- Epicondylitis

FEATURES AND BENEFITS

- Focused pneumatic compression
- Conforms and comfortably cushions sensitive injured anatomy
- Easy to apply, simple design that fits most, reduces inventory needs



REQUIRED MEASUREMENTS

- Universal: Fits elbow circumference 20 – 36 cm / 8-14"

ORDERING INFORMATION

Description	Size	Part#
AirForm Tennis Elbow Support	Universal (fits elbow circumference 8-14" / 20 - 36 cm)	320000

FORMFIT® PRO WRIST

Össur Formfit Pro Wrist is constructed using our superior signature 3D MotionTech knit, and ideal for injuries that may benefit from compression. Malleable aluminum stays can be adjusted to conform to an individual's anatomy. CoolVent™ technology in the basal thumb area enhances breathability and comfort.



MotionTech™



CoolVent™



StayBilizer™



FoldLock™

INDICATIONS FOR USE

For wrist conditions that may benefit from compression, such as:

- Rheumatic Arthritis (RA) and other chronic pain conditions
- Wrist sprain and strains
- Overuse and pain conditions
- Wrist tendinitis

CONTRAINDICATIONS

- Must not be used by individuals for

whom compression is contraindicated.

FEATURES AND BENEFITS

- MotionTech™ for compression and dynamic fit has been proven to be superior to the leading brands at moisture wicking and breathability, keeping patients cool during wear.
- Malleable aluminum stay for a customized fit.
- CoolVent™ technology in thumb base area for extra breathability and comfort.

REQUIRED MEASUREMENTS



Measure wrist circumference at the level of the radial styloid with the hand outstretched. Wrist in 0 degrees of flexion.



Size	Wrist circumference		Blue		Black	
	cm	inch	Left Part#	Right Part#	Left Part#	Right Part#
XS	14 - 15	5½ - 5¾	FP401201	FP401101	FP400201	FP400101
S	15 - 16	5¾ - 6¼	FP401202	FP401102	FP400202	FP400102
M	16 - 17	6¼ - 6¾	FP401203	FP401103	FP400203	FP400103
L	17 - 18	6¾ - 7¼	FP401204	FP401104	FP400204	FP400104
XL	18 - 19	7¼ - 7½	FP401205	FP401105	FP400205	FP400105
XXL	19 - 20	7½ - 7¾	FP401206	FP401106	FP400206	FP400106

EXOFORM® WRIST

INDICATIONS FOR USE

- Mild to Moderate sprains and strains
- Tendonitis
- Post-cast healing
- Other select soft tissue injuries

FEATURES AND BENEFITS

- Ergonomically contoured padding
- Rigid plastic supports contour to the hand and wrist
- Adjustable aluminum palmar stay allows for fine-tuning when needed

REQUIRED MEASUREMENTS

- Wrist circumference



ORDERING INFORMATION

Size	Wrist Circumference	Left	Right
X-Small	5.25 – 6.25" (13 – 16 cm)	507082	507072
Small	6.25 – 7.25" (16 – 18.5 cm)	507083	507073
Medium	7.25 – 8" (18.5 – 20.5 cm)	507085	507075
Large	8 – 9" (20.5 – 23.5 cm)	507087	507077
X-Large	9 – 10" (23.5 – 25.5 cm)	507088	507078

EXOFORM® CARPAL TUNNEL WRIST

INDICATIONS FOR USE

- Carpal tunnel syndrome
- Tendonitis

FEATURES AND BENEFITS

- Allows for full MCP flexion
- Contoured for anatomical fit
- Provides radial and ulnar styloid relief
- Low profile
- Preshaped into a neutral position for comfortable limited motion/flexibility

REQUIRED MEASUREMENTS

- Wrist and palm circumference



ORDERING INFORMATION

Size	Wrist Circumference	Palm Circumference	Left	Right
Small	5.25 – 6.5" (13 – 16.5 cm)	6.75 – 8" (17 – 20.5 cm)	517083	517073
Medium	6.5 – 7.75" (16.5 – 19 cm)	8 – 9.25" (20.5 – 23.5 cm)	517085	517075
Large	7.75 – 10" (19 – 25.5 cm)	9.25 – 10.5" (23.5 – 26 cm)	517087	517077

FORMFIT® WRIST / FORMFIT® WRIST & FOREARM

Formfit Wrist, and Formfit Wrist & Forearm, provides immobilization and protection for patients suffering from tendonitis, select soft tissue injuries – including mild to moderate sprains and strains, as well as for post-cast healing. The lightweight, breathable material and contour fit make this comfortable product an excellent choice.

INDICATIONS FOR USE

- Mild to moderate sprains and strains
- Tendonitis
- Post-cast healing
- Other selected soft tissue injuries

FEATURES AND BENEFITS

- Breathable material
- Contour fit
- Contact Closure Straps
- Customizable Stays

REQUIRED MEASUREMENTS

- Wrist circumference



ORDERING INFORMATION

FORMFIT WRIST & FOREARM

Size	Wrist Circumference	8"		10"	
		Left	Right	Left	Right
XSmall	5.25 – 6.25" (13 – 16 cm)	317082	317072	617082	617072
Small	6.25 – 7.25" (16 – 18.5 cm)	317083	317073	617083	617073
Medium	7.25 – 8" (18.5 – 20.5 cm)	317085	317075	617085	617075
Large	8 – 9" (20.5 – 23.5 cm)	317087	317077	617087	617077
XLarge	9 – 10 (23.5 – 25.5 cm)	317088	317078	617088	617078

ACCESSORIES

GEL PAD FOR HOT/COLD THERAPY

Description	Left	Right
Formfit Wrist (8" version)*	GB20L	GB20R
Formfit Wrist & Forearm (10" version)	GB22L	GB22R



* Fits XSmall/Small. For Medium-Xlarge, please use Form Forearm 10" part numbers

6" WRIST

INDICATIONS FOR USE

- Mild to moderate sprains and strains
- Tendonitis
- Post-cast healing
- Other selected soft tissue injuries

FEATURES AND BENEFITS

- Contact closure straps
- Contoured fit
- Malleable stays

REQUIRED MEASUREMENTS

- Wrist circumference



ORDERING INFORMATION

Size	Wrist Circumference	Left	Right
Small	6.25 – 7.25" (16 – 18.5 cm)	417083C	417073C
Medium	7.25 – 8" (18.5 – 20.5 cm)	417085C	417075C
Large	8 – 9" (20.5 – 23.5 cm)	417087C	417077C
XLarge	9 – 10 (23.5 – 25.5 cm)	417088C	417078C

FORMFIT® WRIST & THUMB, UNIVERSAL

Formfit Universal Wrist and Thumb products are designed to restrict motion of the wrist and the thumb. Providing immobilization and protection from select soft tissue injuries, mild to moderate sprains and sprains, including gamekeeper's thumb and post cast healing. The lightweight, universal design also reduces inventory needs.



INDICATIONS FOR USE

- Mild to moderate sprains and strains
- Tendonitis
- Post-cast healing
- Other selected soft tissue injuries

FEATURES AND BENEFITS

- Quick-Pull close strap
- Positional dorsal stay
- Customizable thumb stay
- Malleable palmar stay

REQUIRED MEASUREMENTS

- Wrist circumference

ORDERING INFORMATION

FORMFIT WRIST

Size	Wrist Circumference	Left	Right
Pediatric	3.5 – 5" (9 – 13 cm)	B-252603301	B-252603300
8" Universal	4.5 – 9.5" (11 – 24 cm)	B-252603109	B-252503109
8" X-Large	7.5 – 11" (19 – 28 cm)	B-252603106	B-252603105

FORMFIT WRIST & FOREARM UNIVERSAL

Size	Wrist Circumference	Left	Right
10" Universal	4.5 – 9.5" (11 – 24 cm)	B-252603210	B-252603209

FORMFIT THUMB SPICA UNIVERSAL

Size	Wrist Circumference	Left	Right
Pediatric	3.5 – 5" (9 – 13 cm)	B-253603301	B-253603300
8" Universal	4.5 – 9.5" (11 – 24 cm)	B-253603109	B-253503109
8" X-Large	7.5 – 11" (19 – 28 cm)	B-253603106	B-253603105

FORMFIT® THUMB SPICA

The Formfit® Thumb Spica is designed to restrict wrist and thumb motion. It provides immobilization and protection for patients suffering from tendonitis (e.g. DeQuervain's Tendonitis), select soft tissue injuries – including mild to moderate sprains and strains and Gamekeeper's Thumb –, as well as for post-cast healing. The lightweight, breathable material and contour fit make this comfortable product an excellent choice.



INDICATIONS FOR USE

- Mild to moderate sprains and strains
- Tendonitis
 - e.g., DeQuervain's Tendonitis
- Post-cast healing
- Other selected soft tissue injuries
 - e.g., Gamekeeper's Thumb

FEATURES AND BENEFITS

- Breathable material
- Contour fit
- Contact Closure Straps
- Customizable Stays

REQUIRED MEASUREMENTS

- Wrist circumference

ORDERING INFORMATION

Size	Wrist Circumference	Left	Right
XSmall	5 – 6" (13 – 16 cm)	3020	3120
Small	6 – 7" (16 – 18.5 cm)	3030	3130
Medium	7 – 8" (18.5 – 20.5 cm)	3050	3150
Large	8 – 9" (20.5 – 23.5 cm)	3070	3170
XLarge	9 – 10" (23.5 – 25.5 cm)	3080	3180

ACCESSORIES

GEL PAD FOR HOT/COLD THERAPY

Description	Left	Right
Formfit Wrist (8" version)*	GB20L	GB20R
Formfit Wrist & Forearm (10" version)	GB22L	GB22R



* Fits XSmall/Small. For Medium-Xlarge, please use Form Forearm 10" part numbers

6" THUMB

INDICATIONS FOR USE

- Gamekeeper's thumb
- Torn collateral ligaments
- DeQuervain Tendonitis

FEATURES AND BENEFITS

- Contact closure straps
- Contoured fit
- Malleable stays

REQUIRED MEASUREMENTS

- Wrist circumference



ORDERING INFORMATION

Size	Wrist Circumference	Left	Right
Small	6.25 – 7.25" (16 – 18.5 cm)	31304C	30304C
Medium	7.25 – 8" (18.5 – 20.5 cm)	31504C	30504C
Large	8 – 9" (20.5 – 23.5 cm)	31704C	30704C
XLarge	9 – 10 (23.5 – 25.5 cm)	31804C	30804C

UNIVERSAL THUMB SPLINT

INDICATIONS FOR USE

- Gamekeeper's thumb
- Torn collateral ligaments
- DeQuervain Tendonitis

FEATURES AND BENEFITS

- Elastic contact closures
- Restricts movement of the MP/CMC
- Light weight and breathable

REQUIRED MEASUREMENTS

- Universal, one size



ORDERING INFORMATION

Description	Part#
Universal Thumb Splint	30000



HIP



REBOUND® HIP

Össur's Rebound Hip brace has a low-profile, lightweight design that supports patient comfort and compliance. It delivers gentle controlling forces and effective functional ROM restriction in extension/flexion and abduction/adduction to promote optimal post-surgical outcomes. The universal frame and softgood are easy to stock, assemble and quickly fit on patients.



INDICATIONS FOR USE

Hip conditions that could benefit from motion restriction following surgery:

- Femoral Acetabular Impingement arthroscopic surgery
- Gluteus medius repair surgery
- Hamstring repair surgery

FEATURES AND BENEFITS

- Low-profile, lightweight design that supports patient comfort and compliance
- Quick-release pegs easily adjust the desired range-of-motion needed for the patient's optimal rehab
- Quick-fit with sliding mechanism on frame to adjust for patient's height
- Optional, adjustable ROM hinge from -10°, 0° 10°, or 20°

REQUIRED MEASUREMENTS

- Circumference at iliac crest height
- Length from greater trochanter to mid patella

ORDERING INFORMATION

REBOUND HIP

Description	Part#
Rebound Hip	98013XRB
Rebound Hip with adjustable hinge	98012XRB
Rotation Strap	RHK-100024

X = 0 (Preset to left) or 5 (Preset to right)

*Frame: Universal; Softgood: Universal fitting to 48" hip circumference

ACCESSORIES

Description	Part#
OR Kit	980301
Belt Extender (Fitting 49-67" hip circumference)	980121RB
Fitting Tool	B-040100000



KNEE



REBOUND® CARTILAGE

The success of an articular cartilage repair is dependent on patient compliance to the rehabilitation process. The protocol is based on restricted cartilage load combined with gradual increase in activity. Many rehabilitation protocols suggest crutches for partial weight-bearing and it is the standard of care for successful outcomes. However, patient compliance to partial weight-bearing instructions with crutches has been shown to be inconsistent.

The Rebound Cartilage is a protective functional solution designed to support regeneration of the knee cartilage after repair procedures such as Microfracture, OATS and (M)ACI. Based on the clinically-proven 3 Points of Leverage System and patented Cartilage Protection Straps, joint unloading is maintained in flexion for protection of the cartilage during the healing process.

Our expertise in Injury Solutions has driven us to take the next step forward by developing indication-specific products which are designed to optimise healing while maintaining function and mobility.

FUNCTIONAL HEALING – WHERE MOBILITY MEANS RECOVERY

INDICATIONS FOR USE

- Articular cartilage defect repair, such as (M)ACI, OATS and Microfracture
- Cartilage knee trauma
- Meniscal repair
- Avascular necrosis
- Condylar bone marrow lesions (bone bruises)

FEATURES AND BENEFITS

- Dual Cartilage Protection Straps™ (CPS) maintain joint unloading in flexion for protection of the cartilage during the healing process.
- Breathable wrap around sleeve integrated with Sensil® silicone coated calf liners and doeskin thigh liners provide increased comfort while ensuring adequate suspension.
- Color coded Quick Fit buckles enable easy application and removal.
- Optional modular Flexion Control kit available when needed per prescribed rehab protocol for range of movement restriction.

REQUIRED MEASUREMENTS

- Circumference measurement: 6" (15 cm) below mid-patella

CUSTOM STANDARD COLORS



No charge for custom standard colors. Additional fees required for other custom colors, custom painted designs and custom printed graphics.

CUSTOM PAINTED DESIGNS



REBOUND® CARTILAGE

ORDERING INFORMATION



SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.

Description	Part#
SmartMeasure Measuring Kit	B-3 12000005
Constant Force Caliper	B-3 12000000

Scans, CMS (CTI only), CCS and casting also accepted. Knee M/L required

REBOUND CARTILAGE

Size	Circumference	Left medial	Right medial	Left lateral	Right lateral
Small	12.2 – 14" (31 – 35.5 cm)	B-254610002	B-254510002	B-254620002	B-254520002
Medium	14 – 16.75" (35.5 – 42.5 cm)	B-254610003	B-254510003	B-254620003	B-254520003
Large	16.75 – 20" (42.5 – 51 cm)	B-254610004	B-254510004	B-254620004	B-254520004
X-Large	20 – 24" (51 – 61.5 cm)	B-254610005	B-254510005	B-254620005	B-254520005

REBOUND CARTILAGE – WITH FLEXION CONTROL

Size	Circumference	Left medial	Right medial	Left lateral	Right lateral
Small	12.2 – 14" (31 – 35.5 cm)	B-254810002	B-254710002	B-254820002	B-254720002
Medium	14 – 16.75" (35.5 – 42.5 cm)	B-254810003	B-254710003	B-254820003	B-254720003
Large	16.75 – 20" (42.5 – 51 cm)	B-254810004	B-254710004	B-254820004	B-254720004
X-Large	20 – 24" (51 – 61.5 cm)	B-254810005	B-254710005	B-254820005	B-254720005

REBOUND CARTILAGE CUSTOM

Description	Part#
Rebound Cartilage Custom Right	B-154500001
Rebound Cartilage Custom Left	B-154600001

Note: Circumference limitation: Thigh circumferences smaller than 34.3 cm or larger than 73.5 cm. Calf circumferences smaller than 25 cm or larger than 61.5 cm.

REBOUND CARTILAGE CUSTOM – WITH FLEXION CONTROL

Description	Part#
Rebound Cartilage Custom – With flexion control Right	B-154700001
Rebound Cartilage Custom – With flexion control Left	B-154800001

REBOUND® CARTILAGE

ACCESSORIES

WRAPAROUND SLEEVE (WITH INTEGRATED SILICONE CALF LINER AND DOESKIN THIGH LINER)

Size	Right medial/Left lateral	Right lateral/Left medial
Small	B-254900072	B-254900082
Medium	B-254900073	B-254900083
Large	B-254900074	B-254900084
X-Large	B-254900075	B-254900085
Optional sleeveless liner kit	Right medial/Left lateral	Right lateral/Left medial
Silicone calf liner, doeskin thigh liner, strap pads	B-254900070	B-254900080
Doeskin calf and thigh liners, suspension strap kit, strap pads	B-427705000	B-427705001

STRAP KITS

Size	Complete strap kit (incl. 4 straps)	CPS strap kit (incl. 2 straps)
Small	B-701290002	B-700290002
Medium	B-701290003	B-700290003
Large	B-701290004	B-700290004
X-Large	B-701290005	B-700290005

UNIVERSAL FLEXION CONTROL KIT

Description	Part#
Left medial	B-705296100
Right medial	B-705295100
Left lateral	B-705296200
Right lateral	B-705295200

OTHER ACCESSORIES

Description	Part#
Velcro Twist-Fit Attachment Kit	B-705290000
Suspension Strap	UNL00213

REBOUND® PCL

The PCL, compared to the ACL, has the intrinsic ability to heal and regain continuity after an injury. The PCL experiences variable tensile forces during knee flexion. Insufficient healing after a PCL rupture leads to significant pathologies, such as, chronic anterior knee pain or knee OA. Healing of the PCL in an elongated position can lead to chronic instability and disability. Therefore, physiologically correct positioning of the PCL during conservative or post-surgical rehabilitation is important.

The custom made Rebound PCL is the world's first Dynamic Force PCL brace. It applies a physiologically correct dynamic force to the PCL during rehabilitation.

Our expertise in Injury Solutions has driven us to take the next step forward by developing indication-specific products which are designed to optimize healing while maintaining function and mobility.

FUNCTIONAL HEALING – WHERE MOBILITY MEANS RECOVERY

INDICATIONS FOR USE

- Non-surgical treatment of PCL ruptures
- Post-surgical rehabilitation for:
 - Partial PCL tears
 - PCL reconstruction / augmentation
 - PCL revision

FEATURES AND BENEFITS

- Dynamic Tension System (DTS) is designed to increase the load on the tibia as the knee goes through flexion, providing an anterior drawer force necessary to decrease load on the healing PCL to allow it to heal
- 3 options to customize the standardized, biochemically correct posterior tibial load for each patient by using shear knobs. And having the possibility to re-adjust the posterior tibial load, e.g. in case of swelling.
- Lightweight brace with proprietary Sensil® silicone calf liner and doeskin thigh liner reduce migration while providing comfort to the patient
- Custom molded posterior frame for optimized fit and support, with open anterior for ease of donning/doffing
- Available in 12 colors

REQUIRED MEASUREMENTS

- Custom-made to measure

STANDARD COLORS



REBOUND® PCL

ORDERING INFORMATION



SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.

Description	Part#
SmartMeasure Measuring Kit	B-3 12000005
Constant Force Caliper	B-3 12000000

Scans, CMS (CTI only), CCS and casting also accepted. Knee M/L required

REBOUND PCL

Description	Part#
Rebound PCL Left	B-124501155
Rebound PCL Right	B-124501154

ACCESSORIES

STRAPS

Size	Part#
X-Small – Large	PCL00061
X-Large – XX-Large	PCL00062

OTHER ACCESSORIES

Description	Part#
Extension Stop Kit	B-705 130042
Flexion Stop Kit	B-705 130043
Shear Knob Kit	B-465350087
Condyle Pad Kit	B-705 130046

CTI® CUSTOM

Introduced in 1983, CTi remains the cornerstone of our ligament bracing product line today. While the original CTi has evolved over time, our basic philosophy remains unchanged – a static support system consisting of a rigid frame will hold the bones in proper alignment and provide high level of control to prevent displacement. The constant development and improvement of this classic design is why the CTi remains the top choice of athletes and leading physicians around the world.

The rigid, carbon composite CTi knee brace frame captures the tibia, providing the brace-to-bone contact needed for maximum stabilization of the knee joint. The non-elastic cruciate strapping combined with a rigid frame that truly captures the tibia, results in a static support system.

The CTi Custom is generally used for high levels of support where comfort and strength are not compromised with weight and bulkiness. It can also offer greater advantages for accurate fitting where leg shape abnormalities present a problem.



INDICATIONS FOR USE

- For knee conditions that may benefit from increased AP and ML stability around the knee, such as: ACL, MCL, LCL, PCL, rotary and combined instabilities

FEATURES AND BENEFITS

- Anatomically-correct Accutrac® hinges with extension stops
- Proprietary Sensil® silicone liners designed to reduce migration
- Total Support System™ for ACL, MCL, LCL, PCL, rotary and combined instabilities
- Pre-preg carbon fiber, rigid frame construction
- Non-corrosive materials, ideal for water sports
- Custom paint option available
- Standard model for moderate to high impact levels
- Pro Sport model for the highest impact levels
- Vapor model for low to moderate impact levels
- Available in 12 colors

REQUIRED MEASUREMENTS

- Custom-made to measure

CUSTOM STANDARD COLORS



CUSTOM PAINTED DESIGNS



No charge for custom standard colors. Additional fees required for other custom colors, custom painted designs and custom printed graphics.

CTI® CUSTOM

ORDERING INFORMATION

CTI CUSTOM

Description	Left	Right
Standard	B-116600010	B-116500010
Pro Sport	B-117600010	B-117500010
Vapor	B-118600010	B-118500010



SMARTMEASURE®

A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.

Description	Part#
SmartMeasure Measuring Kit	B-312000005
Constant Force Caliper	B-312000000

Scans, CMS (CTI only), CCS and casting also accepted. Knee M/L required

ACCESSORIES

CALF LINERS

Description	Part#
CTi Custom liner	B-411500115

STRAPS KITS

Size	Part#
Small, Medium, Large	B-701130104
X-Large, XX-Large, XXX-Large, XXXX-Large	B-701130108

PCL KITS

Size	Part#
Small, Medium	B-705000132
Large, X-Large	B-705000235

CTI® CUSTOM

ACCESSORIES

STRAP PADS*

Size	Part#
Small	B-304010057
Medium, Large	B-304010058
X-Large, XX-Large	B-304010059

*Must also order hook-and-loop B-309020017

AMS (ANTI-MIGRATION SYSTEM) WRAPS

Description	Part#
AMS Wrap EVA	22095
AMS Wrap Neoprene SM-LG (up to 17" in circumference at 3" below knee)	22094
AMS Wrap Neoprene XL (above 17" in circumference at 3" below knee)	22152

MX GEAR

Description	Part#
Patella Cup	B-711130011
Gear Guards	B-705000003
MX Kit S (patella cup, gear guards, undersleeve, stickers)	B-711140002
MX Kit M (patella cup, gear guards, undersleeve, stickers)	B-711140003
MX Kit L (patella cup, gear guards, undersleeve, stickers)	B-711140004
MX Kit XL (patella cup, gear guards, undersleeve, stickers)	B-711140005

SLEEVES

Size	Circumference 15 cm above mid patella	Softsleeve	Neoprene Undersleeve	Neoprene Oversleeve	Sportsleeve	Padded Sports Oversleeve
Small	<15" (38 cm)	20643	20251	20051	20682	20448
Medium	15 – 18" (38 – 45 cm)	20644	20252	20052	20683	20448
Large	18 – 19.5" (45 – 49.5 cm)	20644	20253	20053	20684	20449
X-Large	19.5 – 21.5" (49.5 – 54.5 cm)	20645	20254	20054	20685	20449
XX-Large	21.5 – 24.5" (54.5 – 62 cm)	20645	20255	20055	N/A	
XXX-Large	24.5 – 28" (62 – 71 cm)	20646	20256	20056	N/A	
XXXX-Large	28 – 34" (71 – 86 cm)	20647	20257	20057	N/A	

OTHER ACCESSORIES

Description	Part#
CTi D-Ring Kit Standard 2x 1" and 2x 1.5" D-rings	B-704132600
CTi D-Ring Kit 2" 2x 2" D-rings	B-704132602
CTi Condyle Pad Kit	B-704130150
Add ACL Repair kit	XXXXX
Extension Stop Kit	22130
Flexion Stop Kit (Factory Installed)	19254
Extension strap Kit (Factory Installed)	21144

CTi®

The proven technology of the CTi knee brace is available with off-the-shelf convenience, providing adjustability and an excellent out-of-the-box fit while maintaining the level of support and protection that you expect from a CTi knee brace.

Outfitted with anatomically-correct Accutrac® hinges and flexible subshells and d-rings, the rock-solid CTi carbon frame provides a rigid exoskeleton that stabilizes the knee joint. Breathable liners coated with bio-inert Össur Sensil® Silicone ensure CTi OTS remains properly positioned on the leg, for optimal functional benefit.



INDICATIONS FOR USE

- For knee conditions that may benefit from increased AP and ML stability around the knee, such as: ACL, MCL, LCL, PCL, rotary and combined instabilities

FEATURES AND BENEFITS

- Anatomically-correct Accutrac hinges with extension stops to glide with the knee
- Flexible subshell and Adjustable Tibial Strap capture the tibia and allow adjustments in hinge depth for a more customized fit
- Breathable liners coated with bio-inert Össur Sensil Silicone to reduce migration and skin irritation
- Total Support System™ for ACL, MCL, LCL, PCL, rotary and combined instabilities
- Non-corrosive materials, ideal for water sports
- Available in 2 models: Standard model for moderate to high impact levels, Pro Sport model for the highest impact levels

REQUIRED MEASUREMENTS

- Caliper measurement: Medio-laterally at joint width in a standard, weight-bearing position
- Circumference: Just below the distal border of the patella



CTI®

ORDERING INFORMATION

CTI

Size	Caliper	Circumference	Standard Left	Standard Right	Pro Sport Left	Pro Sport Right
Small	8.9 – 10.2 cm	28 – 34 cm	B-238600112	B-238500112	B-239600112	B-239500112
Medium	10.2 - 11.4 cm	34 – 38 cm	B-238600113	B-238500113	B-239600113	B-239500113
Large	11.4 – 12.7 cm	38 – 41 cm	B-238600114	B-238500114	B-239600114	B-239500114
X-Large	12.7 – 14.0 cm	41 – 44 cm	B-238600115	B-238500115	B-239600115	B-239500115
XX-Large	14.0 – 15.2 cm	44 – 48 cm	B-238600116	B-238500116	B-239600116	B-239500116

Note: M/L measurements are the preferred and most accurate form of measurement for the CTI. Circumferential measurements are approximate.

ACCESSORIES

LINERS

Size	Silicone tibia cuff	Silicone calf	Silicone thigh	Lycra thigh
Small, Medium	B-304030030	B-425580002	B-425571002	B-425570002
Large		B-425580004	B-425571004	B-425570004
X-Large, XX-Large		B-425580005	B-425571004	B-425570004

SILICONE LINER KITS

Size	Description	Part#
Small, Medium	(tibia, calf, thigh and liners and strap pad)	B-703260012
Large	(tibia, calf, thigh and liners and strap pad)	B-703260014
X-Large, XX-Large, XXX-Large, XXXX-Large	(tibia, calf, thigh and liners and strap pad)	B-703260015

STRAP KITS*

Size	Part#
Small, Medium	B-701140103
Large, X-Large, XX-Large	B-701140104

*Must also order hook-and-loop B-309020017

PCL KITS

Size	Part#
Small, Medium	B-705000132
Large, X-Large	B-705000235

SUBSHELLS

Size	Description	Part#
Small	(thigh and calf subshells, screws and washers)	B-713260002
Medium	(thigh and calf subshells, screws and washers)	B-713260003
Large, X-Large	(thigh and calf subshells, screws and washers)	B-713260004
XX-Large	(thigh and calf subshells, screws and washers)	B-713260006
Subshell Screw		B-301060096
Subshell Washer		FR16000

CTI®

ACCESSORIES

AMS (ANTI-MIGRATION SYSTEM) WRAPS

Description	Part#
Silicone AMS Wrap Small, Medium	B-715260002
Silicone AMS Wrap Large, X-Large, XX-Large	B-715260004
Neoprene AMS Wrap Small, Medium	B-410000001
Neoprene AMS Wrap Large, X-Large	B-410000002

MX GEAR

Description	Part#
Patella Cup	B-711130011
Gear Guards	B-705000003
MX Kit S (patella cup, gear guards, undersleeve, stickers)	B-711140002
MX Kit M (patella cup, gear guards, undersleeve, stickers)	B-711140003
MX Kit L (patella cup, gear guards, undersleeve, stickers)	B-711140004
MX Kit XL (patella cup, gear guards, undersleeve, stickers)	B-711140005

SLEEVES

Circumference 6" (15 cm) above mid patella	Softsleeve	Neoprene Undersleeve	Neoprene Oversleeve	Sportsleeve	Padded Sports Oversleeve
<15" (38 cm)	20643	20251	20051	20682	20448
15 – 18" (38 – 45 cm)	20644	20252	20052	20683	20448
18 – 19.5" (45 – 49.5 cm)	20644	20253	20053	20684	20449
19.5 – 21.5" (49.5 – 54.5 cm)	20645	20254	20054	20685	20449
21.5 – 24.5" (54.5 – 62 cm)	20645	20255	20055	N/A	
24.5 – 28" (62 – 71 cm)	20646	20256	20056	N/A	
28 – 34" (71 – 86 cm)	20647	20257	20057	N/A	

STRAP PADS

Size	Part#
Small	B-304010057
Medium, Large	B-304010058
X-Large, XX-Large	B-304010059

OTHER ACCESSORIES

Description	Part#
CTi D-Ring Kit Standard 2x 1" and 2x 1.5" D-rings	B-704132600
CTi D-Ring Kit 2" 2x 2" D-rings	B-704132602
CTi Condyle Pad Kit	B-704130150
Extension Stop Kit	22130
Flexion Stop Kit	19254
Brace Bag	B-307020028

REBOUND® DUAL

The Rebound DUAL knee brace delivers functional support for ligament instabilities, including ligament instability with osteoarthritis, in a single, versatile aluminum brace. Suitable for low-to-high impact activities as well as sliding sports, Rebound DUAL is your new, go-to knee brace.



INDICATIONS FOR USE

- Supports ACL, MCL, LCL, PCL or combined instabilities
- Ligament instabilities with osteoarthritis


FEATURES AND BENEFITS

- Polycentric hinge for stability, durability, and low-profile design
- Adjustable tibia strap captures tibia for optimal stability
- Optional modular SmartDosing® kit available when needed for fine-tuned pain management and offloading
- Breathable ActiveGrip™ calf liner with anti-migration capabilities

REQUIRED MEASUREMENTS

- Caliper measurement: Medio-laterally at joint width in a standard, weight-bearing position
- Circumference: Just below the distal border of the patella

CUSTOM BRACE COLORS



CUSTOM LINER COLORS



REBOUND® DUAL

ORDERING INFORMATION

REBOUND DUAL

Size	Caliper	Circumference	Right	Left
XSmall	3-3.5" (7.6-8.9 cm)	15-16.5" (38.1-41.9 cm)	B-242510001	B-242610001
Small	3.5-4" (8.9-10.2 cm)	16.5-18.5" (41.9-47 cm)	B-242510002	B-242610002
Medium	4-4.5" (10.2-11.4 cm)	18.5-20.5" (47-52.1 cm)	B-242510003	B-242610003
Large	4.5-5" (11.4-12.7cm)	20.5"-22.5" (52.1-57.2 cm)	B-242510004	B-242610004
X-Large	5-5.5" (12.7-14cm)	22.5"-24.5" (57.2-62.2 cm)	B-242510005	B-242610005
XX-Large	5.5-6" (14-15.2cm)	24.5"-27.5" (62.2-69.9 cm)	B-242510006	B-242610006
XXX-Large	6-6.5" (15.2 – 16.5 cm)	27.5"-30" (69.9-76 cm)	B-242510007	B-242610007

M/L at joint width in a standard, weight-bearing position 2 Thigh Circumference: 6" above mid-patella

REBOUND DUAL CUSTOM

Measuring	Right	Left
SmartMeasure, DMS, Coordinate Cast System, Computer Aided Design	B-155510001	B-155610001

REBOUND DUAL ACCESSORIES

Description	Part #
PCL Strap Kit XS-L	B-701090057
PCL Strap Kit XL-XXXL	B-701090058
Field Serviceable D-Ring with Rivets Kit	B-701090059
Smart Dosing Kit RM/LL	B-701090060
Smart Dosing Kit RL/LM	B-701090061
Rebound DUAL Flexion Stop Kit - 45°, 60°, 70° and 90°	B-701090004
Rebound DUAL Extension Stop Kit - 5°, 10°, 20°, 30° and 40°	B-701090008

* Rebound DUAL is supplied with 0° and 10° extension stops. Additional stops are available to order.

SMARTMEASURE®



A digital measuring system for custom knee braces on the iPhone or iPad. The most intuitive, efficient and accurate way to measure for any Össur custom knee brace. Utilizes a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.

Description	Part#
SmartMeasure Measuring Kit	B-312000005
Constant Force Caliper	B-312000000

Scans, CMS (CTI only), CCS and casting also accepted. Knee M/L required

REBOUND® DUAL ST

Evolving from the functional support and versatility of the Rebound® DUAL, Össur offers a knee ligament brace for everyone. The Rebound DUAL ST, a sleek step-through design built with the strength of aircraft aluminum. Leveraging a light-weight profile and patented hi-tech ActiveGrip™ liners, the Rebound DUAL ST is designed for exceptional comfort and fitting with ease.



INDICATIONS FOR USE

For knee conditions that may benefit from increased AP and ML stability around the knee, such as: ACL, MCL, LCL, PCL, and combined instabilities

FEATURES AND BENEFITS

- Step-Through Design: advantage in comfort & fitting with ease
- Sleek, Strong Frame: offering the strength of aircraft aluminum with a light-weight design, the 14" (35 cm) frame length keeps a low-profile
- Anti-Gravity: Activegrip™ liners with the added advantage of Össur's anti-migration strap (AMS) system
- Easy for Everyone: simple strapping system and field-serviceable technology for adjustments on the go
- Versatile, comfortable, and reliable support for improved ligament stability
- Suitable for low-to-high impact activities as well as sliding sports
- Optional PCL kit accessory for support of the PCL
- Available in pre-fabricated sizes XS-3XL and also custom-made

REQUIRED MEASUREMENTS



Measure knee caliper medio-laterally at joint width in standing, weight-bearing position.



Measure thigh circumference 15 cm (6") above mid-patella. Leg in extension.

Size	Knee caliper		Thigh circumference 15 cm (6") above mid patella	
	cm	inch	cm	inch
XS	7.6 - 8.9	3 - 3½	38.1 - 41.9	15 - 16½
S	8.9 - 10.2	3½ - 4	41.9 - 47	16½ - 18½
M	10.2 - 11.4	4 - 4½	47 - 52.1	18½ - 20½
L	11.4 - 12.7	4½ - 5	52.1 - 57.2	20½ - 22½
XL	12.7 - 14	5 - 5½	57.2 - 62.2	22½ - 24½
XXL	14 - 15.2	5½ - 6	62.2 - 69.9	24½ - 27½
3XL	15.2 - 16.5	6 - 6½	69.9 - 76	27½ - 30

REBOUND® DUAL ST

ORDERING INFORMATION

PRE-FABRICATED

Size	Left	Right
XS	B-242620001	B-242520001
S	B-242620002	B-242520002
M	B-242620003	B-242520003
L	B-242620004	B-242520004
XL	B-242620005	B-242520005
XXL	B-242620006	B-242520006
3XL	B-242620007	B-242520007

ACCESSORIES

PCL STRAP KIT

Size	Part#
XS-L	B-701090065
XL-3XL	B-701090066

STRAP KIT

Size	Part#
XS-LG STRAP KIT	B-701090005
XL-XXXL STRAP KIT	B-701090006

CALF LINER

Size	Part#
XS	B-701090074
S	B-701090075
M	B-701090076
L	B-701090077
XL	B-701090078
XXL	B-701090079
3XL	B-701090080

THIGH LINER

Size	Part#
XS	B-701090067
S	B-701090068
M	B-701090069
L	B-701090070
XL	B-701090071
XXL	B-701090072
3XL	B-701090073

OTHER ACCESSORIES

Description	Part#
Extension stop kit	B-701090008
Flexion stop kit	B-701090004
Condyle pad kit	B-701090009
Field servicable D-ring kit	B-701090059

REBOUND® DUAL ST CUSTOM

Evolving from the functional support and versatility of the Rebound® DUAL, Össur offers a knee ligament brace for everyone. The Rebound DUAL ST, a sleek step-through design built with the strength of aircraft aluminum. Leveraging a light-weight profile and patented hi-tech ActiveGrip™ liners, the Rebound DUAL ST is designed for exceptional comfort and fitting with ease.



25 Color Options*

Customize to match your team colors



*Custom availability. Colors are approximate

INDICATIONS FOR USE

For knee conditions that may benefit from increased AP and ML stability around the knee, such as: ACL, MCL, LCL, PCL, rotary and combined instabilities

FEATURES AND BENEFITS

- Step-Through Design: advantage in comfort & fitting with ease
- Sleek, Strong Frame: offering the strength of aircraft aluminum with a light-weight design, the 14" (35 cm) frame length keeps a low-profile
- Anti-Gravity: Activegrip™ liners with the added advantage of Össur's anti-migration strap (AMS) system
- Easy for Everyone: simple strapping system and field-serviceable technology for adjustments on the go
- Versatile, comfortable, and reliable support for improved ligament stability
- Suitable for low-to-high impact activities as well as sliding sports
- Optional PCL kit accessory for support of the PCL
- Custom-made (also available in pre-fabricated sizes XS-3XL)

REQUIRED MEASUREMENTS



Measure calf circumference 15 cm (6") below mid-patella. Leg in extension.



Measure calf circumference 7.5 cm (3") below mid-patella. Leg in extension.



Measure knee caliper medio-laterally at joint width in standing, weight-bearing position.



Measure thigh circumference 15 cm (6") above mid-patella. Leg in extension.



Measure thigh circumference 7.5 cm (3") above mid-patella. Leg in extension.

REBOUND® DUAL ST CUSTOM

ORDERING INFORMATION

CUSTOM

Left	Right
B-155620002	B-155520001

ACCESSORIES

PCL STRAP KIT

Size	Part#
XS-L	B-701090065
XL-3XL	B-701090066

SPORT UNDERSLEEVE

Size	Part#
S	20682
M	20683
L	20684
XL	20685

NEOPRENE UNDERSLEEVE

Size	Part#
S	20251
M	20252
L	20253
XL	20254
XXL	20255
3XL	20256
4XL	20257

NEOPRENE OVERSLEEVE

Size	Part#
S	20051
M	20052
L	20053
XL	20054
XXL	20055
3XL	20056
4XL	20057

PADDED SPORTS OVERSLEEVE

Size	Part#
S	20448
L	20449

CALF LINER

Size	Part#
XS	B-701090074
S	B-701090075
M	B-701090076
L	B-701090077
XL	B-701090078
XXL	B-701090079
3XL	B-701090080

THIGH LINER

Size	Part#
XS	B-701090067
S	B-701090068
M	B-701090069
L	B-701090070
XL	B-701090071
XXL	B-701090072
3XL	B-701090073

STRAP KIT

Size	Part#
XS-L	B-701090005
XL-3XL	B-701090006

OTHER ACCESSORIES

Description	Part#
Extension stop kit	B-701090008
Flexion stop kit	B-701090004
Condyle pad kit	B-701090009
Field servicable D-ring kit	B-701090059

PARADIGM®

INDICATIONS FOR USE

- ACL, MCL, LCL, rotary and combined instabilities
- Protection and stabilization of surgical grafts
- Recommended for low to medium contact / impact

FEATURES AND BENEFITS

- Carbon fiber Triax® shells with flexible subshell allows for anatomical differences while providing firm control
- Hard-anodised, polycentric hinges
- Adjustable Tibial Strap captures the tibia and allows adjustments in hinge depth for a more customized fit
- Full ROM kit is available, enabling post surgical use
- Large, "open window" frame design for comfort and compliance



REQUIRED MEASUREMENTS

- Circumference: 15 cm above mid-patella

CUSTOM STANDARD COLORS



No charge for custom standard colors. Additional fees required for other custom colors.

ORDERING INFORMATION

Size	Circumference	Left	Right
X-Small	13.5 – 15.5" (34.3 – 39.4 cm)	B-210600001	B-210500001
Small	15.5 – 17.5" (39.4 – 44.5 cm)	B-210600002	B-210500002
Medium	17.5 – 20.5" (44.5 – 52.1 cm)	B-210600003	B-210500003
Large	20.5 – 23.5" (52.1 – 59.7 cm)	B-210600004	B-210500004
X-Large	23.5 – 26.5" (59.7 – 67.3 cm)	B-210600005	B-210500005
XX-Large	26.5 – 29.5" (67.3 – 74.9 cm)	B-210600006	B-210500006
Custom	N/A	B-110600020	B-110500020

ACCESSORIES

Size	Strap and Pad Kit	Liner Kit Left	Liner Kit Right
X-Small	B-701060001	B-703066001	B-703065001
Small	B-701060002	B-703066002	B-703065002
Medium	B-701060003	B-703066003	B-703065003
Large	B-701060004	B-703066004	B-703065004
X-Large	B-701060005	B-703066005	B-703065005
XX-Large	B-701060006	B-703066006	B-703065006

Description	Part#
Flexion / Extension Stop Kit	B-705000504

EXTREME®

INDICATIONS FOR USE

- Mild to severe ACL, MCL, LCL, PCL, rotary, combined instabilities, laxity or sprains
- Protection and stabilization of surgical grafts

FEATURES AND BENEFITS

- The proven single-upright design provides dynamic stability, offering maximal support when the knee is under stress
- The unique Dynamic Force Strap® minimizes bulk and maximizes agility
- Custom fabricated brace
- Available in standard and in SmartDosing® version. Tool-free SmartDosing dial option enables fine-tuning and can be adjusted even when wearing the brace under clothing.
- SmartDosing version available with Sensil® silicone calf liner to reduce migration
- Available in matte and glossy paint finishing and in 12 standard colors. Also optional is to have the brace custom painted.



REQUIRED MEASUREMENTS

- Custom-made to measure

ORDERING INFORMATION

Description	Part#
Extreme Left with DFS strap	B-108600010
Extreme Right with DFS strap	B-108500010
Extreme with SmartDosing Left	B-108600020
Extreme with SmartDosing Right	B-108500020

ACCESSORIES

Description	Part#
Patella cup	B-250110001
Hinge cover	B-251020000
AMS Wrap (custom made)	B-426000150
AMS Sleeve (custom made)	B-426000250

SLEEVES

Size	Össur knee undersleeve	Sofsleeve	Neoprene undersleeve	Neoprene oversleeve
X-Small	B-250150101			
Small	B-250150102	20643	20251	20051
Medium	B-250150103	20644	20252	20052
Large	B-250150104	20644	20253	20053
X-Large	B-250150105	20645	20254	20054
XX-Large	B-250150106	20645	20255	20055
XXX-Large		20646	20256	20056
XXXX-Large		20647	20257	20057

FLEX™

INDICATIONS FOR USE

- ACL, MCL, LCL, PCL, rotary and functional support
- Long-term use for post-op, rehab and functional support
- Recommended for low to moderate activity

FEATURES AND BENEFITS

- Accutrac® hinges with easy to use flexion/extension stops
- Heat-moldable, carbon composite frame accommodates edema and/or changes in leg shape during rehabilitation
- Gel Fit™ Condylar pad for improved suspension and comfort
- Total Support System™ for ACL, MCL, LCL, PCL, rotary and combined instabilities
- Adjustable Tibial Strap captured the tibia and allows adjustments in hinge depth for a customized fit
- Large, "open window" frame design accommodates post-op dressing and edema
- Non-corrosive materials, ideal for water sports



REQUIRED MEASUREMENTS

- Caliper measurement: Medial-lateral at joint width in a standard, weight-bearing position
- Circumference: Just below the distal border of the patella

ORDERING INFORMATION

Size	Caliper	Circumference	Standard		PCL	
			Left	Right	Left	Right
Small	3.5 – 4" (8.9 – 10.2 cm)	11 – 12.625" (27.9 – 32.1 cm)	FLEX-SL	FLEX-SR	FLEXPCL-SL	FLEXPCL-SR
Medium	4 – 4.625" (10.2 – 11.4 cm)	12.625 – 14.5" (32.1 – 36.2 cm)	FLEX-ML	FLEX-MR	FLEXPCL-ML	FLEXPCL-MR
Large	4.625 – 5.125" (11.4 – 12.4 cm)	14.5 – 16" (36.2 – 39.1 cm)	FLEX-LL	FLEX-LR	FLEXPCL-LL	FLEXPCL-LR
X-Large	5.125 – 5.625" (12.4 – 13.3 cm)	16 – 17.625" (39.1 – 42.2 cm)	FLEX-XLL	FLEX-XLR	FLEXPCL-XLL	FLEXPCL-XLR

ACCESSORIES

CALF & THIGH PADS

Size	Calf Pad Left	Calf Pad Right	Thigh Pad Left	Thigh Pad Right
Small	44048	44049	44046	44047
Medium	44040	44041	44042	44043
Large	44044	44045	44038	44039
X-Large	44052	44053	44050	44051

FLEX™

ACCESSORIES

STRAPS KITS

Size	Part#
Small	44071
Medium	44069
Large	44070
X-Large	44068

OTHER ACCESSORIES

Description	Part#
Flexion / Extension Stop Kit	44035
Hinge Cover Kit	44083
Condyle Pad Gel (1/4in)	22183
Condyle Pad Kit Neoprene (1/8in + 1/16in + cover)	44066
AMS (Anti-Migration System) Wrap S/M	23550
AMS (Anti-Migration System) Wrap L/XL	23551

FORMFIT® MCL

INDICATIONS FOR USE

- Mild to moderate Medial Collateral Ligament (MCL) sprains (Grade I and Grade II)

FEATURES AND BENEFITS

- Single, steel ROM lateral hinge (0° or 20° extension stop, 70° flexion stop)
- Dynamic Force Strap resist valgus forces
- Hybrid wrap/sleeve for easy application
- Adjustable tibial strap for improved suspension

REQUIRED MEASUREMENTS

- Circumference: 6" (15 cm) above mid-patella



ORDERING INFORMATION

Size	Circumference	Left	Right
X-Small	13 – 15.5" (34.3 – 39.4 cm)	B-237600001	B-237500001
Small	15.5 – 17.5" (39.4 – 44.5 cm)	B-237600002	B-237500002
Medium	17.5 – 20.5" (44.5 – 52.1 cm)	B-237600003	B-237500003
Large	20.5 – 23.5" (52.1 – 59.7 cm)	B-237600004	B-237500004
X-Large	23.5 – 26.5" (59.7 – 67.3 cm)	B-237600005	B-237500005

C180™ ROCKET

INDICATIONS FOR USE

- ACL, MCL, LCL, PCL, rotary and combined instabilities
- Recommended for low to medium contact/impact
- Recommended for pediatric usage

FEATURES AND BENEFITS

- Anatomically-correct Accutrac® hinges with easy-to-use extension stops
- Heat-Moldable frame for customization
- Adjustable Tibial Strap captures the tibia and allows adjustments in hinge depth for a more customized fit
- Available with an optional patella cup
- Total Support System™ for ACL, MCL, LCL, PCL, rotary and combined instabilities
- Semi-rigid frame consisting of woven carbon fiber infused with acrylic resin
- Lightweight, steamlined design for patient comfort and compliance
- Non-corrosive materials, ideal for water sports



REQUIRED MEASUREMENTS

- Caliper measurement: Medial-lateral at joint width in a standard, weight-bearing position
- Circumference: Just below the distal border of the patella

ORDERING INFORMATION

Size	Caliper	Circumference	Left	Right
Small	3 – 3.5" (7.6 – 8.9 cm)	9.5 – 11" (24.1 – 27.9 cm)	C180KSL	C180KSR
Large	3.3 – 3.75" (8.3 – 9.5 cm)	10 – 12" (26 – 29.8 cm)	C180KLL	C180KLR
X-Large	3.5 – 4" (8.9 – 10.2 cm)	11 – 12.5" (27.9 – 32.1 cm)	C1802XSL	C1802XSR

ACCESSORIES

Description	Part#	
Motocross kit Small / Large	10833	
Motocross kit X-Large	B-705130072	
Patella Cup Small / Large	25029	
Patella Cup X-Large	B-711130000	
Calf Liner Small	25061 (L)	26062 (R)
Calf Liner Large	25057 (L)	26058 (R)
Calf Liner X-Large	25106 (L)	25107 (R)
Thigh Liner Small	25060 (L)	26063 (R)
Thigh Liner Large	25056 (L)	26059 (R)
Thigh Liner X-Large	25104 (L)	25105 (R)

REBOUND® POST OP KNEE

The simple, versatile post-operative knee bracing solution that is easy for everyone. Whether adjusting ROM, setting the drop lock, or altering the brace length, you'll appreciate the simplicity of the Rebound Post-op knee brace.

INDICATIONS FOR USE

- ACL, PCL, MCL, LCL repairs/injuries
- Meniscal repairs
- Tibial plateau fractures
- Patellar tendon repairs
- Osteochondral repairs
- Condylar fractures
- Sprains/strains of the knee
- HTOs (High Tibial Osteotomy)

FEATURES AND BENEFITS

- Tool-free design and multi-purpose Smart-Fit™ packaging
- Available in universal, left or right contour and pediatric
- Telescoping frame from 18.5"-27" (Adult) and 15"-21" (Pediatric)
- OneTouch™ easy to adjust ROM hinge with lock-out capability
- Color-coding and numbered frame for patient education
- Quick-fit, gripping buckles for easy application

REQUIRED MEASUREMENTS

- Circumference: 6" (15 cm) above mid-patella



ORDERING INFORMATION

Description	Thigh circumference	Cool	Full Foam
Pediatric	12 – 25" (30.5 cm – 63.5 cm)	B-8381470	
Universal, Left contoured	16 – 32" (41 cm – 81 cm)	B-8381680	B-8382680
Universal, Right contoured	16 – 32" (41 cm – 81 cm)	B-8381580	B-8382580
Universal	16 – 32" (41 cm – 81 cm)	B-8381480	B-8382480
Universal XL	16 – 38" (41 cm – 96.5 cm)	B-8381490	B-8382490

ACCESSORIES

Description	Part#
ROM stop locks	B-8380005
Alligator hook and loop kit	B-8380006
Adhesive hook and loop kit pediatric	B-8380008
Adhesive hook and loop kit universal / universal XL	B-8380007

REPLACEMENT SOFTGOODS

Description	Part#
ROM stop locks	B-8380005
Alligator hook and loop kit	B-8380006
Adhesive hook and loop kit pediatric	B-8380008
Adhesive hook and loop kit universal / universal XL	B-8380007

FORMFIT® POST-OP KNEE

Formfit Post-Op Knee provides range of motion control and immobilization following injury or surgery. Length adjustment is tool free and easy. The low profile ROM hinge can be adjusted with a simple pull, rotate and release motion. Wide range of flexion and extension settings allow you to choose the optimal range of motion.

FEATURES AND BENEFITS

- Breakaway struts for easy, tool free length adjustment
- OneTouch™ easy to adjust ROM hinge with lock-out capability from -10° to 30° (10° incremental)
- Extension control from -10° to 90° (10° incremental)
- Flexion control from 0° to 120° (10° incremental)
- Available in 3 models, 3 lengths, and cool or full foam



INDICATIONS FOR USE

- Post-op/post-injury ROM control or immobilization. Indicated for ligament and fracture-related surgeries and when controlled range of motion is necessary

REQUIRED MEASUREMENTS

- Thigh circumference 15 cm (6") above mid patella

ORDERING INFORMATION

Size	Thigh circumference	Cool	Full Foam
Knee PO 24in DL	UP TO 27" (69cm)	FPOK403100	FPOK403200
Knee PO XL 24in DL	27" – 35" (69-89cm)	FPOK403102	FPOK403202
Knee PO 28in DL	UP TO 27" (69cm)	FPOK403101	FPOK403201
Knee PO XL 28in DL	27" – 35" (69-89cm)	FPOK403103	POK403203
Knee PO 16in DL	UP TO 27" (69cm)	-	FPOK403300
Knee PO XL 16in DL	27" – 35" (69-89cm)	-	FPOK403301

SOFT GOODS

Description	Cool	Full Foam
Replacement Softgood, Regular	403175	403275
Replacement Softgood, XL	403176	403277
Replacement Softgood, Regular Long	403276	403278

REPLACEMENT SOFTGOODS

Size	
ROM stop locks (set of 20 lock clips)	B-8380005
Strap Reduction Kit Cool	403199
Strap Reduction Kit Full Foam	403299

EXOFORM[®] KNEE IMMOBILIZER

INDICATIONS FOR USE

- Post-injury
- Post-trauma

FEATURES AND BENEFITS

- Slide-to-size adjustment enables you to provide the right fit for the largest variety of patient lengths and widths
- provides superior fit and limits migration
- Cool foam version can be worn directly against the skin for greater comfort on warm days
- Full foam version provides warmth and compression

REQUIRED MEASUREMENTS

- Universal – fits leg inseams of 18 – 24" (46-61 cm)



ORDERING INFORMATION

Description	Leg circumference	Part#
Cool	up to 28" (71.1 cm)	222100
Full foam	up to 28" (71.1 cm)	222000
Full foam XXXL	up to 38" (96.5 cm)	222050

ACCESSORIES

Description	Part#
Replacement softgood, Cool	222101
Replacement softgood, Full Foam	222001
Replacement softgood, Full Foam, 3XL	222051

FORMFIT® KNEE IMMOBILIZER

INDICATIONS FOR USE

- Emergency immobilization
- Pre and post-surgical knee immobilization

FEATURES AND BENEFITS

- Adjustable medial/lateral stays and rigid posterior stay provide stability
- Soft breathable foam keeps patients cool
- Universal 3-piece design simplifies inventory needs

REQUIRED MEASUREMENTS

- Universal



ORDERING INFORMATION

Length	Part#
12"	200009
14"	200010
16"	200016
18"	200118
20"	200020
22"	200022
24"	200030

ACCESSORIES

Description	Part#
Hot/Cold Therapy Gel Pad	GB7

FX™ PATELLA STABILIZER

INDICATIONS FOR USE

- Ideal for Patellofemoral Pain Syndrome
- Patellofemoral Stress Syndrome (Runner's Knee)
- Patellar Subluxation/Dislocation
- Chondromalacia
- Patellar Tendinopathy (Tendonitis or Jumper's Knee)
- Arthritis and Osgood-Schlatter disease.

FEATURES AND BENEFITS

- Offers superior patellar tracking through a combination of maximum coronal plane control and minimal flexion and extension resistance
- The detachable thigh strap provides a customizable fit. Additionally, the coronal adjuster allows for further specific fit
- The patented stabilizer offers superior control and assists in restoring normal kinematics to the joint
- Breathable Coolflex sleeve minimizes heat and perspiration



REQUIRED MEASUREMENTS

- Knee circumference at mid-patella

ORDERING INFORMATION

Size	Knee circumference	Part#
X-Small	10.5 - 12"	B-218000001
Small	12 - 13.5"	B-218000002
Medium	13.5 - 15"	B-218000003
Large	15 - 17"	B-218000004
X-Large	17 - 19"	B-218000005
XX-Large	19 - 21"	B-218000006
XXX-Large	21 - 23"	B-218000007

ACCESSORIES

Size	Knee circumference	Part#
XSmall	Osgood-Schlatter Buttress	B-304040037
Medium	Osgood-Schlatter Buttress	B-304040038
XLarge	Osgood-Schlatter Buttress	B-304040039

FORMFIT® HINGED LATERAL J

INDICATIONS FOR USE

- Lateral patellar subluxation and dislocations
- Lateral patellofemoral malalignment
- Patellofemoral maltracking and pain
- Chondromalacia due to lateral maltracking

FEATURES AND BENEFITS

- Breathable material
- Removable hinges
- Easy pull-on
- Lateral J buttress
- Open popliteal



REQUIRED MEASUREMENTS

- Thigh Circumference: 6" above mid-patella
- Calf Circumference: 6" below mid-patella

ORDERING INFORMATION

Size	Thigh circumference	Knee circumference	Calf circumference	Right leg	Left leg
X-Small	13 - 15.5" (33 - 39.4cm)	12-13" (30.5-33cm)	10-12" (25.4-30.5cm)	709152	709162
Small	15.5-18.5" (39.4-47cm)	13-14" (33-35.6cm)	12-14" (30.5-35.6cm)	709154	709164
Medium	18.5-21" (47-53.3cm)	14-15" (35.6-38.1cm)	14-16" (35.6-40.6cm)	709155	709165
Large	21-23.5" (53.3-59.7cm)	15-17" (38.1-43.2cm)	16-18" (40.6-45.7cm)	709157	709167
X-Large	23.5-26.5" (59.7-67.3cm)	17-19" (43.2-48.3cm)	18-20" (45.7-50.8cm)	709158	709168
XX-Large	26.5-29.5" (67.3-74.9cm)	19-21" (48.3-53.3cm)	20-22" (50.8-55.9cm)	709159	709169

FORMFIT® KNEE, HINGED

The Formfit® Knee brace is for patients suffering from mild to moderate strains, sprains or ligament tears (ACL/PCL/MCL/LCL). It provides excellent stability, protection and compression during injury recovery. Made from strong and breathable materials it offers comfort while remaining highly durable. The open popliteal area reduces material bulging and potential irritation, while the malleable rigid hinge arms provide an intimate fit.

The Formfit Knee brace comes in either a wrap or a sleeve version.

INDICATIONS FOR USE

- Mild to moderate ACL and PCL tears
- MCL and LCL tears
- Combined mild knee sprains and strains
- Instabilities of the knee

FEATURES AND BENEFITS

- Breathable material
- Easy pull-on
- Malleable hinges
- Universal buttress
- Open popliteal

REQUIRED MEASUREMENTS

- Measure thigh circumference 6" above mid-patella



ORDERING INFORMATION

Size	Knee Circumference	Sleeve	Wrap
XSmall	13" – 15.5" (33 – 39 cm)	504252	503252
Small	15.5 – 18.5" (39 – 47 cm)	504253	503253
Medium	18.5 – 21" (47 – 53 cm)	504255	503255
Large	21 – 23.5" (53 – 60 cm)	504257	503257
XLarge	23.5 – 26.5" (60 – 67 cm)	504258	503258
XXLarge	26.5 – 29.5" (67 – 75 cm)	504259	503259
XXXLarge	29.5 – 32" (75 – 81 cm)		503260

Accessories

Description	Part#
Formfit Knee Brace Buttress (Sleeve only)	503000

FORMFIT® KNEE ROM

The Formfit® Knee brace is for patients suffering from mild to moderate strains, sprains or ligament tears (ACL/PCL/MCL/LCL). It provides excellent stability, protection and compression during injury recovery. Made from strong and breathable materials it offers comfort while remaining highly durable. The open popliteal area reduces material bulging and potential irritation, while the malleable rigid hinge arms provide an intimate fit.



The Formfit Knee brace comes in either a wrap or a sleeve version.

INDICATIONS FOR USE

- Mild to moderate ACL and PCL tears
- MCL and LCL tears
- Combined mild knee sprains and strains
- Instabilities of the knee

FEATURES AND BENEFITS

- Breathable material
- Easy pull-on
- Malleable hinges
- Universal buttress
- Open popliteal

REQUIRED MEASUREMENTS

- Measure thigh circumference 6" above mid-patella

ORDERING INFORMATION

Size	Knee Circumference	Sleeve		Wrap	
		Long ROM Sleeve	Short ROM Sleeve	Long ROM Wrap	Short ROM Wrap
XSmall	13" – 15.5" (33 – 39 cm)	508352	508252	507352	507252
Small	15.5 – 18.5" (39 – 47 cm)	508353	508253	507353	507253
Medium	18.5 – 21" (47 – 53 cm)	508355	508255	507355	507255
Large	21 – 23.5" (53 – 60 cm)	508357	508257	507357	507257
XLarge	23.5 – 26.5" (60 – 67 cm)	508358	508258	507358	507258
XXLarge	26.5 – 29.5" (67 – 75 cm)	508359	508259	507359	507259
XXXLarge	29.5 – 32" (75 – 81 cm)			507360	507260

REBOUND® KNEE BRACE

With six variations to choose from, Össur's Rebound hinged knee brace product line is a comprehensive portfolio of high-quality braces designed to serve the needs of a broad range of patients.

INDICATIONS FOR USE

- Mild to moderate knee sprains and strains
- Mild to moderate Anterior Cruciate Ligament (ACL) and Posterior Cruciate Ligament (PCL) tears
- Medial Collateral Ligament (MCL) and Lateral Collateral Ligament (LCL) tears

FEATURES AND BENEFITS

- Featuring our unique, 3-layer Cooltech™ fabric, the Rebound Knee Brace is highly breathable, reduces heat retention and dries ultrafast.
- Low profile hinges with extension and flexion stop kit options (ROM version) and malleable arms to provide a custom fit and knee movement control.
- Extension Stops: 0, 5, 10 (pre-installed), 20, 30, 40
- Flexion stops: 45, 60, 75, 90
- Contoured pre-tensed inset fabric that reduces irritation and material bulging in the popliteal area as well as inventory
- Removable straps allow application from right or left side and strap pull tab allow for consistent donning by patient
- Universal buttresses are included with all sleeve versions for customizable patella support



REQUIRED MEASUREMENTS

- Thigh Circumference: 6" (15 cm) above mid-patella

ORDERING INFORMATION

REBOUND ROM KNEE

Size	Circumference	Sleeve Long	Sleeve Short	Wrap Long	Wrap Short
X-Small	11.5 – 13.75" (29 – 35 cm)	708152	708052	707152	707052
Small	13.75 – 16" (35 – 41 cm)	708153	708053	707153	707053
Medium	16 – 18" (41 – 46 cm)	708155	708055	707155	707055
Large	18 – 20.5" (46 – 52 cm)	708157	708057	707157	707057
X-Large	20.5 – 22.5" (52 – 57 cm)	708158	708058	707158	707058
XX-Large	22.5 – 24.75" (57 – 63 cm)	708159	708059	707159	707059
XXX-Large	24.75 – 28.5" (63 – 72 cm)			707160	707060
XXXX-Large	28.5 – 32" (72 – 81 cm)			707161	707061

REBOUND KNEE

Size	Circumference	Sleeve Short	Wrap Short
X-Small	11.5 – 13.75" (29 – 35 cm)	704052	703052
Small	13.75 – 16" (35 – 41 cm)	704053	703053
Medium	16 – 18" (41 – 46 cm)	704055	703055
Large	18 – 20.5" (46 – 52 cm)	704057	703057
X-Large	20.5 – 22.5" (52 – 57 cm)	704058	703058
XX-Large	22.5 – 24.75" (57 – 63 cm)	704059	703059
XXX-Large	24.75 – 28.5" (63 – 72 cm)	704060	703060
XXXX-Large	28.5 – 32" (72 – 81 cm)	704061	

NEOPRENE HINGED KNEE SUPPORT

INDICATIONS FOR USE

- Moderate ligament instabilities

FEATURES AND BENEFITS

- Accordion molding in the popliteal area reduces bunching behind the knee
- Contoured fit for additional comfort and protection
- The hinges are molded into the brace to maintain a low profile

REQUIRED MEASUREMENTS

- Knee circumference at mid-patella



ORDERING INFORMATION

Size	Knee circumference	Part#
Small	11.5 – 13" (29 – 33 cm)	302153BLK
Medium	13 – 14.5" (33 – 37 cm)	302155BLK
Large	14.5 – 16" (37 – 41 cm)	302157BLK
X-Large	16 – 17.5" (41 – 44 cm)	302158BLK
XX-Large	17.5 – 19" (44 – 48 cm)	302159BLK

NEOPRENE WRAPAROUND HINGED KNEE SUPPORT

INDICATIONS FOR USE

- For patients with unstable or injured knees that may benefit from additional knee support, such as post-operation rehabilitation, tendonitis, bursitis, arthritis; and sprains and strains.

FEATURES AND BENEFITS

- Dual action hinges incorporate our patented three dimensional moulding technology
- The hinges are molded into the brace to maintain a low profile

REQUIRED MEASUREMENTS

- Knee circumference at mid-patella



ORDERING INFORMATION

Size	Knee circumference	Part#
Small	13 – 14"	212003
Medium	14 – 15"	212005
Large	15 – 16"	212007
X-Large	16 – 17"	212008
XX-Large	17 – 18"	212009
3X-Large	18 – 23"	212010
4X-Large	19 – 25"	212011

NEOPRENE KNEE SUPPORT WITH STABILIZED PATELLA

INDICATIONS FOR USE

- For patients suffering from knee conditions that may benefit from additional knee and patella support such as patellar tendonitis, subluxation dislocation, post-operation rehabilitation, post arthroscopic surgery, condromalacia, bursitis, arthritis and tendonitis.

FEATURES AND BENEFITS

- Molded hardware to maintain a low profile
- Accordion molding in the popliteal area eliminates bunching
- Soft cushioned liner to offer unequaled comfort
- Varied areas of thickness give additional protection
- Lycra® liner offers a breathable contact area for comfortable and extensive use



REQUIRED MEASUREMENTS

- Knee circumference at mid-patella

ORDERING INFORMATION

Size	Knee circumference	Part#
Small	11.5-13"	302563BLK
Medium	13-14.5"	302565BLK
Large	14.5-16"	302567BLK
X-Large	16-17.5"	302568BLK
XX-Large	17.5-19"	302569BLK

FORMFIT® PRO KNEE

Össur Formfit® Pro Knee is a compressive knee sleeve that features a SmartPad™ buttress for targeted support and cushioning. Our patented design realigns the patella, suitable for conditions such as Patella Tendinitis (Jumper's Knee) and Chondromalacia Patella. The dynamic, 3-D MotionTech™ knit is supportive, meets certified medical grade compression requirements and has been proven to be superior to the leading brands at moisture wicking and breathability, keeping your knee cool during wear. Lightweight, removable spiral stays ensure optimum shape and fit, while distinct patella graphics make donning intuitive.



MotionTech™



CoolVent™



StayBilizer™



SmartPad™

INDICATIONS FOR USE

For knee conditions that may benefit from compression, such as:

- Patella tendinitis
- Chondromalacia Patella
- Mild sprains and strains
- Knee pain
- Feeling of instability


CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated.

FEATURES AND BENEFITS

- MotionTech™ for compression and dynamic fit with superior moisture wicking and breathability, keeping your knee cool during wear.
- Certified medical grade compression (20-30 mmHg).
- Lightweight CoolVent™ knit in popliteal area for added breathability and comfort.
- Patented SmartPad™ for targeted support and cushioning.
- Patella graphics for intuitive placement.
- Patented FoldLock™ system for easier and more secure application.
- Lightweight, removable spiral steel stays.

REQUIRED MEASUREMENTS



Measure calf circumference
15 cm (6") below mid-patella.
Leg in extension.



Measure thigh circumference
15 cm (6") above mid-patella.
Leg in extension.



Size	Calf circumference 15 cm (6") below mid patella		Thigh circumference 15 cm (6") above mid patella		Blue	Black
	cm	inch	cm	inch	Part#	Part#
XS	26 - 31	10¼ - 12¼	36 - 41	14⅞ - 16⅞	FP101001	FP100001
S	31 - 36	12¼ - 14⅞	41 - 46	16⅞ - 18⅞	FP101002	FP100002
M	36 - 41	14⅞ - 16⅞	46 - 51	18⅞ - 20⅞	FP101003	FP100003
L	41 - 46	16⅞ - 18⅞	51 - 56	20⅞ - 22	FP101004	FP100004
XL	46 - 50	18⅞ - 19¾	56 - 60	22 - 23⅝	FP101005	FP100005

NEOPRENE KNEE SLEEVE

INDICATIONS FOR USE

- For patients suffering from weak and/or injured knees that may benefit from additional knee support.

FEATURES AND BENEFITS

- Unique accordion design enhances comfort by eliminating bunching in the popliteal area
- Open or closed patella

REQUIRED MEASUREMENTS

- Knee circumference at mid-patella



ORDERING INFORMATION

		0.3 cm Neoprene		0.6 cm Neoprene	
Size	Circumference	Open Patella	Closed Patella	Open Patella	Closed Patella
Small	11.5 – 13" (29 – 33 cm)	402013BLK	412013BLK	302013BLK	312013BLK
Medium	13 – 14.5" (33 – 37 cm)	402015BLK	412015BLK	302015BLK	312015BLK
Large	14.5 – 16" (37 – 41 cm)	402017BLK	412017BLK	302017BLK	312017BLK
X-Large	16 – 17.5" (41 – 44 cm)	402018BLK	412018BLK	302018BLK	312018BLK
XX-Large	17.5 – 19" (44 – 48 cm)	402019BLK	412019BLK	302019BLK	312019BLK

PREMIUM ELASTIC KNEE SUPPORT

INDICATIONS FOR USE

- Knee conditions that may benefit from an elastic support, such as swelling, feeling of instability, and/or stability around the knee cap.

FEATURES AND BENEFITS

- Cool, contoured, comfortable fit
- Made with high-quality, long-lasting, breathable materials
- Flexible support and compression
- Reduces swelling from minor injuries
- Adds stability to the knee cap

REQUIRED MEASUREMENTS

- Knee circumference at mid-patella



ORDERING INFORMATION

Size	Knee circumference	Part#
Small	13 – 15" (33 – 38 cm)	129603
Medium	15 – 18" (38 – 46 cm)	129595
Large	18 – 21" (46 – 53 cm)	129587
X-Large	21 – 24" (53 – 61 cm)	129578

ELASTIC KNEE SUPPORT

INDICATIONS FOR USE

- Knee conditions that may benefit from an elastic support, such as swelling and feeling of instability.

FEATURES AND BENEFITS

- Cool, contoured, comfortable fit
- Made with high-quality, long-lasting, breathable materials
- Flexible support and compression
- Reduces swelling from minor injuries

REQUIRED MEASUREMENTS

- Knee circumference at mid-patella



ORDERING INFORMATION

Size	Knee circumference	Part#
Small	13 – 15" (33 – 38 cm)	13350
Medium	15 – 18" (38 – 46 cm)	13351
Large	18 – 21" (46 – 53 cm)	13352
X-Large	21 – 24" (53 – 61 cm)	13353

FOOT AND ANKLE



REBOUND® DIABETIC WALKER

Offloading has proven to be a key factor in foot ulcer treatment and entails reducing pressure of the plantar area in general and the ulcer site in particular. The Rebound Diabetic Walker is uniquely designed for effective treatment of plantar foot ulcers, with easy access and re-fitting during wound care, while improving the patient's mobility during the healing process.

The Rebound Diabetic Walker has major advantages over diabetic shoes and shoe inserts in its healing effectiveness and over total contact casting due to its ease of use, speed of application and low cost.

Our expertise in Injury Solutions has driven us to take the next step forward by developing indication-specific products which are designed to optimise healing while maintaining function and mobility.



FUNCTIONAL HEALING – WHERE MOBILITY MEANS RECOVERY

INDICATIONS FOR USE

- Plantar ulcers needing offloading

FEATURES AND BENEFITS

- Effective tri-laminate insole with removable shock-absorbing hexagons that promote ulcer offloading and a Plastazote, heat-moldable top layer that evenly distributes the pressure and minimizes friction
- Unique pressure relief valve for preventing over-inflation, while balancing compression with customized fit
- Positioning Portals allow the clinician to check the fit of the walker at the heel bottom, in order to ensure appropriate positioning of the foot in the walker
- Anatomically designed shell that immobilizes the foot and ankle at 90 degrees while providing protection
- Flex Edge® technology for relief of pressure points
- Rocker Bottom outsole for a smooth roll over while allowing for individual natural toe out
- Soft and breathable, low friction, machine washable liner for hygiene promotion and prevention of infection
- Tamper evident strapping for compliance monitoring

REQUIRED MEASUREMENTS

- Shoe size



REBOUND® DIABETIC WALKER

ORDERING INFORMATION

REBOUND DIABETIC WALKER

Size	Men's Shoe Size US	Women's Shoe Size US	Shoe size UK	Shoe size EU	Part#
Small	4.5-7	5.5-8.5	4 – 6.5	36 – 39	B-242500002
Medium	7-10.5	8.5-12	7 – 10	39.5 – 44	B-242500003
Large	10.5-12.5	12-14	10.5 – 12	44.5 – 46.5	B-242500004
X-Large	12.5-15	14-16.5	12.5 – 14.5	47 – 50	B-242500005

ACCESSORIES

REPLACEMENT SOFTGOODS

Size	Liners	Insoles
Small	B-418500102	B-418500202
Medium	B-418500103	B-418500203
Large	B-418500104	B-418500204
X-Large	B-418500105	B-418500205

WALKER COVERS

Size	Hygiene covers (2 pcs)	Durable weather cover	Disposable weather cover (20 pcs)
Small	B-418500602	B-418500502	B-418500512
Medium	B-418500603	B-418500503	
Large	B-418500604	B-418500504	B-418500514
X-Large	B-418500605	B-418500505	

DIABETIC SOCKS (2 PCS)

Size	Part#
Small, Medium	B-418500703
Large, X-Large	B-418500704

COMPLIANCE STRAPS (5 PCS)

Size	Part#
Universal	B-418500300



HYGIENE COVER



DURABLE WEATHER COVER



DISPOSABLE WEATHER COVER



ACHILLES WEDGE KIT



INSOLE

REBOUND® AIR WALKER

The Rebound Air Walker features a host of details designed to provide added protection, superior compression, simplified application and usage, improved comfort and, ultimately, successful patient outcomes.



INDICATIONS FOR USE

REBOUND AIR WALKER

Foot and ankle conditions that may benefit from motion restriction, such as:

- Grade 2 and 3 ankle sprains
- Stable fractures
- Tendon or ligamentous sprain
- During postoperative or rehab phase. If used with the Achilles Wedge Kit, Achilles ruptures can be treated as well.

REBOUND AIR WALKER LOW TOP

Foot and ankle conditions that may benefit from motion restriction, such as:

- Stable foot fractures
- Tendon or ligamentous sprain
- During postoperative or rehab phase.

FEATURES AND BENEFITS

- Full circumferential shell with additional semi-soft toe bumper for frontal toe protection
- Soft Flex Edge™ reduces pressure points
- Full rocker bottom is especially engineered to be low and wide to help promote a natural gait
- Stretch-to-fit shell for optimal fit
- Shell can be trimmed for additional stretch in the calf
- Cushioning and breathable softgood provides comfort
- Linked medial/lateral bladders for easy inflation and equal compression
- Integrated pump and air release valve for easy access

REQUIRED MEASUREMENTS

- Shoe size



REBOUND® AIR WALKER

ORDERING INFORMATION

REBOUND AIR WALKER

Size	Men's Shoe Size US	Women's Shoe Size US	Standard	Low Top
X-Small	2 – 4	3.5 – 5.5	B-242900001	B-242900061
Small	4 – 7	5.5 – 8.5	B-242900002	B-242900062
Medium	7 – 10.5	8.5 – 12	B-242900003	B-242900063
Large	10.5 – 12.5	12 – 14	B-242900004	B-242900064
X-Large	12.5 – 15	14 – 16.5	B-242900005	B-242900065

ACCESSORIES

REPLACEMENT SOFTGOODS

Size	Slim Liner / High Top	Slim Liner / Low Top
X-Small	B-418030101	B-418030161
Small	B-418030102	B-418030162
Medium	B-418030103	B-418030163
Large	B-418030104	B-418030164
X-Large	B-418030105	B-418030165

REBOUND ACHILLES WEDGE KIT

Size	Part#
Small	AWK0001
Medium	AWK0005
Large	AWK0009
Universal	AWK00015

WALKER COVERS

Size	Hygiene covers (2 pcs)	Durable weather cover	Disposable weather cover (20 pcs)
Small	B-418500602	B-418500502	B-418500512
Medium	B-418500603	B-418500503	
Large	B-418500604	B-418500504	B-418500514
X-Large	B-418500605	B-418500505	

OTHER ACCESSORIES

Description	Part#
Adhesive Insole for Wedge	WHW-W0107
Strap Extension	B-502100101
Rebound Hightop Socks SM/MD	B-418500703
Rebound Hightop Sock LG/XL	B-418500704

AFO DYNAMIC

The AFO Dynamic provides a strong, lightweight solution for patients with varying degrees of drop foot secondary to a variety of neurological conditions. Össur has taken the carbon fiber technology from its Flex-Foot® product and incorporated it into the new AFO Dynamic to provide increased support over traditional ankle foot orthoses. The energy storing properties of carbon fiber make it an ideal material for this dynamic ankle foot orthosis because it provides real energy return. Best for those people leading busy lives, the AFO Dynamic will ensure a more efficient walk for the user.



INDICATIONS FOR USE

Drop foot secondary to:

- Cerebro Vascular Accident (CVA)
- Mild knee instability
- Other neurological conditions
- Charcot Marie Tooth disease (HMSN)
- Multiple Sclerosis

Contraindications:


- Severe ankle – foot deformities
- Fluctuating edema
- Severe spasticity

FEATURES AND BENEFITS

- Carbon fiber design provides excellent strength-to-weight ratio and dynamic motion
- The tapered carbon heel stores energy and manages forces at heel strike
- Full-length toe lever provides full support while foot is loaded, as well as reduced shock to the contralateral limb
- Creates medial-lateral stability at mid-stance
- Reduces energy consumption by reducing hip hike
- Medial strut makes it less visible in shoe
- Open anterior tibia shell provides superior anatomical fit, and reduces heat and perspiration
- Weight limit: 120kg (265lbs)

REQUIRED MEASUREMENTS

- Shoe size
- Calf height at thickest part of the calf



AFO DYNAMIC

ORDERING INFORMATION

AFO DYNAMIC

Size	Calf Height	Shoe Size US		Shoe Size UK		Shoe Size EU cm		Left	Right
		Male	Female	Male	Female	Male	Female		
X-Small	30 – 34 cm	3 – 5	4 – 6.5	2.5 – 4.5	< 4	34 – 37	< 37	O-102021	O-102011
Small	32 – 36 cm	5.5 – 7.5	7 – 9	5 – 7	4.5 – 6.5	37.5 – 40	37.5 – 40	O-102022	O-102012
Medium	34 – 38 cm	8 – 10.5	9.5 – 12	7.5 – 10	7 – 9.5	40.5 – 44.5	40.5 – 44	O-102024	O-102014
Large	36 – 40 cm	11 – 14		10.5 – 13.5		45 – 48.5		O-102026	O-102016

ACCESSORIES

ACCESSORY KITS

Size	Part#
X-Small	O-102101
Small	O-102102
Medium	O-102104
Large	O-102106

Note: Accessory kits come with foam padding, upper and lower Velcro straps, and Velcro dots.

AFO LEAF SPRING

AFO Leaf Spring is a pre-fabricated polypropylene ankle-foot orthosis designed to support flaccid drop foot. It provides a semirigid section for toe clearance and support, yet has a thin, flexible foot plate for easy trimming. The absence of a heel section makes the Leaf Spring more comfortable to wear and provides a better fit in shoes.

INDICATIONS FOR USE

- Drop foot
- Mild drop foot secondary to other neurological pathologies

FEATURES AND BENEFITS

- Lightweight injection-molded polypropylene
- Thin, flexible foot part which is easy to trim with scissors
- Absence of heel section for improved comfort



REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

Size	Shoe Size US	Shoe Size UK	Shoe Size EU	Height	Left	Right	Foot Plate Length
Small – Medium	7.5 – 10.5	3.5 – 6.5	36 – 39	10.75" (27 cm)	O-110227	O-110127	24 cm
Large	11 – 14	7 – 10.5	39.5 – 45	12" (30 cm)	O-110230	O-110130	28 cm

AFO LIGHT

AFO Light provides support only where it is needed, eliminating unnecessary bulk and providing a cosmetically-appealing look. Using the best of Flex-Foot® technology and the energy storing properties of carbon fiber, AFO Light offers a strong and lightweight solution for people with mild to moderate Drop Foot.



INDICATIONS FOR USE

- Mild to moderate drop foot
- Contraindications:
 - Severe ankle – foot deformities
 - Medial-lateral instability
 - Severe spasticity
 - Extreme activity
 - Fluctuating edema

FEATURES AND BENEFITS

- Energy return and propulsion assistance of the limb at toe off
- Supports the foot for ground clearance during swing phase
- Carbon fiber design provides excellent strength-to-weight ratio and dynamic motion
- Gradual plantar flexion at initial contact
- Stiffness gradient at heel and the toe correspond to ground reaction forces as they are applied to the foot
- Easy to don and doff
- Weight limit: 100kg (220lbs)

REQUIRED MEASUREMENTS

- Shoe size
- Calf height at thickest part of the calf

ORDERING INFORMATION

AFO LIGHT

Size	Calf height	Shoe Size US		Shoe Size UK		Shoe Size EU		Left	Right
		Male	Female	Male	Female	Male	Female		
Small	32 – 36 cm	Up to 7.5	7 – 9	< 7	4.5 – 6.5	< 40	37.5 – 40	O-121022	O-121012
Medium	34 – 38 cm	8 – 10.5	9.5 – 12	7.5 – 10	7 – 9.5	40.5 – 44.5	40.5 – 44	O-121024	O-121014
Large	36 – 40 cm	11 – 14		10.5 – 13.5		45 – 48.5		O-121026	O-121016
X-Large	38 – 42 cm	12 – 15		11.5 – 14.5		46 – 50		O-121028	O-121018

ACCESSORIES

Description	Part#
Replacement Kit	O-122101

REBOUND® FOOT-UP®

Foot-Up is a lightweight ankle orthosis that offers dynamic support for Drop Foot or similar complaints. The revolutionary orthosis provides visible improvement in gait by supporting the foot the moment it is raised.



INDICATIONS FOR USE

- Drop foot
- Support of dorsiflexion

FEATURES AND BENEFITS

- Streamlined
- Lightweight
- Invisible when worn
- Comfortable to wear
- Breathable material
- Padded
- Easy to fit
- Separate inlays available
- Available in black and beige



REQUIRED MEASUREMENTS

- Basic version: Circumference of ankle, measured around the malleoli
- Shoeless version: Circumference of foot immediately behind the ball of the foot

ORDERING INFORMATION

FOOT-UP

Size	Ankle Circumference	Beige	Black
Medium	7 – 8.5" (18 - 21 cm)	07811 -1	07810-1
Large	8.5 – 10" (22 - 26 cm)	07811 -2	07810-2
X-Large	10 – 13" (27 - 33 cm)	07811 -3	07810-3

SHOELESS VERSION*

Size	Foot Circumference	Beige	Black
Small	< 8" (20 cm)	27811-1	27810-1
Medium	8 – 9.5" (20 - 24 cm)	27811-2	27810-2
Large	9.5 – 11" (24 - 28 cm)	27811-3	27810-3
X-Large	> 11" (28 cm)	27811-4	27810-4

*Shoeless version item comprises of foot wrap only and must be ordered with the Foot-Up unless the patient already has one.

ACCESSORIES

Description	Beige	Black
Foot-Up plastic insert	17811	17810

FORMFIT® WALKER, AIR

INDICATIONS FOR USE

- Soft tissue injuries
- Grade 2 and 3 sprains
- Stable fractures
- Post-operative, trauma or rehab use

FEATURES AND BENEFITS

- High and low top versions available
- Integrated pneumatic system
- Low, wide rocker bottom promotes a stable, natural gait
- Shock absorbing sole for reduced impact at heel strike
- Breathable, seamless liner



REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

Size	Men's Shoe Size US	Women's Shoe Size US	Shoe Size UK	Shoe Size EU	Standard	Low Top
X-Small	2 – 4	3.5 – 5.5	1.5 – 3.5	33 – 35.5	A-W0200BLK	A-W0300BLK
Small	4.5 – 7	6 – 8	4 – 6.5	36 – 39	A-W0400BLK	A-W0500BLK
Medium	7.5 – 10.5	8.5 – 11.5	7 – 10	39.5 – 44	A-W0600BLK	A-W0700BLK
Large	10.5 – 12.5	11.5 – 13.5	10.5 – 12	44.5 – 46.5	A-W0800BLK	A-W0900BLK
X-Large	12.5 +	13.5 +	12.5 – 14.5	47 – 50	A-W1000BLK	A-W1100BLK

ACCESSORIES

REPLACEMENT SOFTGOODS

Size	Standard	Low Top
X-Small	12349	12356
Small	10675	12357
Medium	10678	12358
Large	10681	12359
X-Large	12347	12360

FORMFIT® WALKER

INDICATIONS FOR USE

- Soft tissue injuries
- Grade 2 and 3 sprains
- Stable fractures
- Post-operative, trauma or rehab use

FEATURES AND BENEFITS

- High and low top versions available
- Shock absorbing sole reduces impact during heel strike
- Low, wide rocker bottom promotes a stable, natural gait
- Breathable, seamless liner



REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

Size	Men's Shoe Size US	Women's Shoe Size US	Shoe Size UK	Shoe Size EU	Standard	Low Top
X-Small	2 – 4	3.5 – 5.5	1.5 – 3.5	33 – 35.5	W0200BLK	W0300BLK
Small	4.5 – 7	6 – 8	4 – 6.5	36 – 39	W0400BLK	W0500BLK
Medium	7.5 – 10.5	8.5 – 11.5	7 – 10	39.5 – 44	W0600BLK	W0700BLK
Large	10.5 – 12.5	11.5 – 13.5	10.5 – 12	44.5 – 46.5	W0800BLK	W0900BLK
X-Large	12.5 +	13.5 +	12.5 – 14.5	47 – 50	W1000BLK	W1100BLK

ACCESSORIES

REPLACEMENT SOFTGOODS

Size	Standard	Low Top
X-Small	WKSG200BLK	WKSG0300BLK
Small	WKSG400BLK	WKSG500BLK
Medium	WKSG600BLK	WKSG700BLK
Large	WKSG800BLK	WKSG900BLK
X-Large	WKSG1000BLK	WKSG1100BLK

OTHER ACCESSORIES

Description	Part#
10 cm softgood width extension	WKSG-1BLK
Tibia and calf plates for use when additional splinting is required	WKSP100
Hot/Cold therapy insert for safe and reusable cryotherapy	WKGEL100
Heel wedge, Small	AWK0001
Heel wedge, Medium	AWK0005
Heel wedge, Large	AWK0009
Adhesive insole for heel wedge	WHW-W0107
Walker Strap Extender	10882

ROM WALKER

INDICATIONS FOR USE

- Foot and ankle conditions that may benefit from motion restriction or controlled range of motion, such as:
- Grade 2 and 3 ankle sprains
- Stable fractures
- Achilles ruptures
- During postoperative or rehab phase

FEATURES AND BENEFITS

- Lightweight and low profile aluminum side struts
- Rocker bottom design assists gait
- Decreased heel-strike force on calcaneus
- Simple to apply and adjust by circumferential straps
- Dial lock hinges with adjustable ROM limitation
- ROM control from 0-45° of dorsi/plantar flexion (7.5° interval)
- Available in Standard and Air version



ORDERING INFORMATION



Shoe size.

Size	Shoe size				Standard	Air
	Men's Shoe Size US	Women's Shoe Size US	Shoe Size UK	Shoe Size EU	Part#	Part#
XS	up to 4	up to 5.5	up to 3.5	up to 35.5	RY-WK002-XS	RY-WK003-XS
S	4.5 - 7	6 - 8.5	4 - 6.5	36 - 39	RY-WK002-SM	RY-WK003-SM
M	7.3 - 9	9 - 10.5	7 - 9	39.5 - 43	RY-WK002-MD	RY-WK003-MD
L	9.5 - 11.5	10.5 - 13.5	9.5 - 10.5	43.5 - 45.5	RY-WK002-LG	RY-WK003-LG
XL	12 and up	14 and up	11 and up	46 and up	RY-WK002-XL	RY-WK003-XL

DH OFFLOADING WALKER®

INDICATIONS FOR USE

- Treatment of plantar ulcers

FEATURES AND BENEFITS

- Shear force reduction by stabilizing foot and ankle at 90°
- Wide low profile rocker bottom provides a smooth roll-through and increased stability for altered gaits
- The softgood's patented inner bootie prevents foot from migrating on the insole
- Designed around Össur's unique Active Offloading Insole
- Provides a safe environment for the simultaneous use of wound care material
- Easy to apply and modify

REQUIRED MEASUREMENTS

- Shoe size



ORDERING INFORMATION

DH OFFLOADING WALKER

Size	Shoe Size US (Men)	Shoe Size US (Women)	Shoe Size EU	Part#
X-Small	2 – 4	3.5 – 5.5	33 – 35.5	DH0200BLK
Small	4.5 – 7	6 – 8	36 – 39	DH0400BLK
Medium	7.5 – 10.5	8.5 – 11.5	39.5 – 44	DH0600BLK
Large	10.5 – 12.5	11.5 – 13.5	44.5 – 46.5	DH0800BLK
X-Large	12.5+	13.5+	47 – 50	DH1000BLK

ACCESSORIES

REPLACEMENT PARTS

Size	Softgoods	Insoles
X-Small	Call for info	11742C
Small	10958	PP04C
Medium	10959	PP06C
Large	10960	11740C
X-Large	Call for info	11741C

DH OFFLOADING POST-OP STYLE SHOE

INDICATIONS FOR USE

- Plantar wounds

FEATURES AND BENEFITS

- Designed to accommodate either the left or the right foot
- Designed around Össur's patented Active Offloading Insole
- Designed to accommodate bulky dressings
- Super-soft nylon upper with reinforced heel counter provides extra comfort and stability
- Patented removable toe piece aides in wound accessibility for daily wound site inspection

REQUIRED MEASUREMENTS

- Shoe size



ORDERING INFORMATION

DH POST-OP SHOE

Size	Shoe Size US (Men)	Shoe Size US (Women)	Shoe Size EU	Part#
Small	4 – 6	6.5 – 8.5	36 – 39	10340
Medium	6.5 – 8.5	9 – 11	39.5 – 44	10341
Large	9 – 11	12+	44.5 – 46.5	10342
X-Large	11.5 - 14	13+	47 – 50	10343

ACCESSORIES

REPLACEMENT PARTS

Size	Insoles
Small	12365
Medium	12366
Large	12367
X-Large	12368

PEDIATRIC WALKER

INDICATIONS FOR USE

Foot and ankle conditions that may benefit from motion restriction, such as:

- Grade 2 and 3 ankle sprains
- Stable fractures
- During postoperative or rehab phase.

FEATURES AND BENEFITS

- High quality rigid outer skeleton for protection
- Dense foam lining within the sole to provide comfort and support
- Shock absorbing sole to improve comfort and protects the user from slips, providing additional security in each step
- Uprights give additional stability by restricting motion of the leg
- Lightweight design helps in keeping kids active

REQUIRED MEASUREMENTS

- Shoe size



ORDERING INFORMATION

PEDIATRIC WALKER

Size	Shoe Size US	Part#
Small	Child 6.5 – 9.5	PEW0400
Medium	Child 10 – Youth 1	PEW0600
Large	Youth 1.5 – 2.5	PEW0800

REPLACEMENT LINERS

Size	Shoe Size US	Part#
Small	Child 6.5 – 9.5	PEW0001
Medium	Child 10 – Youth	PEW0002
Large	Youth 1.5 – 2.5	PEW0003

METAL UPRIGHT AIR WALKER

INDICATIONS FOR USE

Foot and ankle conditions that may benefit from motion restriction, such as:

- Grade 2 and 3 ankle sprains
- Stable fractures
- During postoperative or rehab phase

FEATURES AND BENEFITS

- Malleable struts allow you to customized fit
- Rocker bottom sole promotes a more natural gait
- Pneumatic pump allows air compression to be adjusted as needed for a more customized fit.



REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

Size	Men's shoe size US	Women's shoe size US	Part#
X-Small	up to 4	up to 5.5	MAW0200
Small	4.5 – 7	6 – 8	MAW0400
Medium	7.5 – 10.5	8.5 – 11.5	MAW0600
Large	10.5 – 12.5	11.5 – 13.5	MAW0800
X-Large	12.5 & up	13.5 & up	MAW1000

FORMFIT PRO ANKLE

Össur Formfit® Pro Ankle is a compressive ankle sleeve with 2 integrated malleoli pads, designed to worn after mild ankle injuries. The dynamic 3-D MotionTech™ knit is supportive and meets certified medical grade compression requirements. Proven superior to leading brands for moisture wicking and breathability, your ankle stays cool during wear. CoolVent™ technology in the area in front of the ankle further enhances breathability and comfort.



MotionTech™



CoolVent™



SmartPad™



FoldLock™

INDICATIONS FOR USE

For Ankle conditions that may benefit from compression, such as:

- Ankle injury rehabilitation
- Mild sprains and strains
- Mild instability


CONTRAINDICATIONS

- Must not be used by individuals for whom compression is contraindicated

FEATURES AND BENEFITS

- MotionTech™ for compression and dynamic fit with superior moisture wicking and breathability, keeping your ankle cool during wear
- Certified medical grade compression (20-30 mmHg)
- CoolVent™ technology in the area in front of the ankle for added breathability and comfort
- 2 integrated malleoli SmartPads™ for targeted support and cushioning
- Patented FoldLock™ system for easier and more secure application

REQUIRED MEASUREMENTS



Measure ankle circumference diagonally over the base of the heel and across the ankle bone. Ankle in 90 degrees.



BLUE

Size	Ankle circumference		Blue		Black	
	cm	inch	Left	Right	Left	Right
XS	18 - 20	7 1/8 - 7 7/8	FP201201	FP201101	FP200201	FP200101
S	20 - 23	7 7/8 - 9	FP201202	FP201102	FP200202	FP200102
M	23 - 26	9 - 10 1/4	FP201203	FP201103	FP200203	FP200103
L	26 - 29	10 1/4 - 11 3/8	FP201204	FP201104	FP200204	FP200104
XL	29 - 32	11 3/8 - 12 5/8	FP201205	FP201105	FP200205	FP200105

REBOUND® ANKLE BRACE

INDICATIONS FOR USE

- Grade 1, 2, 3 and high ankle sprains
- Transition out of a walker and/or cast after fracture or surgery
- Chronic ankle instability
- Prophylactic use

FEATURES AND BENEFITS

- Adjustable calf band accommodates a broad range of leg circumferences and creates a more custom fit
- Open front, crossover strap design allows simple and secure application
- Version with stability strap offers numerous strapping configurations for a range of stability levels. Can be stepped down or removed as the patient progresses through rehabilitation



REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

Size	Men's Shoe Size US	Women's Shoe size US	Shoe Size UK	Shoe Size EU	Left	Right
Small	6 – 7.5	7.5 – 9	2 – 7.5	35 – 39.5	B-23 1602512	B-23 1502512
Medium	8 – 12	9.5 – 13.5	8 – 10.5	39.5 – 42	B-23 1602513	B-23 1502513
Large	12.5 – 16	N/A	11+	42+	B-23 1602514	B-23 1502514

ACCESSORIES

Description	Part#
Width Extender Kit	B-23 1900000
Stability Strap Kit	B-701 160080

FORMFIT® ANKLE BRACE

The Formfit® Ankle brace provides stability and protection for the ankle – whether used prophylactically, for chronic instability or following an injury. A hybrid between a rigid ankle stirrup and a soft ankle support, the Formfit Ankle ensures a close fit which limits inversion, eversion, flexion and extension movements. The low profile design fits easily into shoes and the quick lace up feature makes it easy to don and doff. Breathable material ensures it can be comfortably worn all day.

The Formfit Ankle with Figure 8 strap version, offers an even greater level of support, by further limiting the range of motion of the ankle.



INDICATIONS FOR USE

- Grade I-II sprains and strains
- Chronic instability
- Prophylactic use

FEATURES AND BENEFITS

- Breathable material
- Quick lace Up
- Inversion/Eversion Control
- Available with Figure 8 strap



REQUIRED MEASUREMENTS

- Measure ankle circumference diagonally over the base of the heel and across the ankle bone.

ORDERING INFORMATION

Size	Ankle Circ.*	Formfit Ankle Brace	Formfit Ankle with Figure 8
X-Small	11 – 12" (28 – 30.5 cm)	B-212010001	B-212000001
Small	12 – 13" (30.5 – 33 cm)	B-212010002	B-212000002
Medium	13 – 14" (33 – 35.5 cm)	B-212010003	B-212000003
Large	14 – 15" (35.5 – 38 cm)	B-212010004	B-212000004
X-Large	15 – 16" (38 – 40.5 cm)	B-212010005	B-212000005

FORMFIT® ANKLE WITH SPEEDLACE

The Formfit® Ankle Brace with Speedlace provides stability and protection for the ankle – whether used prophylactically, for chronic instability or following an injury. Made from durable materials, its low profile, quick-lace design is comfortable and fits easily into shoes. Removable semi-rigid stays provide additional support and the Figure 8 straps are designed for inversion/eversion control.

INDICATIONS FOR USE

- Grade I-II sprains and strains
- Chronic instability
- Prophylactic use

FEATURES AND BENEFITS

- Durable materials
- Quick-Pull closure straps
- Figure 8 straps
- Removable stays
- Circumferential cuff

REQUIRED MEASUREMENTS

- Measure ankle circumference diagonally over the base of the heel and across the ankle bone.



ORDERING INFORMATION

Size	Ankle Circ.*	Part#
XXSmall	9 – 10" (23 cm – 25.5 cm)	W-10621
XSmall	10 – 11" (25.5 cm – 28 cm)	W-10622
Small	11 – 12" (28 cm – 30.5 cm)	W-10623
Medium	12 – 13" (30.5 cm – 33 cm)	W-10625
Large	13 – 14" (33 cm – 35.5 cm)	W-10627
XLarge	14 – 15" (35.5 cm – 38 cm)	W-10628
XXLarge	15 – 16" (38 cm – 40.5 cm)	W-10629

*Measure ankle circumference diagonally over the base of the heel and across the ankle bone.

AIRFORM[®] UNIVERSAL INFLATABLE ANKLE STIRRUP

INDICATIONS FOR USE

- Acute Ankle Injuries
- Post-cast healing and post-walker support
- Chronic Instability

FEATURES AND BENEFITS

- Universal design fits both left and right ankles, reducing inventory needs
- Soft Flex-Edge molding provides a comfort barrier between the shell and ankle
- Inflatable air bladders provides adjustable pneumatic compression
- Adjustable heel strap allow for easy width adjustments



ORDERING INFORMATION

Size	White Single	White 6 pack	Black Single
Youth/Pony	8025 1	8025 16	
Adult	80250	802506	80350

AIRFORM® INFLATABLE ANKLE STIRRUP

INDICATIONS FOR USE

- Ankle injury rehabilitation/ immobilization for grades I & II strains and sprains
- Chronic Instability
- Post Cast healing

FEATURES AND BENEFITS

- Inflatable bladder provides adjustable pneumatic compression
- Side specific right or left sizing



ORDERING INFORMATION

Size	Left	Right
Youth/Pony	70200	70201
Adult	80200	80201

AIRFORM® PRE-INFLATED ANKLE STIRRUP

INDICATIONS FOR USE

- Ankle injury rehabilitation/ immobilization for grades I & II strains and sprains
- Chronic Instability
- Post Cast healing

FEATURES AND BENEFITS

- Pre-inflated, pulsating air bladder with regulated chamber reduces bottoming out and decreases edema
- No straws or over inflation
- Adjustable pivoting heel strap increases comfort and allows for full range of motion
- Side specific right or left sizing



ORDERING INFORMATION

Size	Left	Right
Youth/Pony	70100	70101
Adult	80100	80101

FORMFIT® HONEYCOMB ANKLE STIRRUP

INDICATIONS FOR USE

- Ankle injury rehabilitation/ immobilization for grades I & II strains and sprains
- Chronic Instability
- Post Cast healing

FEATURES AND BENEFITS

- Patented Honeycomb pads
- Flex-Edges and adjustable pivoting heel straps



ORDERING INFORMATION

Size	Left	Right
Youth/Pony	40100	40101
Adult	30100	30101

FORMFIT® FOAM ANKLE BRACE

INDICATIONS FOR USE

- Ankle injury rehabilitation/ immobilization for grades I & II strains and sprains
- Chronic Instability
- Post Cast healing

FEATURES AND BENEFITS

- Patented Honeycomb pads
- Flex-Edges and adjustable pivoting heel straps



ORDERING INFORMATION

Size	Left	Right
Youth/Pony	20100	20101
Adult	50100	50101

GEL ANKLE STIRRUP

INDICATIONS FOR USE

- Ankle injury rehabilitation / immobilization for grades I & II strains and sprains
- Chronic Instability
- Post Cast healing

FEATURES AND BENEFITS

- Safe hot/cold therapy
- Patented lace-lock feature for added stability
- Available 50% surface or 80%



ORDERING INFORMATION

Size	50% Gel Surface		80% Gel Surface	
	Left	Right	Left	Right
Youth/Pony	E200PL	E200PR		
Adult	E5050L	E5050R	E8020L	E8020R

EXOFORM® DORSAL NIGHT SPLINT

INDICATIONS FOR USE

- Plantar fasciitis
- Post-static pain
- General heel pain
- Drop foot
- Achilles tendonitis
- Lower extremity overuse injuries

FEATURES AND BENEFITS

- Adjustable aluminum stays can be positioned to maintain flexion between 85° – 90°
- Soft, Flex edge over molding improves patient comfort
- Full calf and heel strapping prevent brace from rotating
- Slip resistant padding



REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

Size	Men's Shoe Size US	Women's Shoe Size US	Shoe Size UK	Shoe Size EU	Part#
Medium	5 – 9	6 - 10	4.5 – 8	37.5 – 42.5	W-50085
Large	9.5 – 14	10.5 – 14	8.5 – 13.5	43 – 48.5	W-50087

FORMFIT® NIGHT SPLINT

INDICATIONS FOR USE

- Plantar fasciitis
- Lower extremity overuse injuries
- Achilles tendinopathy

FEATURES AND BENEFITS

- Straps adjust from 0° – 10° of dorsiflexion
- Contoured for a comfortable fit
- Soft foam interior is washable and breathable
- quick release buckles for easy application
- Slip resistant tread



REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

Size	Shoe Size US		Shoe Size UK		Shoe Size EU		Part#
	Male	Female	Male	Female	Male	Female	
Small	Up to 7	Up to 7.5	Up to 6.5	Up to 5	Up to 39	Up to 38	50023
Medium	7.5-10	8-10.5	7 – 9.5	5.5 – 8	39.5 – 44	38.5 – 43	50025
Large	10.5 and up	11 and up	10 and up	8.5 and up	44.5 and up	43.5 and up	50027

ACCESSORIES

Description	Size	Part#
Foam Wedge		10700
Replacement Liner	Small	11760
Replacement Liner	Medium	11761
Replacement Liner	Large	11762

AIRFORM[®] NIGHT SPLINT

INDICATIONS FOR USE

- Plantar Fasciitis
- Achilles Tendonitis
- Drop foot
- Post-static pain

FEATURES AND BENEFITS

- Light weight, low profile design
- Easy to adjust
- Available in singles or 6 packs

REQUIRED MEASUREMENTS

- Shoe size




ORDERING INFORMATION

Size	Shoe Size US		Shoe Size UK		Shoe Size EU		Part#	
	Male	Female	Male	Female	Male	Female	Standard	6 Pack
Medium	5 – 9	6 – 10	4.5 – 8.5	3.5 – 7.5	37 – 43	36 – 42	50045	500456
Large	9.5 – 14	10.5 – 15	9 – 13.5	8 – 12	43.5 – 48.5	42.5 – 48	50047	500476

SQUARE TOE POST-OP SHOE

FEATURES AND BENEFITS

- Square Toe design to increase protection as square acts as bumper.
- Rocker bottom sole reduces pressure on plantar surface resulting in increased comfort and compliance.
- Strapless closure eliminates buckle pressure on the forefoot and can expand to accommodate dressing up to 17" (43cm) in circumference.



Shoe size.



ORDERING INFORMATION

Size	Shoe Size US		Part#	
	Male	Female	Male	Female
Small	6 - 8	4 - 6	18113	18103
Medium	8.5 - 10	6.5 - 8	18115	18105
Large	10.5 - 12	8.5 - 10	18117	18107
X-Large	12.5 - 14	10.5 - 15	18118	-

SOFT TOP POST-OP SHOE

INDICATIONS FOR USE

- For recovery after soft tissue procedures and post trauma

FEATURES AND BENEFITS

- The cotton-poly blend upper aids in breathability and comfort
- The contact closure system is designed for ease of use.

REQUIRED MEASUREMENTS

- Shoe size



ORDERING INFORMATION

Size	Shoe Size US		Shoe Size UK		Shoe Size EU		Part#	
	Male	Female	Male	Female	Male	Female	Male	Female
Pediatric	Preschool 10 - Youth 1						18002	
Small	4 - 6	5 - 6	3.5 - 5.5	2.5 - 3.5	35.5 - 38	35 - 36	18013	18003
Medium	6.5 - 8.5	7 - 8	6 - 8	4 - 5.5	38.5 - 42	36.5 - 38.5	18015	18005
Large	9 - 11	9 +	8.5 - 10.5	6 +	42.5 - 45	39 +	18017	18007
X-Large	11.5 - 14		11 - 13.5		45.5 - 48.5		18018	

MESH TOP POST-OP SHOE

INDICATIONS FOR USE

- For recovery after soft tissue procedures and post trauma

FEATURES AND BENEFITS

- Breathable mesh upper with reinforced heel counter and padded heel collar, more stability overall to the rear of the foot while reducing slippage
- Semi-rigid sole provides superior protection and control after osteotomies
- Extended strap and dual side closure system accommodates extra bulky dressings



REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

Size	Shoe Size US		Shoe Size UK		Shoe Size EU		Male	Female
	Male	Female	Male	Female	Male	Female		
Small	4 – 6	5 – 6	3.5 – 5.5	2.5 – 3.5	35.5 – 38	35 – 36	17013	17003
Medium	6.5 – 8.5	7 – 8	6 – 8	4 – 5.5	38.5 – 42	36.5 – 38.5	17015	17005
Large	9 – 11	9 +	8.5 – 10.5	6 +	42.5 – 45	39 +	17017	17007
X-Large	11.5 – 14		11 – 13.5		45.5 – 48.5		17018	

CANVAS ROCKER BOTTOM CAST SHOE

INDICATIONS FOR USE

- For use with a short leg cast

FEATURES AND BENEFITS

- Cushioned rocker bottom cast sandal
- Contact closure, open toe

REQUIRED MEASUREMENTS

- Insole length



ORDERING INFORMATION

Size	Insole length	Blue	Black
X-Small	7" (18 cm)	308XS	
Small	7.5 – 9" (19 – 23 cm)	308SM	B308SM
S/Medium	9 – 9.75" (23 – 25 cm)	308SM/MED	B308S/M
Medium	10 – 10.5" (25.5 – 26.5 cm)	308MED	B308MED
M/Large	10.5 – 11" (26.5 – 28 cm)	308MLG	B308MLG
Large	11 – 12" (28 – 30.5 cm)	308LG	B308LG
X-Large	12 – 13.75" (30.5 – 35 cm)	308XLG	B308XLG



SPECIALITY CARE,
CASTING & TRACTION



COLD RUSH®

INDICATIONS FOR USE

- For cold therapy pain relief, as needed

FEATURES AND BENEFITS

- Consistent cooling and excellent compression, even at high elevation.
- Brushless motor produces a whisper quiet unit that allows patients to get the rest they need to properly rehabilitate.
- Highest ice-to-water ratio in the industry and a design that produces less heat transfer from the motor allows therapy pain relief for up to 6 continuous hours.
- Unique locking lid prevents spills
- Insulated walls improve cooling



COLD RUSH® COMPACT

INDICATIONS FOR USE

- For anything that may benefit from cold therapy

FEATURES AND BENEFITS

- Convenient, quite and easy-to-use
- Single flow setting
- Easy fill instructions



ORDERING INFORMATION

Description	Cold Rush Part #	Cold Rush Compact Part #
Device Only USA	B-232000010	B-232100010
Device w/ Knee Pad	B-232000011	B-232100011
Device w/ Unv Pad	B-232000012	B-232100012
Device w/ Hip Pad LT	B-232000013	B-232100013
Device w/ Hip Pad RT	B-232000014	B-232100014
Device w/ Shoulder Pad SM	B-232000015	B-232100015
Device w/ Shoulder Pad LG	B-232000016	B-232100016
Device w/ Ankle Pad	B-232000017	B-232100017
Device w/ Rectangular Pad SM	B-232000018	B-232100018
Device w/ Lumbar Pad	B-232000019	B-232100019


ACCESSORY

Description	Part #
Replacement Power Cord	B-301090001


COLD RUSH® PADS

Compatible with Cold Rush and Cold Rush Compact.


Description	Part #
Cold Rush knee pad	B-232002000
Cold Rush hip pad, left	B-232602100
Cold Rush hip pad, right	B-232502100
Cold Rush shoulder pad, small	B-232002202
Cold Rush shoulder pad, large	B-232002204
Cold Rush universal pad	B-232002300
Cold Rush rectangular pad, 5" x 10" (12.7 cm x 25.4 cm)	B-232002402
Cold Rush lumbar pad, large	B-232002404
Cold Rush foot & ankle pad	B-232002500




Knee Pad




Hip Pad, Left




Hip Pad, Right




Shoulder Pad,
Small




Shoulder Pad,
Large




Universal Pad



Lumbar Pad



Rectangular Pad



Foot & Ankle Pad

ELASTIC BANDAGE®

ORDERING INFORMATION

Description	Size	U/M	Part #
2" Metal Tab	2" x 5 yds	10 ea	2020
3" Metal Tab	3" x 5 yds	10 ea	2021
4" Metal Tab	4" x 5 yds	10 ea	2022
6" Metal Tab	6" x 5 yds	10 ea	2023
2" Velcro	2" x 5 yds	10 ea	C5920
4" Velcro	4" x 5 yds	10 ea	C5930
3" Velcro	3" x 5 yds	10 ea	C5940
6" Velcro	6" x 5 yds	10 ea	C5960

CAST PADDING

ORDERING INFORMATION

Width	U/M	Synthetic	Cotton
2"	12 rolls	V-101003	V-110002
3"	12 rolls	V-101003	V-110003
4"	12 rolls	V-101004	V-110004
6"	6 rolls	V-101006	V-110006

STOCKINETTE

- 100% synthetic stockinette
- Tubular stockinette for use with fiberglass, 25 yard dispenser box

ORDERING INFORMATION

Width	U/M	Synthetic	Cotton
1"	1 ea	2101	2000A
2"	1 ea	2102	2001
3"	1 ea	2103	2002
4"	1 ea	2104	2003
6"	1 ea	2106	2005

ALOETOUCH® EXAM GLOVES

FEATURES AND BENEFITS

- Certified by the International Aloe Science Council. This certification guarantees that the aloe vera gel in Aloetouch exam gloves meets the highest standards for product purity, thereby providing 100% of aloe vera's natural benefits.
- Designed for people who frequently change gloves.



ORDERING INFORMATION

Size	Fits	U/M	Part #
Small	6 – 6.5	100/bx, 10 bx/cs	MDS195074
Medium	7 – 7.5	100/bx, 10 bx/cs	MDS195075
Large	8 – 8.5	100/bx, 10 bx/cs	MDS195076
XLarge	9 – 9.5	100/bx, 10 bx/cs	MDS195077

Note: Packed in boxes of 10

TECHFORM® PREMIUM

INDICATIONS FOR USE

- Immobilization

FEATURES AND BENEFITS

- Glides off the roll for easier application
- Soft, conformable rolls allow for easy handling
- “Rub-responsive“ resin formula adapts to the technician’s technique – the more it is rubbed, the smoother the cast gets
- Excellent lamination ensures a strong and durable cast
- Smoother finished cast is non-abrasive and will not snag on clothing
- Precision welded ends do not fray and lay down easily
- Finger holes in the pouch allow for easy opening



ORDERING INFORMATION

WIDTH	WH	BK	OR	RD	PK	HP	LB	DB	PR	GN	RB
1.5"	154WH										
2"	204WH	204BK	204OR	204RD	204PK	204HP	204LB	204DB	204PR	204GN	204RB
3"	304WH	304BK	304OR	304RD	304PK	304HP	304LB	304DB	304PR	304GN	304RB
4"	404WH	404BK	404OR	404RD	404PK	404HP	404LB	404DB	404PR	404GN	404RB
5"	504WH										

Note: Packed in boxes of 10



TECHFORM® ORTHO-ROLL®

Developed in collaboration with orthopaedic professionals, Techform Ortho-Roll represents the first major breakthrough in orthopaedic splinting technology in decades.

In focus groups, Techform Ortho-Roll rated from "Good" to "Excellent" in every key category, including handling, conformability, padding, water retention, dispensing method, ease-of-use, overall appearance, and overall satisfaction. Techform Ortho-Roll simply performs better.

FEATURES AND BENEFITS

- Quickdry is a patented, revolutionary orthopaedic padding material. Quickdry padding never really stays wet. Water simply passes straight through it. The result is less risk of skin maceration. Plus, it saves you application time.
- Instructions are right on the box
- Box opens easily along perforation
- Closure rod creates an air tight seal



ORDERING INFORMATION

TECHFORM ORTHO-ROLL

Part#	Size	Part#	Size
TR0115	1" x 15ft.	TR0415	4" x 15ft.
TR0215	2" x 15ft.	TR0515	5" x 15ft.
TR0315	3" x 15ft.	TR0615	6" x 15ft.

ACCESSORIES

Part#	Size
13302	Positioning Strap (12 pack)
13332	Closure Rod
P-550606	Splinting Guide Poster

TECHFORM® PRE-CUT SPLINTS

Techform Pre-Cut Splints have the same key features as our patented Ortho-Roll splinting material, including our patented Quickdry padding. Each splint comes pre-packaged in an individual foil pouch, perfect for those who prefer to eliminate waste of splint material.

FEATURES AND BENEFITS

- Our selection includes 5 sizes, a full variety of lengths and widths to accommodate all your splinting needs.
- Pre-cut splints come in various lengths so you can choose the appropriate size to match your specific needs



ORDERING INFORMATION

Part#	Size	Quantity
TC0312	3" x 12"	10/Box
TC0335	3" x 35"	5/Box
TC0415	4" x 15"	10/Box
TC0430	4" x 30"	5/Box
TC0530	5" x 30"	10/Box



FOLLOW ÖSSUR ON



USA (800) 233-6263
WWW.OSSUR.COM

CANADA (800) 663-5982
WWW.OSSUR.CA

