

FACILITY INFORMATION (Billing)			
Account #	Shipping same as billing address		
Clinician Name	Ship to Address		
Company Name			
Company Address			
City	State	Zip	
Telephone			
Purchase Order #			
Patient Identifier	Ship Via	Next Day	2 Day
Please provide email address or fax number for order confirmation			



DIGITS MODEL (Please indicate type)		AFFECTED SIDE(S)	
i-Digits Quantum	i-Digits Access	Left	Right
1 digit TBX50008	1 digit TBX50006		
2 digits TBX50016	2 digits TBX50014		
3 digits TBX50024	3 digits TBX50022		
4 digits TBX50032	4 digits TBX50030		
5 digits TBX50040	5 digits TBX50038		



**Extra Small**  
Extra Small mounting plate with a 5, 5, 3, 2 digit configuration will be supplied as standard. If an alternative configuration is required, then please note digit strength may be compromised.



**Medium**  
Medium mounting plate with a 6, 6, 5, 3 digit configuration will be supplied as standard. If an alternative configuration is required, then please note digit strength may be compromised.



**INDICATE DIGITS AND SIZE**

**Left Prosthesis Configuration**

- thumb
- index
- middle
- ring
- pinky

**Right Prosthesis Configuration**

- thumb
- index
- middle
- ring
- pinky

INPUT OPTIONS	
50Hz Remote	
50Hz Compact	
Large FSR	
Small FSR	

ACCESSORIES		
iPod: (included with i-digits quantum)	Yes	No
Covers: Black Clear		

CHOOSE WARRANTY		
Year 3	Year 4	Year 5

**NOTES:**

## FABRICATION PROCESS

Please take impressions using the instructions below. Watch the video at [youtu.be/kfevZeYOvCo](https://youtu.be/kfevZeYOvCo) or by reading the QR-code further down on this page.

Note: It is recommended that you refrigerate the silicone 24 hours ahead of taking the impressions.

### Step 1:

Apply Vaseline to the affected hand and arm up to the elbow. If patient has hair, pat-down with Vaseline and then wipe the excessive Vaseline off with a paper towel.




### Step 2:

Take impression up to the elbow crease on the affected side.



### Step 3:

When the silicone is completely cured, remove the mold. Have the patient wiggle their digit(s) or wrist gently while you pull the mold off. The silicone can be ripped slightly at the top for ease of removal.



SAVE | PRINT | CLEAR FORM | PREVIOUS | NEXT

page 2 of 2

FOLLOW ÖSSUR ON



ÖSSUR NORDIC  
+46 18 18 22 00  
[info@ossur.com](mailto:info@ossur.com)

[WWW.OSSUR.SE](http://WWW.OSSUR.SE)

