



UNLOADER ONE® X

Clinically proven
osteoarthritis pain relief

There's only one Unloader®. With more than 15 clinical studies* demonstrating their efficacy, Unloader OA knee braces have been the industry standard for close to 4 decades. And they just got even better with the latest addition to the Unloader family, Unloader One® X.

Designed with patient compliance in mind. Like our popular Unloader One® brace, Unloader One X features the Össur 3-Point Leverage System. This patented system is clinically proven to improve mobility by providing pain relief from unicompartmental knee OA or degenerative meniscal tears. Equally important, Unloader One X incorporates a host of new features (see below) designed to improve patient comfort and compliance. Because even the best brace isn't effective if it isn't worn.

PRODUCT HIGHLIGHTS

- Stylish, low-profile frame with soft edges
- Breathable and skin friendly soft goods
- Machine washable for easy cleaning
- Streamlined, integrated SmartDosing® system for easy self-dosing of pain relief
- Color-coded patient touch points for intuitive use
- Redesigned Quick-Fit buckles with lever arms for easy closure, even by users with limited strength


*data on file at Ossur

The Unloader One X is easily adaptable to individual needs and anatomical shapes, with the help of various accessories and replacement parts. For ordering information, consult the Unloader One X catalog page or visit ossur.com/unloader-one-x

UNLOADER ONE® X SIZING

Size	Measurements (in)	Measurements (cm)
XSmall	9.8-12.2"	25-31cm
Small	12.2-14"	31-36cm
Medium	14-16.75"	36-43cm
Large	16.75-20"	43-51cm
XLarge	20-24"	51-62cm

Circumference: 6"/15cm below mid-patella



ORDERING INFORMATION

STANDARD

Size	Right Medial	Left Medial	Right Lateral	Left Lateral
XSmall	UOX0518711	UOX0618711	UOX0528711	UOX0628711
Small	UOX0518712	UOX0618712	UOX0528712	UOX0628712
Medium	UOX0518713	UOX0618713	UOX0528713	UOX0628713
Large	UOX0518714	UOX0618714	UOX0528714	UOX0628714
XLarge	UOX0518715	UOX0618715	UOX0528715	UOX0628715

Recommended for patients between 5'5" - 6'4" (1.65m - 1.92m).
Hinge is placed on the side of the affected compartment.

SHORT

Size	Right Medial	Left Medial	Right Lateral	Left Lateral
XSmall	UOX1518711	UOX1618711	UOX1528711	UOX1628711
Small	UOX1518712	UOX1618712	UOX1528712	UOX1628712
Medium	UOX1518713	UOX1618713	UOX1528713	UOX1628713
Large	UOX1518714	UOX1618714	UOX1528714	UOX1628714
XLarge	UOX1518715	UOX1618715	UOX1528715	UOX1628715

Recommended for patients under 5'5" (1.65m).
Hinge is placed on the side of the affected compartment.

ACCESSORIES

Description	Part #	Size	Configuration
ActiveGrip Calf Liner	UOX516153	XS/S	LL/RM
ActiveGrip Calf Liner	UOX515153	XS/S	RL/LM
ActiveGrip Calf Liner	UOX516163	M/XL	LL/RM
ActiveGrip Calf Liner	UOX515163	M/XL	RL/LM
Thigh Extension Liner	UOX545101	NA	RL/LM
Thigh Extension Liner	UOX546101	NA	LL/RM

EASE OF USE

Description	Part #
Sundial + Finger Loops	UOX1200001
Sundial	UOX800001



Sundial (for easier grip on SmartDosing dial)



Finger Loops (for an even easier opening of buckles)

