

NEW

Electronically rotating thumb glides
effortlessly into grip position

18 different grip options available




touch
bionics BY ÖSSUR

i-limb® ultra

An externally powered, multi-articulating prosthetic hand with five individually powered digits and powered thumb rotation, with manual override.

- Individually motorized digits with compliant grip and stall-out ability
- Speed boost:
Adjustable open/close speed for digits
- Auto-grasp:
Anti drop safety feature for added reassurance
- Vari-grip:
Feature for enhanced grip strength
- Multiple covering options available to suit individual lifestyles
- Titanium digits
Check out our titanium strengthened digits which increase the maximum carry load by 50%. Available in S/M/L.

METHODS OF CONTROL:





ÖSSUR
TOUCH
SOLUTIONS

Össur Touch Solutions are designed to ensure the best possible outcomes for people with upper-limb deficiencies.

Technical Information				
Voltage	7.4 V (nominal)			
Max. Current	5 A			
Battery Capacity	Rechargeable lithium polymer; 7.4 V (nominal); 2000 mAh capacity; 1300 mAh capacity			
Max. hand load limit (static limit)	40kg/88lb (Extra Small) 90kg/198lb (Small/Medium/Large)			
Finger carry load (static limit)	20kg/44lbs (Extra Small) 32kg/71lbs (Small/Medium/Large)			
Time from open position to full power grip	0.8 seconds			
Device Weight		Extra Small	Small	Medium/Large
	QWD	472g/1.04lb	512g/1.13lbs	528g/1.16lbs
	WD	432g/0.96lbs	472g/1.04lbs	488g/1.08lbs
	Flexion Wrist	572g/1.26lbs	612g/1.35lbs	628g/1.38lbs
	Friction Wrist	467g/1.03lbs	507g/1.12lbs	523g/1.15lbs

Note: Titanium digits add an additional 30g/0.07lbs per hand

Standard Package	
Part #	Description
TBX50185	PK1 i-limb Ultra QWD/WD/FRIC
TBX50186	PK2 i-limb Ultra Flexion

Titanium Package	
Part #	Description
TBX50482	PK1 i-limb Ti Ultra QWD/WD/Fr
TBX50483	PK2 i-limb Ti Ultra Flexion

FOLLOW ÖSSUR ON
f t i YouTube in

WWW.OSSUR.ASIA

© ÖSSUR, P-600084, 1. 2018

Össur APAC
2F, W16 B
No. 1801 Hongmei Road
200233, Shanghai, China

TEL: 86 21 6127 1707
FAX: 86 21 6127 1799
asia@ossur.com

 **ÖSSUR**
LIFE WITHOUT LIMITATIONS