

touch  
bionics BY ÖSSUR

## i-digits™ access

design. dexterity.  
intelligent motion.

Fully customized, multi-articulating partial hand prosthesis.

- Anywhere from 1-5 digits can be replaced
- Up to 12 grips available
- Form fitting anatomical design
- Self-contained wristband houses batteries & microprocessor
- Digits move independently and work in conjunction with remaining fingers
- Full wrist movement




ÖSSUR


## TOUCH SOLUTIONS

Össur Touch Solutions are designed to ensure the best possible outcomes for people with upper-limb deficiencies.

i-digits™ access is appropriate for partial hand absence where the level of loss or deficiency is distal to the wrist and proximal of the metacarpophalangeal joint. Anywhere from one to five digits can be replaced.



### Digit positionig options:



### Digit size options:



FOLLOW ÖSSUR ON



WWW.OSSUR.ASIA

Össur APAC  
2F, W16 B  
No. 1801 Hongmei Road  
200233, Shanghai, China

TEL: 86 21 6127 1707  
FAX: 86 21 6127 1799  
asia@ossur.com

