


Össur Custom Solutions Ohio 6640 Riverside Drive, Suite 360 Dublin, OH 43017 P: (949) 382-3775 F: (949) 271-2390

FACILITY INFORMATION (Billing)			
Account #	Shipping same as billing address		
Clinician Name	Ship to Address		
Company Name			
Company Address			
City	State	Zip	
Telephone			
Purchase Order #			
Sales Quotation#	Ship Via	Next Day	2 Day
Please provide email address or fax number for order confirmation			

PATIENT INFORMATION	
Patient Identifier	Insurance Company



DIGITS MODEL (Please indicate type)		AFFECTED SIDE(S)	
Quantum	Access	Left	Right



Extra Small
Extra Small mounting plate with a 5, 5, 3, 2 digit configuration will be supplied as standard. If an alternative configuration is required, then please note digit strength may be compromised.



Medium
Medium mounting plate with a 6, 6, 5, 3 digit configuration will be supplied as standard. If an alternative configuration is required, then please note digit strength may be compromised.



INDICATE DIGITS AND SIZE	
Left Prosthesis Configuration	Right Prosthesis Configuration
thumb	thumb
index	index
middle	middle
ring	ring
pinky	pinky

INPUT OPTIONS	
60Hz Remote	
60Hz Compact	
Large FSR	
Small FSR	

ACCESSORIES		
iPod: (included with i-digits quantum)	Yes	No
Covers:	Black	Clear
Add Ons:	Livingskin bundle (passive functional prosthesis)	
Extended Touch Care:	3	4 5

NOTES:

FABRICATION PROCESS

Please take impressions using the instructions below. Watch the video at youtu.be/kfevZeYOvCo

Note: It is recommended that you refrigerate the silicone 24 hours ahead of taking the impressions.

Step 1:

Apply Vaseline to the affected hand and arm up to the elbow. If patient has hair, pat-down with Vaseline and then wipe the excessive Vaseline off with a paper towel.



Step 2:

Take impression up to the elbow crease on the affected side.



Step 3:

When the silicone is completely cured, remove the mold. Have the patient wiggle their digit(s) or wrist gently while you pull the mold off. The silicone can be ripped slightly at the top for ease of removal.



The impressions will be used to create a clear silicone socket with PETG frame.

Össur Academy in Ohio Central Fabrication

Any motion limitations or requests for diagnostic fabrication:



Ship impressions to: Ossur Custom Solutions Ohio 6640 Riverside Drive, Suite 360 Dublin, OH 43017
P: (949) 382-3775 F: (949) 271-2390

SAVE | PRINT | CLEAR FORM | PREVIOUS | NEXT

page 2 of 2

FOLLOW ÖSSUR ON



USA (800) 233-6263
WWW.OSSUR.COM

CANADA (800) 663-5982
WWW.OSSUR.CA

