

ÖSSUR ARMS - i-DIGITS FABRICATION ORDER FORM

Return shipments to: Össur Custom Solutions, 6640 Riverside Dr, Suite 360, Dublin, OH 43017
 For other inquiries: touch.fabrication@ossur.com



FACILITY INFORMATION (Billing)

Account #	Shipping same as billing address		
Clinician Name	Ship to Address		
Company Name			
Company Address			
City	State	Zip	
Telephone			
Purchase Order #			
Sales Quotation#	Ship Via	Next Day	2 Day
Please provide email address or fax number for order confirmation			

REQUIRED PATIENT INFORMATION SIDE(S)

Last Name:	First Name:	LEFT	RIGHT
------------	-------------	------	-------

INCLUDED ITEMS	ALTERNATIVE PARTS USED								
Kit includes: 1. Fitting instructions manual 2. Silicone diagnostic socket and frame 3. Plaster model to modify 4. Digit(s) with plate for alignment 5. Alternative Parts Kit *All items in the kit must be returned with the diagnostic package or additional fees will result.	CHECK USED PARTS	DESCRIPTION	QUANTITY IN KIT	CHECK USED PARTS	DESCRIPTION	QUANTITY IN KIT	CHECK USED PARTS	DESCRIPTION	QUANTITY IN KIT
		Nothing Used	X		Thumb Size 5	1		Flat Mounting Plate R	1
		Size 6 digit	1		Bellows	1		Flat Mounting Plate L	1
		Size 5 digit	1		Wristband	1		LP Remote Electrodes	2
		Size 3 digit	1	L R	Stepped Mounting Plate Medium	1/ea		Compact Electrodes	2
		Size 2 digit	1	L R	Stepped Mounting Plate XS	1/ea			

NEXT STEP

DIAGNOSTIC PROSTHESIS	REQUESTS	NOTES
	Bilateral user modifications to improve independence Position anchor points for radial suspension strap Other	

DEFINITIVE PROSTHESIS	SOCKET COLOR	FRAME COLOR	FRAME MATERIAL
	Black Gray Battery Color Other such as Red or Blue <i>(Please include example print-out or swatch of color)</i> Hand-mixed Skin Tone based on Swatch # by	Color (i.e. Black, Blue, Red, etc.) Otto Back Caucasian Freds Legs Style bit.ly/SleeveArt + \$40 Other <i>If using a swatch provider other than Touch Solutions, please send us a physical swatch</i> Hand-mixed Skin Tone based on Swatch # by	Fiberglass Lamination Other <i>Note: Carbon fiber is conductive and cannot be used with myoelectric devices.</i> <b style="background-color: #0056b3; color: white;">FRAME FINISH Smooth Finish Weave Featured (Carbon Look) Finish Notes:
	OTHER		
	Bilateral user modifications to improve independence Position anchor points for radial suspension strap	Zipper modification to prevent skin pinching Material in palm for increased friction during grasp	
	NOTES		

CHECKLIST

INITIAL TO CONFIRM	REQUIREMENTS (See associated "Fitting Instructions" pages)
	I mounted the electrodes in the silicone socket (pg2)
	I aligned and mounted the digits to the diagnostic frame (pg3)
	I modified the positive model (pg4)
	I marked the trimlines and electrode locations on the positive model (pg4)
	I packaged the positive model, socket, frame, componentry, and paperwork for return without damage