


Sri Ganesha Healthcare

FULL HD CAMERA

Model: SGH-Full HD



Intelligent Digital Full HD Camera

- Resolution: 1920 (H) × 1080 (V) effective pixels
- Sensor: 1/2.8" CMOS
- TV Lines: Over 1100

Full Intelligent automatic white balance, Image-freeze technology, brightness adjustment and surgery modes for recommending settings for various surgeries with one click on camera head.

- Aspect ratio 16:9 with auto shutter and brightness control.
- Water proof camera head IP8X to prevent the water during surgery.
- Progressive scan with brightness adjustment.
- In-built USB continues Recording / photos .

Outputs: HDMI, DVI, BNC and SDI

Note: Technical specification may/may-not change without any prior information.

