



LUCAS[®] 2

CHEST COMPRESSION SYSTEM

Genuine Accessories from Physio-Control



**Thank you for choosing Physio-Control
as your partner in helping save lives and
improve patient care.**

The LUCAS 2 chest compression system is shipped with one battery, patient straps, three suction cups, a carrying bag and the Instructions for Use. In this catalog you will find the overview of accessories designed to meet your needs and optional power solutions. If you have any questions about accessories or disposables, please contact your Physio-Control representative.

GENERAL ACCESSORIES AND DISPOSABLES



LUCAS Carrying Bag

LUCAS 2 is stored and carried in a padded carrying bag complete with storage pouches for a spare battery and external power supply.
11576-000038



LUCAS 2 Disposable Suction Cup

LUCAS 2 has a disposable suction cup which can be easily changed after each use.
11576-000046 (3-Pack)
11576-000047 (12-Pack)



LUCAS Stabilization Strap

The cushion and straps are used to provide further stabilization relative to the patient and to prevent downward movement of the LUCAS device during operation.
21996-000064
11576-000037 (4-Pack)



LUCAS 2 Rubber Bumper (Pair)

The rubber bumper accessory is available to protect the hinge end caps.
11576-000070



LUCAS Patient Straps

The LUCAS patient straps (two required per device) are used to secure the patients arms to the support legs of LUCAS during transportation.
11576-000050 (Pair)
11576-000051 (Pair, 3-Pack)



LUCAS Back Plate Grip Tape

The LUCAS back plate grip tape is designed to increase the friction of the back plate and to reduce the risk of the device slipping when in contact with other hard surfaces.
11576-000052
11576-000053 (3-Pack)



Standard LUCAS Back Plate (Replacement)

Replacement LUCAS Back Plate is available in the event of damage or contamination of original.
21996-000044

BATTERIES & POWER ACCESSORIES



LUCAS 2 Battery

The LUCAS 2 battery is a rechargeable Lithium Polymer (LiPo) battery which typically last for 45 minutes of operation. You can charge the battery in less than 4 hours either in the external battery charger or while installed in the LUCAS device and connected to an external power supply.
11576-000039
11576-000040 (4-Pack)



LUCAS 2 Battery Charger

The external battery charger is a stand-alone unit intended to charge one LUCAS 2 Battery. It typically takes less than 4 hours to charge the battery. A mounting bracket is included for space saving installation.
11576-000060 (U.S., Canada, Japan)
11576-000061 (EU)
11576-000062 (UK)
11576-000063 (Australia)
11576-000068 (Japan)



LUCAS 2 Aux Power Supply (S/N 3009 0181 and above)

The external power supply (100-240VAC, 50/60 Hz) can be connected to the LUCAS 2 device and supports prolonged operation as well as charging of the LUCAS 2 Battery while in the device.
11576-000055 (U.S.)
11576-000056 (EU)
11576-000057 (UK)
11576-000058 (Japan)
11576-000059 (Australia)



LUCAS 2 Car Cable

The car cable (12-24V DC) can be connected to the LUCAS 2 device and supports prolonged operation as well as charging of the LUCAS 2 Battery while in the device.
11576-000048

PCI APPLICATIONS




LUCAS PCI Back Plate

The carbon fiber LUCAS PCI Back Plate is intended specifically for use in the cath lab. It is fully radiotranslucent, with minimum shadows and no metal screws.

11576-000064

SOFTWARE UPGRADES




LUCAS 2, Software Version 2.1 Upgrade Kit

Available from a certified Physio-Control Service Technician, the software version 2.1 upgrade includes a quick fit feature which allows for optimal height placement of the pressure pad and earlier recognition and notification if a patient is too small.

43576-000001 (U.S.)

43576-000002 (International English)

INSTRUCTIONS FOR USE AND TRAINING



Instructions for Use - Replacement (Software version 2.1)

- 26500-003434 English (EN)
- 26500-003433 International English (INTL EN)
- 26500-003430 Croatian (HR)
- 26500-003431 Chinese (ZH)
- 26500-003435 German (DE)
- 26500-003436 French (FR)
- 26500-003437 Dutch (NL)
- 26500-003438 Danish (DA)
- 26500-003439 Swedish (SV)
- 26500-003440 Norwegian (NO)
- 26500-003442 Italian (IT)
- 26500-003441 Spanish (ES)
- 26500-003443 Czech (CS)
- 26500-003444 Polish (PL)
- 26500-003445 Finnish (FI)
- 26500-003446 Greek (EL)
- 26500-003447 Portuguese (PT)
- 26500-003448 Hungarian (HU)
- 26500-003449 Russian (RU)
- 26500-003450 Turkish (TU)
- 26500-003455 Japanese (JA)

Instructions for Use - Replacement (Software version 2.0)

- 26500-003084 English (EN)
- 26500-003085 International English (INTL EN)
- 26500-003086 German (DE)
- 26500-003087 French (FR)
- 26500-003088 Dutch (NL)
- 26500-003090 Danish (DA)
- 26500-003089 Swedish (SV)
- 26500-003091 Norwegian (NO)
- 26500-003092 Italian (IT)
- 26500-003093 Spanish (ES)
- 26500-003094 Czech (CS)
- 26500-003095 Polish (PL)
- 26500-003108 Finnish (FI)
- 26500-003109 Greek (EL)
- 26500-003110 Portuguese (PT)
- 26500-003111 Hungarian (HU)
- 26500-003184 Japanese (JA)
- 26500-003428 Russian (RU)
- 26500-003429 Turkish (TU)



Product DVD

The LUCAS product DVD features an introduction on the physiology of CPR with Tim Phalen followed by a LUCAS 1 Inservice video and a LUCAS 2 Inservice video.

A great DVD to have for new or refresher training.

This product must be ordered through your Physio-Control Sales Representative (reference GDR 3304502)

For further information, please contact Physio-Control at 800.442.1142 (U.S.), 800.895.5896 (Canada) or visit our website at www.physio-control.com.



Physio-Control Headquarters
Redmond, WA 98052
www.physio-control.com

Customer Support
P. O. Box 97006
Redmond, WA 98073
Toll free 800 442 1142
Fax 800 426 8049

Physio-Control Canada
Physio-Control Canada Sales, Ltd.
Toll free 800 895 5896
Fax 866 430 6115