

(ROM) walker

S-FIT-91-006 & 007 | Walking boot

Indications

- Post ligament, soft tissue and tendon repairs
- Conservative or post-operative treatment of fractures

Features & benefits

- Can block the ankle or limit range of motion
- Plantar flexion: 0°, 10°, 20°, 30° and 40°
- Dorsal flexion: 20°, 10° and 0°
- Special sole design for natural movement
- Set of three wedges available for Achilles tendon relief



This brace is part of our own FITLINE brand.

Orthopaedic braces and shoes, designed to obtain a maximum of comfort and to maintain a normal lifestyle.

Move - feel - live !

 **Fitline®**



Reference
A. ROM Walker

Reference
B. Walker

Size

Shoe Size

S-FIT-91-006-U2

S-FIT-91-007-U2

S

32-37

S-FIT-91-006-U3

S-FIT-91-007-U3

M

38-42

S-FIT-91-006-U4

S-FIT-91-007-U4

L

43-48



Moovartes bv - info@moovartes.com – www.moovartes.com

Tel. BE: +32 (0)89 710 400 – Tel. NL: +31 (0)46 4369110