

ES-EX022180 | Diabetic shoe with wide fit.

Indications

- Post-operative treatment of **hallux valgus** (minimal-invasive protocol)
- Protection of bone structures and surrounding tissues of the forefoot
- Post-operative treatment after a **Morton's neuroma**
- **Swelling and injury of the forefoot** that make wearing a shoe impossible
- **Protection and support** of foot, toes and forefoot in an adequate position

Features & benefits

- Lightweight shoe with an **ergonomic wide fit**
- Made to be easily worn with a cast or bandage
- **Highly adjustable** thanks to the velcro strap
- **Ambidextrous** design: can be worn left and right



Also available
with the Plantar
S-Cloud insole



Ordering information

Reference Biomed with basic S-Air insole	Reference Biomed with S-Cloud insole	Reference Toe cap	Shoe size	Size
ES-EX022180-XS	ES-EX022181-XS	ES-CT-BIOMED-XS	35 - 36	XS
ES-EX022180-S	ES-EX022181-S	ES-CT-BIOMED-S	37 - 38	S
ES-EX022180-M	ES-EX022181-M	ES-CT-BIOMED-M	39 - 40	M
ES-EX022180-L	ES-EX022181-L	ES-CT-BIOMED-L	41 - 42	L
ES-EX022180-XL	ES-EX022181-XL	ES-CT-BIOMED-XL	43 - 44	XL
ES-EX022180-XXL	ES-EX022181-XXL	ES-CT-BIOMED-XXL	45 - 46	XXL

