

# tecnicca heel

## ES-EX023000 | Diabetic shoe with movable heel

### Indications

- Post-surgery rehabilitation and/or wounds or ulcers of the heel area
- Post-surgery fractures and/or calcaneus ulcers (second treatment) and ulcers at the insertion of the Achilles tendon
- Bursitis

### Features & benefits

The balance of the weight is moved to the front of the foot thanks to the **unloading insole** and the outsole, which is designed to totally release the heel from the pressure load.

- **Square shape** offers space for bandages
- **Movable rear-cover** to preserve the heel from impacts
- Non-slip outsole to ensure safety
- Total opening with velcro-straps
- Removable insole



**Movable  
rear-cover**



### Ordering information

Reference	Shoe size	Size
ES-EX023000-S	34 - 36	S
ES-EX023000-M	37 - 39	M
ES-EX023000-L	40 - 42	L
ES-EX023000-XL	43 - 45	XL

