

Telescopic post-op elbow brace with graduated movement

Indications

- Compound fractures of elbow, distal humerus, proximal radius and ulna
- Post-operative immobilisation following tendon and ligament repairs
- Useful in the rehabilitation phase following surgical and conservative treatment of elbow fractures
- Elbow subluxations and dislocations
- Chronic elbow injuries
- Elbow hyperextension
- Post-traumatic treatment



Flexion from -10° to $+120^{\circ}$
with intervals of 10°

Extension from -10° to $+90^{\circ}$
with intervals of 10°

Fast lock joint from -10° to $+90^{\circ}$
with intervals of 10°



Features & benefits

- **Telescopic and easy to adapt** to the patient thanks to the easy-to-use buttons on the valves
- **Lateral articulated hinge** in a light aluminium alloy
- New **"PULL&SET"** range of motion control for easy adjustment
- **Maximum comfort** thanks to the heat-sealed padding in soft lycra that comes into contact with the skin



Ordering information

Ref. M24-017 - choose left or right side

Moovartes bv - info@moovartes.com - www.moovartes.com

Tel. BE: +32 (0)89 710 400 – Tel. NL: +31 (0)46 4369110



**ORTHOSERVICE
RO+TEN**

Take care feel better