

# malleorapid

## Ankle brace with steel stays and figure-8 elastic straps.

Malleorapid is recommended for **II degree sprain trauma** and as an adjuvant in the resumption of **sporting activity**. The combined action of the steel stays and the 8-figure elastic straps limits inversion/eversion of the ankle.

### Indications

- Conservative treatment following a grade 2 sprain trauma
- Acute tendinopathies
- Post-surgical rehabilitation
- Aid to resuming sporting activity after grade 2 and 3 sprain trauma
- Chronic joint pain caused by tibiotarsal arthrosis



### Features & benefits

- Easy to wear thanks to open heel closure
- Simple and quick drawstring adjustment system
- Anatomical and lowprofile design for use with most shoes
- Multifunctional use, even during sports activities

### Ordering information

Fits right and left

Reference	Size	Ankle circumference
OS-24081-S	S	19 - 22 cm
OS-24081-M	M	22 - 26 cm
OS-24081-L	L	26 - 30 cm

