

ROM walker

Walking boot with range of motion

Indications

- Malleolar fractures conservative treatment
- Grade 3 ankle sprains
- Post-surgical care for osteosynthesis of malleolar fractures, ankle ligament reconstruction and Achilles tendon tenorrhaphy
- Walking recovery in consolidation delay of tibial fractures
- Charcot foot

Features & Benefits

- Structure made of durable, light plastic
- Non-slip rubber sole
- Inner padding of the sole made from non-toxic Evazote®
- Struts with Velcro for padding grip
- Soft foam pad with double lining in brushed nylon for Velcro® grip
- Velcro® fastening tapes
- Washable padding

Range of Motion

Plantar & dorsal flexion:

0° - 7,5° - 15° - 22,5° - 30° - 37.5° - 45°

Ordering information

Fits right and left

Reference	Size	Shoe Size
OS-M600-S	S	36 - 39
OS-M600-M	M	39 - 44
OS-M600-L	L	44 - 47



Moovartes bv - info@moovartes.com - www.moovartes.com

Tel. BE: +32 (0)89 710 400 – Tel. NL: +31 (0)46 4369110



ORTHOSERVICE
RO+TEN

Take care feel better