

Airstep Walker

Walker with inflatable air padding

Indications

- Conservative treatment of of malleolar fractures
- Severe sprains (grade 3)
- Post-surgical phase of osteosynthesis of malleolar fractures
- Ligamentous reconstructions of the ankle and of the tenorrhaphy of the Achilles tendon
- Stable fractures of the foot, ankle, tibia and fibula distal

Contra-indications

The application of this brace is advised only when it is documented that the fracture is stable and there are acceptable limits of angular and rotational misalignments

Features & benefits

- **Reinforcing shells** on the front and back, for increased stability
- **Inflatable air pads** for optimum fit to the limb and improved stability
- Structure in lightweight, sturdy plastic
- Fits right and left

Ordering information

Reference	Size	Shoe size
OS-M500-S	S	35 - 39
OS-M500-M	M	39 - 44
OS-M500-L	L	44 - 47
OS-M500-XL	XL	> 47

WEDGES

Reference	Side
OS-M501L	Left
OS-M501R	Right

DIABETIC SETS

Reference	Size	Shoe size
OS-M505-S	S	35 - 39
OS-M505-M	M	39 - 44
OS-M505-L	L	> 44

