

# Airstep Tight Walker

Walker with **two** inflatable air pads

**NEW**

## Indications

- Conservative and post-surgical treatment of foot/ankle and leg fractures
- Serious bone distortion grade 1 to 3
- Conservative treatment after a tear in the Achilles tendon
- Post-surgical treatment after tenorrhaphy of the Achilles tendon
- Pre- and post-operative/post-traumatic immobilisation

## Contra-indications

The application of this brace is advised only when it is documented that the fracture is stable and there are acceptable limits of angular and rotational misalignments

## Features & benefits

- **Reinforcing shells** on the front and back, for increased stability
- **Two inflatable air pads** for optimum fit to the limb and improved stability
- Structure in lightweight, sturdy plastic; terminal part in rubbery material
- **Integrated pump** to inflate/deflate the cells
- Sole in non-slip material
- Optional padding for additional adaptation to the leg
- Fits right and left

## Ordering information

Reference long	Reference short	Size	Shoe size
OS-M550N-S	OS-M560N-S	S	35 - 39
OS-M550N-M	OS-M560N-M	M	39 - 44
OS-M550N-L	OS-M560N-L	L	44 - 47
OS-M550N-XL	OS-M560N-XL	XL	> 47

## WEDGES

Reference	Side
OS-M501L	Left
OS-M501R	Right

## DIABETIC SETS

Reference	Size	Shoe size
OS-M505-S	S	35 - 39
OS-M505-M	M	39 - 44
OS-M505-L	L	> 44



**SHORT version**



**Moovartes bv** - [info@moovartes.com](mailto:info@moovartes.com) - [www.moovartes.com](http://www.moovartes.com)

Tel. BE: +32 (0)89 710 400 – Tel. NL: +31 (0)46 4369110



**ORTHOSERVICE  
RO+TEN**

Take care feel better