

GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis

New
BX Stent
Graft




Flexible Strength. Proven Success.




BALLOON EXPANDABLE
ENDOPROSTHESIS

100% RESTORATION OF LUMEN DIAMETER




Before



After

≤ 30% residual stenosis due to high radial strength,
even in highly calcified and non-compliant lesions*

Flexible Strength. Proven Success.



VIABAHN®
VBX


BALLOON EXPANDABLE
ENDOPROSTHESIS

100% STENT DELIVERY


-
- 100% delivery to target lesion
 - 100% stent retention
 - 100% deployment at the target site
 - 234 devices delivered
 - 50% bilateral treatment
 - 42% kissing stents
 - 32% TASC II C & D including occlusions
 - 18% contralateral deliveries
 - predilatation not required

Flexible Strength. Proven Success.

96.9% PRIMARY PATENCY



95.3% Primary Patency in
TASC II C & D lesions



Expand to every demand

The GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis (VBX Stent Graft) is designed to offer a unique combination of technical and clinical benefits in a broad range of peripheral vessels. This advanced BX stent graft, an evolution of proven vascular solutions, provides precise delivery and positive outcomes that physicians can rely on for complex cases.*

Flexible Strength. Proven Success.

Advanced technology and unique design



Stainless steel independent rings enhance flexibility

- Minimizes foreshortening
- Provides high radial strength




Semi-compliant covered balloon

- Enables diameter customization
- Improves device retention and trackability in tortuous anatomies



CBAS Heparin Surface for lasting thromboresistance



Expanded range of diameters and lengths



Proven leader with 20 years of peripheral stent graft clinical experience

- Employs the covered stent technology of GORE® VIABAHN® Endoprosthesis



GORE® VIABAHN® VBX Balloon Expandable Endoprosthesis

Flexible Strength. Proven Success.



See flexible strength in action:
goremedical.com/vbx-expands



W. L. GORE & ASSOCIATES, INC.

Flagstaff, AZ 86004

+65.67332882 (Asia Pacific)

00800.6334.4673 (Europe)

800.437.8181 (United States)

928.779.2771 (United States)

goremedical.com

* Delivery and outcomes data derived specifically from iliac occlusive disease clinical studies.
Bismuth J, Gray BH, Holden A, Metzger C, Panneton J; VBX FLEX Study Investigators. Pivotal study of a next-generation balloon-expandable stent-graft for treatment of iliac occlusive disease. Journal of Endovascular Therapy 2017;24(5):629-637.

Products listed may not be available in all markets.

CBAS is a trademark of Carmeda AB, a wholly owned subsidiary of W. L. Gore & Associates, Inc.
GORE®, VBX, VIABAHN®, and designs are trademarks of W. L. Gore & Associates.
© 2017 W. L. Gore & Associates GmbH AW3505-EN1 OCTOBER 2017