

FM 150 Vein trainer





Continuously adjustable massage intensity



Improves the venous blood circulation

Revitalising pressure massage

Promotes blood circulation and relieves heavy and tired legs

Looses tension

Air pressure massage with inflating and deflating air cushion

Continuously adjustable massage intensity

Alternate inflation of the leg cuffs

Incl. Safety STOP button

With deflation function - for emptying the air cushions quickly

Incl. timer function 10/20/30 min.

Easy operation with manual control switch

Versatile application – whether you're at home or on the go

Individual adjustment of the leg cuff using hook-and-loop fasteners

Storage pouch included

Medical device

Battery or mains adaptor driven

Mains operation:

Input: 100-240 V~; 50/60 Hz; 200 mA

Output: 6 V === 500 mA; 3 W

Incl. 4 x 1.5V battery AA

Product measurements (LxWxH): Leg cuff: 61 x 67 x 0.5 cm

Product weight: approx. 1.4 kg

3 year guarantee

Sales unit: 2 / Shipping carton: -

EAN no.: 4211125 64325 6

Item no.: 643.25



german

engineering

Physio Line

643.25_2020-01-14_03_DS_EN_BEU Subject to errors and changes