

INTERMEDIATE DRESSAGE TEST 140 (2020)

Interval between horses - 6 mins 30 Secs Arena 20m x 40m

To be ridden in a snaffle or simple double bridle		Max. Marks	
1.	Enter in Working Trot		
X	Halt Immobility Salute Proceed Collected Trot		
C	Track Right	10	
2.	MXK	Change the Rein in Medium Trot (Sitting or Rising).....	10
3.	KA	Collected Trot (Mark for Transitions at M & K).....	10
4.	After A	Turn onto the 3/4 line and from the FDK Line Shoulder In Left To the MGH Line then Track Left.....	10
5.	CH	Collected Trot	
	HXF	Change the Rein in Medium Trot (Sitting or Rising).....	10
6.	FA	Collected Trot (Mark for Transitions at H&F).....	10
7.	After A	Turn onto the 3/4 Line and from the FDK Line Shoulder In Right To the MGH line then track Right	10
8.	C	Collected Canter Right	
	CM	Collected Canter Right	
	MBF	Medium Canter.....	10
9.	F	Collected Canter Right and 1/2 10 Metre Circle return to the Track at B	10
10.	Between		
B & M	1/2 20 Metre Circle Left in Counter Canter		
EK	Counter Canter	10	
11.	K	Medium Walk	
	KAF	Medium Walk	10
12.	F	1/2 20 Metre Circle Left in Extended Walk.....	10
13.	KA	Collected Walk	
	A	Collected Canter Left Direct from Walk	10
14.	AF	Collected Canter	
	FBM	Medium Canter Left.....	10
15.	M	Collected Canter and 1/2 10 Metre Circle Return to the Track at B	10
16.	Between		
B & F	1/2 20 Metre Circle in Counter Canter		
EH	Counter Canter	10	
17.	HCMB	Collected Trot.....	10
18.	B	1/2 10 Metre Circle to X	
	XG	Collected Trot	
	G	Halt, Immobility Salute	10
		Leave the Arena at Free Walk at an appropriate place	
Collective Marks			
19.	Paces (freedom and regularity)	10	
20.	Impulsion (desire to move forward, elasticity of the steps, suppleness of the back, and engagement of the hind quarters)	10	
21.	Submission (attention and confidence; harmony, lightness and ease of the movements; acceptance of the bridle and lightness of the forehand)	10	
22.	Rider (position and seat of the rider, correct use of the aids and effectiveness of the aids.....	10	
		TOTAL 220	

N.B.In INTERMEDIATE TESTS, all trot work must be executed “sitting”