HAVEI™



interwoven



DESIGN STORY

BY QDESIGN

The Havei series was born from a simple but powerful question: What if healthcare spaces felt more like home?

In environments where patients often feel most vulnerable, the furniture around them plays a critical role in how safe, supported, and comforted they feel. Too often, these spaces default to sterile, utilitarian forms — prioritizing function at the expense of humanity.

We set out to change that.

Informed by site visits, caregiver interviews, and patient interactions, the design was shaped by direct observations of how people sit, move, and care within healthcare spaces. We sought to better understand health trends and learn how Havei can improve the delivery of care. We all agree that patient experience matters. And in healthcare, a recliner isn't just another place to sit. It's often holding someone in a fragile state, at one of the most vulnerable moments of their life. So, Havei needed to be a part of the care experience — to offer support, be safe and welcoming, and also ready for the unexpected.

Havei reimagines the healthcare recliner as more than a clinical tool. It's designed with intentional curves, soft transitions, and residential cues that invite rest without sacrificing performance. Every detail, from the silhouette to the mechanisms, was refined through several years of research and collaboration with nurses and facility planners. The result: a series that elevates both patient experience and caregiver support.



Design is the first caregiver.

PATIENT

Havei is designed to support patient comfort and dignity in every detail. Soft-touch arm caps, quiet recline, and a wide, gently contoured seat create a calming experience, whether for recovery, rest, or waiting. With carefully considered proportions, including seat height, recline angle, and reach zones, Havei helps patients feel secure and supported. Its approachable design brings warmth and ease to a variety of care settings.













CAREGIVER

Havei was built to accommodate frontline needs. We believe that when caregivers have the right tools, patients get better care. With a focus on access, safety, and ease of operation, Havei provides a supportive environment where comfort meets confidence. Specific consideration was taken to simplify functions such as transfers, mobility, and positioning, while keeping patient comfort in the forefront.



FAMILY & FACILITY

From family members, who may spend extended hours visiting or sleeping in the chair, to facility staff, Havei is engineered to be easy to move and operate. Smooth mobility and intuitive controls make Havei simple to reposition and adjust, whether it's for patient need, guest use, or it's purely operational. Clearly labeled controls and a QR code link to a video ensure ease of use for anyone.







24"W Patient Recliner Shown with Solid Surface Arm Caps, Fold-Down Tray, and Central Locking Casters



24"W Patient Recliner Shown with Infusion Arm Caps, Headrest, and Tech Arm



24"W Patient Recliner Shown with Transfer Arms



27"W Patient Recliner Shown with Egress Arm Caps and Central Locking Casters



27"W Patient Recliner Shown with Loop Arms



3-Position Recliner, High Back Shown with Egress Arm Caps and Metal Base



3-Position Recliner, Mid Back Shown with Casters and Fold-Down Tray



Shown with Solid Surface Arm Caps and Metal Base

Highlighted Features Please reference the price list for all features.

- 1. Contrasting upholstery available on the back and arms, which are the highest wear points, so that those areas can be upholstered with more durable materials.
- 2. The 90° back angle allows patients to sit upright for feeding and visiting without needing pillows placed behind them, which provides occupant comfort and supports caregiver access.
- 3. Dual transfer arms are designed to avoid interference with recline or footrest controls. They provide easy access for patient transfers (assisted or independent).
- 4. Optional heated seat and back soothes discomfort and supports healing.
- 5. The infinite-position recline helps patients stay comfortable for longer periods of time, while caregivers also benefit from dual recline controls and smooth mobility.

- 6. The utility hook is located at a convenient height and location for Foley bags.
- 7. The Trendelenburg function tilts the chair so that the feet are higher than the head, which is ideal for fast-paced environments where quick response is essential.
- 8. The central locking caster system features dual-wheel casters, steer assist, and full mobility modes, making transport easy and secure, even in tight spaces.
- 9. The pull out footrest provides additional leg support that is easy to stow away when not in use.
- 10. The USB port provides unobstructed access to charge devices while seated.



To learn more about the design of Havei, watch our interview with Qdesign.



