

Reference strains for Antibiotic Susceptibility Testing

Recommended strains by EUCAST / CA-SFM / VetCAST committees

CIP Number	Species
CIP 103429	<i>Staphylococcus aureus</i>
CIP 104340	<i>Streptococcus pneumoniae</i>
CIP 103214	<i>Enterococcus faecalis</i>
CIP 103570	<i>Haemophilus influenzae</i> II
CIP 70.2T	<i>Campylobacter jejuni</i>
CIP 76.24	<i>Escherichia coli</i> O6
CIP 76.110	<i>Pseudomonas aeruginosa</i>
CIP 102181	<i>Escherichia coli</i>
CIP 110856	<i>Staphylococcus aureus</i>
CIP 104676	<i>Enterococcus faecalis</i>
CIP 104604	<i>Haemophilus influenzae</i> I
CIP 103426T	<i>Mannheimia haemolytica</i>
CIP 110855	<i>Klebsiella quasipneumoniae similipneumoniae</i> (ex <i>K. pneumoniae</i>)

Reference strain for MIC determination of anti-tuberculous agents

CIP Number	Species
CIP 64.31T	<i>Mycobacterium tuberculosis</i>

Anaerobes

CIP Number	Species
CIP 77.16T	<i>Bacteroides fragilis</i>
CIP 103409T	<i>Clostridium perfringens</i> A