

705419 - Taille S

705429 - Taille M

705439 - Taille L

705449 - Taille XL



# **Technical product specification**

Product name	sempertip	Version / Index no:
Spec code	LLT-260YW-N-8CE / LLT-260RB-N-8CE	somportin Version E Echruny 2020 EN
Date of issue	03.02.2020	sempertip_Version F_February 2020_EN

#### **General information**

Type

protective glove for general use

Labelling

information printed on packaging

Shape

anatomical

Material

Natural Rubber Latex (NRL)

Colour

yellow or blue

Inside

flock lined

Outside

chlorinated

Cuff / surface

wave cut cuff / fish scale pattern

Shelf life

3 years

Available sizes

S (6-7) M (7-8) L (8-9) XL (9-10)

# Dimensions, physical properties and biocompatibility

Glove length

median 300 ± 10 mm

Minimum wall

not specified at finger

thickness

at palm  $0.35 \pm 0.05$  mm (single measured)

at cuff

not specified

Glove width

median S 100  $\pm$  5 mm, M 108  $\pm$  5 mm, L 115  $\pm$  5 mm, XL 123  $\pm$  5 mm

**Force at Break** 

not specified

**Tensile Strength** 

≥ 10 MPa (before aging) / ≥ 7.5 MPa (after aging) Elongation at Break ≥ 500% (before aging) / ≥ 375% (after aging)

Residual powder /

not specified

Powder content

## Performance requirements and inspection levels

Freedom from holes (Barrier)

**AQL 4.0** 

(sampling in accordance with ISO 2859-1, G-1)

Dimensions and physical properties

not specified

### Standards, guidelines & quality certificates

**Quality certification** 

ISO 9001

Conformity to regulations

- PPE Regulation (EU) 2016/425: Category I

- Regulation (EC) 1935/2004 on Food Contact Materials

Conformity to standards

EN 420, ISO 2859-1



# **Technical product specification**

Product name	sempertip	Version / Index no:
Spec code	LLT-260YW-N-8CE / LLT-260RB-N-8CE	sempertip_Version F_February 2020_EN
Date of issue	03.02.2020	Semperup_version i _i ebidary 2020_Liv

### Instructions and additional statements

Storage instruction

Store in original packaging in a dry and dark place at 10 °C to 30 °C. Refer to guidelines of storage of rubber products as described in ISO 2230:2002. Ensure that storage area is kept cool, dry and dust free, avoid ventilation and storage close to photocopy equipment. Copper ions discolour the glove. Protect gloves against ultraviolet light sources, such as sunlight and oxidizing agents. Storage above 30 °C will lead to accelerated aging and should be avoided.

Cautionary statement and ingredient information

This product contains natural rubber latex which may cause allergic reactions, including anaphylactic responses.

This product contains accelerators (Dithiocarbamate types, Zinc-mercaptobenzothiazol) not to be used in a hypersensitivity of these

substances.

For further information, a list of substances contained in the glove is available upon request.

**Reporting system** 

Medical device vigilance and reporting system

not applicable

A. Wöss Director J. Glantschnig | Head of Regulatory Affairs and Contract Management

Sempermed

L. Rieger

Head of Product Management

Remark

Replaces all previous versions.

All standards references refer to the date of document issue.