

Routine Cleaning

Remove ordinary dirt and smudges with mild soap and water. A 5:1 ratio of water to bleach solution may be used as a disinfectant. Rinse the surface with clean water after disinfecting. Dry with a soft, lint-free cloth or towel. The use of conditioners or protectants is not required nor recommended for use on Sta-Kleen PC upholstery – its cleanability is permanent and will not wear out.

Stain Removal

Upholstery protected with Sta-Kleen PC is resistant to most common stains. To keep your furniture looking new, stains such as ballpoint pen or permanent marker can be dry-erased with a clean, lint-free cloth. Gently rub the area until the stain has been removed. Wet or gooey stains such as food stains (e.g., ketchup or jelly) or topical stains (e.g., antiseptics, lotions and cream) should first be wiped off with a clean cloth or sponge, then follow the instructions above.

Stubborn Stains

If a ghost stain remains, apply a small amount of household rubbing alcohol (isopropyl alcohol) to a clean, lint-free cloth and rub the stain until it has been removed. Rinse with a clean, damp cloth.