

# Mini BlotBoy™ and BlotBoy™

## Benchmark 3-Dimensional Rockers



**BlotBoy™**  
12 x 12 in. platform



**Mini BlotBoy™**  
10.5 x 7.5 in. platform



*B3D1008  
shown with stacking platform*

Modeled after our popular BioMixer series of “nutating” shakers, these platform rockers provide the perfect 3-dimensional motion for all blot related applications. Their tilt angle and speed have been optimally set for gentle, but thorough mixing in gel trays, boxes and other flat vessels.

**Mini BlotBoy** and **BlotBoy** come standard with autoclavable flat mats and are ideal for blot washes and gel staining/destaining. Both models have been designed to conserve bench space. The **Mini BlotBoy** is supplied with a 10.5 x 7.5 inch platform while the larger **BlotBoy**'s 12 x 12 inch platform expands the useable work surface by 80%.

The **Mini BlotBoy** and **BlotBoy** can accept stackable platforms and handle loads up to 2kg. Both have maintenance free brushless motors and are safe for use in cold rooms and incubators.

- Optimized speed & tilt for gel blotting, washing & staining
- Tri-directional motion for thorough, gentle mixing
- Stacking platforms available
- Incubator & cold room safe

### Technical Data

Speed:	12 rpm*
Tilt Angle:	+/- 8 degrees
Platform Size:	
Mini BlotBoy:	10.5x7.5 in. / 27x19 cm
BlotBoy:	12x12 in. / 30x30 cm
Load Capacity:	1.6 kg / 3.5 lbs.
Operating Temp. Range:	+4°C to +65°C
Dimensions (WxDxH):	
Mini BlotBoy:	10.5 x 7.5 x 8 in.
BlotBoy:	12 x 12 x 8.3 in.
Weight:	2.0 kg / 4.4 lb
Warranty:	2 Years
Electrical:	115V AC, 60 Hz, 0.4A 230V AC, 50 Hz, 0.2A

\*Speed of 230V model is 10 rpm

### Ordering Information

<b>B3D1008*</b>	<b>Mini BlotBoy, 10.5 x 7.5 in. platform with flat mat</b>
B3D-STACK	Stacking Platform, flat mat, 10.5 x 7.5 in.
B3D-FLAT	Flat Mat, 10.5 x 7.5 in.
B3D-DIMPLED	Dimpled Mat, 10.5 x 7.5 in.
<b>B3D1308*</b>	<b>BlotBoy, 12 x 12 in. platform with flat mat</b>
BR1000-STACK	Stacking Platform, flat mat, 12 x 12 in.
BR1000-FLAT	Flat Mat, 12 x 12 in.
BR1000-DIMPLED	Dimpled Mat, 12 x 12 in.

\* To order a product in 230V, please add -230 to item number.

**Benchmark**   
Scientific

**Benchmark Scientific Inc.**  
PO Box 709  
Edison, NJ 08817  
Phone: 908-769-5555  
Fax: 908-222-1864  
Web: [www.BenchmarkScientific.com](http://www.BenchmarkScientific.com)  
Email: [Info@BenchmarkScientific.com](mailto:Info@BenchmarkScientific.com)