



# ACCURATE DATA IN LESS TIME

## **BRADYPRINTER I7100** VIAL LABEL APPLICATOR

Ensuring high quality data at higher volumes is a challenge. But, it's also 100% possible thanks to an automated solution like the **BradyPrinter i7100 Vial Label Applicator**. This device works with the BradyPrinter i7100 Industrial Label printer peel model, taking printed labels and applying them to vials and tubes. This is faster and more efficient than writing or hand-applying labels. It also reduces labeling time for lab techs, which saves both time and money. As daily volumes increase, so should your efficiency. And now, with the Vial Label Applicator, it finally can.

- ▣ **Efficient** — automated labeling saves time when used with the i7100 printer
- ▣ **Accurate** — allows precision tracking by eliminating error-prone, handwritten labels
- ▣ **Affordable** — pays for itself in up to 100 hours of use

Catalog #	Description
151156	BradyPrinter i7100 Vial Label Applicator
151752	BradyPrinter i7100 Industrial Label Printer peel version 300 dpi with Vial Label Applicator
151753	BradyPrinter i7100 Industrial Label Printer peel version 600 dpi with Vial Label Applicator

## Label materials

Catalog #	Color	Label width inch (mm)	Label length inch (mm)	Labels per roll	Recommended ribbon	Fits tube size*
<b>B-7425 Polypropylene — Economical polypropylene with moderate chemical resistance; ideal for labeling dry, room temperature specimen containers; must be applied to containers prior to being froze</b>						
THT-235-7425-2-AC	White	1.500 (38.1)	1.000 (25.4)	2000	R4300	1.8 mL-10 mL test tubes
THT-137-7425-2-AC	White	2.000 (50.8)	1.000 (25.4)	2000	R4300	1.8 mL-10 mL test tubes
THT-235-7425-2	White	1.500 (38.1)	1.000 (25.4)	2000	R4302	1.8 mL-10 mL test tubes
THT-179-7425-2	White	1.000 (25.4)	1.000 (25.4)	2000	R4302	0.5 mL- 2 mL microcentrifuge tubes
THT-137-7425-2	White	2.000 (50.8)	1.000 (25.4)	2000	R4302	1.8 mL-10 mL test tubes
THT-258-7425-2	White	1.750 (44.5)	1.500 (38.1)	2000	R4302	1.8 mL-10 mL test tubes
THT-136-7425-2	White	1.500 (38.1)	0.750 (19.1)	2000	R4302	1.8 mL-10 mL test tubes
THT-256-7425-2	White	1.300 (33.0)	0.600 (15.2)	2000	R4302	0.5 mL- 2 mL microcentrifuge tubes
THT-53-7425-2	White	2.000 (50.8)	0.500 (12.7)	2000	R4302	1.8 mL-10 mL test tubes
<b>B-499 Nylon Cloth — Nylon cloth with permanent adhesive for easy handling and conformability</b>						
THT-125-499-3	White	0.900 (22.9)	0.500 (12.7)	2000	R6000	0.5 mL- 2mL microcentrifuge tubes
THT-185-499-3	White	1.750 (44.5)	1.000 (25.4)	3000	R6000	1.8 mL-10 mL test tubes
THT-136-499-3	White	1.500 (38.1)	0.750 (19.1)	3000	R6000	1.8 mL-10 mL test tubes

\* Tube size may vary by brand

 **For more information, visit  
BradyID.com or call (888) 272-3946**

USA  
Customer Service: 1-888-272-3946  
Inside Sales: 1-888-311-0775  
BradyID.com

Canada  
Customer Service: 1-800-263-6179  
BradyCanada.ca

Mexico  
Customer Service: 1-800-262-7777  
Inside Sales: 1-800-262-7777 ext 177  
BradyLatinAmerica.com