



DSA-815 Notes on Memory Utilization

Date:04/09/12

The DSA815(-TG) comes with internal memory and a slot for an external USB memory sticks.

The information below provides an overview of the storage options available.

Internally, the “disk” of the 815(-TG) has two different directories.

C: used to store setups (.set)
for users 0 to user 7
size of 4 K

D: holds the other types of files

state files (* .sta)
numbers 0 – 9
1 K files

Trace files (.trc)
numbers 0 – 6
size 6K

Correction files (.cbl)
numbers 0 -9
size 1 K



Limit Files (.lim)
numbers 0-1
8 K

external Usb Memory(E:)

Mkr Table (.mrk, csv, bin)
1 K

Peak Table (csv, bin)
1 K

When the print button is pressed, a .bmp file (1.1 Meg) file is created

The data of the plot can also be written to an external USB stick as .CSV file (603 rows of data – 2 header lines – 601 data points).

More notes:

Firmware version 00.01.03 allows Mkr and Peak Tables to only be written to the external USB memory Stick.