


Western Blot Boxes

- For wholes, halves, strips or minis
- A perfect box for every blot
- Reduce antibody usage
- Hinged lid reduces evaporation
- Opaque black sizes available



B1200-7



B1200-13

B1200-17

Cat. No.	Description	Application	Qty
B1200-3	Micro-strip, 73 x 30 x 19mm	Various	5
B1200-5	Mini-rectangle, 73 x 51 x 32mm	Half mini-protein gels	5
B1200-7	Rectangle, 89 x 65 x 25mm	Mini-protein gels	5
B1200-7BK	Rectangle, 89 x 65 x 25mm, Black	(Bio-Rad Proteran 3 & similar)	5
B1200-9	Long strip, 95 x 30 x 16mm	Protean 3 (cut up)	5
B1200-11	Medium strip, 111 x 52 x 47mm	Various	5
B1200-13	Large, 117 x 89 x 28mm	Cambrex PAGER Gold, Novex, Might Small, SureLock mini	5
B1200-13BK	Large, 117 x 89 x 28mm, Black		5
B1200-15	Square, 116 x 116 x 32mm	Owl Penguin P9DS/Puffin	5
B1200-15BK	Square, 116 x 116 x 32mm, Black		5
B1200-17	Extra Large, 152 x 102 x 32mm	Bio-Rad Criterion, Amersham	5
B1200-19	XXL, 152 x 102 x 51mm	Bio-Rad Criterion, Amersham	5
B1300-8BK	Square, 85x85x28mm, Pull-Off Lid, Blk	Novex™ Mini Gels	10

No matter what size your blot is, there's an MTC box that fits it perfectly. With boxes that are correctly matched to common blot sizes, you can fit more blots on shakers/rockers while keeping reagent and antibody usage to an absolute minimum. Additionally, all corners are formed at a 90° angle and the bottom surface is uniformly flat, further reducing consumption of costly reagents and antibodies. (For example, popular mini-gels require only about 4mL.)

At a price that makes them truly disposable, all sizes are engineered from crystal clear polystyrene (except B1200-7BK, B1200-13BK, B1200-15BK & B1300-8BK, which are an opaque black that protects light sensitive samples). All sizes except B1300-8BK are also equipped with an easy-to-open hinged lid that reduces evaporation.

Western MultiBoxes™



B1212



B1203



B1204



B1206

New from MTC Bio, the MultiBox™ Western Blotting Box family. The popularity of our 6-place box has proven that researchers love the ease and convenience of using one box to process multiple gels. We've taken some of our most popular box sizes, and multiplied them, allowing end-users to compare results side by side in real time.

Cat. No.	Compartments	Box Dims. (LxWxH)	Compartment Dims. (LxWxH)	Qty
B1212	2 compartments	85 x 174 x 30mm	85 x 85 x 30mm	6
B1203	4 compartments	72 x 115 x 28mm	32 x 53 x 28mm	6
B1204	4 compartments	80 x 142 x 28mm	38 x 71 x 38mm	6
B1206	6 compartments	107 x 210 x 35mm	33 x 103 x 35mm	4