


PowerSpin™ DX Centrifuge

UNICO® is proud to present the PowerSpin™ DX Centrifuge

With 3 rotor options this centrifuge is a great answer for every day laboratory tasks.

- 24 place fixed rotor for 2-15ml tubes (#C8724)
- 6 place horizontal rotor for 2-15ml tubes (#C8706)
- 4 place horizontal rotor for up to 50ml tubes (#C8704)

The PowerSpin™ DX has variable speed 500 – 4000 RPM, to insure a clean and complete separation. Quiet and dependable, with a maintenance free motor, it can be used to process blood, urine, fecal and other specimens.



The all digital front panel with 3 pre-programmed options makes it easy to select most commonly run settings. The oversized display shows status, time, speed (RPM) and RCF (g). Precise speed control enables you to comply with CLIA recommendations for specific specimens.

- Zero RPM Safety Latching System prevents lid from being opened while rotor is spinning
- Quiet running, only 50 decibels at 3400 RPM
- 99 minute digital timer, 1 sec. increment
- 500-4000 RPM, 100 RPM increment

- 3 pre-programmed options for blood, fecal and urine
- Remembers last setting for quick repeat runs
- Adjustable Acceleration and Deceleration
- Stainless steel rotor chamber for easy cleaning
- Available in **110V** and **220V** - Contact **UNICO** for details



UNICO Professional Digital Non-Contact Photo-Tachometer

- Highly accurate, scientific tool
- Easy to operate
- Standard adhesive reflective tape (included)
- Aim and measure--RPMs instantly display on large LCD screen.
- Includes custom protective cushioned storage case
- Order pn# C800-75



Model C8724

- 24 place fixed angle rotor (pn#C8600-03 — rotor only)
- Designed to accept 10ml blood tubes and up to 15ml urine tubes
- 2ml tubes with small tube adaptor (pn# C800-04 — set of 2)
- 5ml tubes with short tube adaptor (pn# C800-03 — set of 10)



Model C8706

- 6 place horizontal rotor (pn#C8600-01 — rotor only)
- Specifically designed to keep tubes vertical at rest
- Designed to accept 10ml blood tubes and up to 15ml urine tubes
- 2ml tubes with small tube adaptor (pn# C800-04 — set of 2)
- 5ml tubes with short tube adaptor (pn# C800-03 — set of 10)



Model C8704

- 4 place horizontal rotor (pn#C8600-02 — rotor only)
- Designed to accept 50ml tubes
- With adaptors can accept most of the tubes used in laboratory
- 6 to 15ml tubes with adaptor (pn# C8604-02 — set of 2)
- 3 to 5ml tubes with adaptor (pn# C8604-03 — set of 2)
- 0.5 to 2ml tubes with adaptor (pn# C8604-04 — set of 2)

UNICO Model	Item	Speed	Max. RCF	Rotor	Tube Size	Tube Holdster 24 place
PowerSpin™ DX	C8724	500 – 4000 RPM	2508	24 place fixed	2-15ml	Optional
PowerSpin™ DX	C8724H	500 – 4000 RPM	2508	24 place fixed	2-15ml	Included
PowerSpin™ DX	C8706	500 – 4000 RPM	2508	6 place horizontal	2-15ml	Optional
PowerSpin™ DX	C8706H	500 – 4000 RPM	2508	6 place horizontal	2-15ml	Included
PowerSpin™ DX	C8704	500 – 4000 RPM	2508	4 place horizontal	50 ml	Optional
PowerSpin™ DX	C8704H	500 – 4000 RPM	2508	4 place horizontal	50 ml	Included