



# STIMULIABILITY APPROACH

## 7 Parts

1. Directly target nonstimulable sounds
2. Make targets the joint focus of attention
3. Associate speech sounds with hand/body motions
4. Associate speech sounds with alliterative characters
5. Encourage vocal practice
6. Ensure early success
7. Ensure successful communicative attempts



## Steps of Therapy

- Elicit ⅓ of stimulability probe: each consonant is elicited in isolation, then with a vowel in CV, VCV, and VC context. Repeat with 2 more vowels and the same consonant
- Review sounds and motions: model sound and associated movement
- Engage in play-based stimulability activities: 3 different 10-minute activities, SLP should model sound and encourage imitation
- Elicit palindrome probe