

S-FIT-80-002 | Functional knee brace

Indications

- Conservative treatment of advanced unilateral osteoarthritis
- Post-operative pressure relief & compartmental protection
- Instability of the cruciate ligaments
- Functional pre- or post-operative care of knee instabilities


Features & benefits

- Extension stops: 0°, 10°, 20°, 30° and 45°
- Flexion stops: 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90° and 120°
- Immobilisation at: 0°, 10°, 20°, 30° and 45°
- High quality, 100% aluminium frame
- Pressure relief of the medial or lateral hinge compartment by the 3-point principle off-loads the knee into varus or valgus position between 5-15°



This brace is part of our own FITLINE brand.


Orthopaedic braces and shoes, designed to obtain a maximum of comfort and to maintain a normal lifestyle. **Move - feel - live !**



Size: Measure circumference 15 cm above the mid-patella.

Size	Circumference
XS	31-37 cm
S	37-45 cm
M	45-51 cm
L	51-57 cm
XL	57-67 cm
XXL	67-74 cm

The correct choice of orthosis



Deformity	Varus deformity (bowed leg)	Valgus deformity (knocked knee)
Affected side	Medial	Lateral
Treatment	Medial relief	Lateral relief
Varus / Valgus adjuster joint	Lateral	Medial
Direction of force applied	Lateral to medial	Medial to lateral
Orthosis	ArtroFit Varus	ArtroFit Valgus

VARUS

Reference Left	Reference Right	Size
S-FIT-80-001-L1	S-FIT-80-001-R1	XS
S-FIT-80-001-L2	S-FIT-80-001-R2	S
S-FIT-80-001-L3	S-FIT-80-001-R3	M
S-FIT-80-001-L4	S-FIT-80-001-R4	L
S-FIT-80-001-L5	S-FIT-80-001-R5	XL
S-FIT-80-001-L6	S-FIT-80-001-R6	XXL

VALGUS

Reference Left	Reference Right	Size
S-FIT-80-002-L1	S-FIT-80-002-R1	XS
S-FIT-80-002-L2	S-FIT-80-002-R2	S
S-FIT-80-002-L3	S-FIT-80-002-R3	M
S-FIT-80-002-L4	S-FIT-80-002-R4	L
S-FIT-80-002-L5	S-FIT-80-002-R5	XL
S-FIT-80-002-L6	S-FIT-80-002-R6	XXL

