

# Introducing

the GORE® TAG® Conformable Thoracic Stent  
Graft with ACTIVE CONTROL System.

Experience new levels of control and achieve predictable  
patient outcomes with precise placement and trusted  
conformability—even in patients with complex  
anatomies, such as acute aortic angles.



**Controlled. Conformable. Predictable.**

Experience new levels of control at [goremedical.com/predictable](http://goremedical.com/predictable)



## GORE® TAG® Conformable Thoracic Stent Graft

CATALOGUE NUMBER	INTENDED AORTIC DIAMETER (mm)	PROXIMAL DIAMETER (mm)	DISTAL DIAMETER (mm)	ENDOPROSTHESIS LENGTH (cm)	RECOMMENDED GORE® DRYSEAL FLEX SHEATH SIZE (Fr)	OVERSIZING RANGE (%)	PARTIALLY UNCOVERED STENT LENGTH (mm)
TGM212110E	16–19.5	21	21	10	18	8–31	3
TGM262110E	19.5–24 / 16–19.5	26	21	10	20	8–33	4
TGM262610E	19.5–24	26	26	10	20	8–33	4
TGM282810E	22–26	28	28	10	20	8–27	4
TGM282815E	22–26	28	28	15	20	8–27	4
TGM312610E	24–29 / 19.5–24	31	26	10	22	7–33	4
TGM313110E	24–29	31	31	10	22	7–29	4
TGM313115E	24–29	31	31	15	22	7–29	4
TGM313120E	24–29	31	31	20	22	7–29	4
TGM343410E	27–32	34	34	10	22	6–26	5
TGM343415E	27–32	34	34	15	22	6–26	5
TGM343420E	27–32	34	34	20	22	6–26	5
TGM373710E	29–34	37	37	10	24	9–28	5
TGM373715E	29–34	37	37	15	24	9–28	5
TGM373720E	29–34	37	37	20	24	9–28	5
TGM404010E	31–37	40	40	10	24	8–29	6
TGM404015E	31–37	40	40	15	24	8–29	6
TGM404020E	31–37	40	40	20	24	8–29	6
TGM454510E	34–42	45	45	10	24	7–32	6.5
TGM454515E	34–42	45	45	15	24	7–32	6.5
TGM454520E	34–42	45	45	20	24	7–32	6.5



**W. L. GORE & ASSOCIATES, INC.**  
Flagstaff, AZ 86004

+65.67332882 (Asia Pacific)  
00800.6334.4673 (Europe)  
800.437.8181 (United States)  
928.779.2771 (United States)

[goremedical.com](http://goremedical.com)

Products listed may not be available in all markets.

GORE®, ACTIVE CONTROL, TAG®, and designs are trademarks of W. L. Gore & Associates.  
© 2017 W. L. Gore & Associates, Inc. AW2850-EN3 AUGUST 2017